EN

Conquest TTC™ Through-the-Channel

Lithotripter Cable

CS

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ENGLISH

INTENDED USE

This device is used to mechanically crush stones in the biliary duct when other methods of endoscopic removal have failed.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Used in conjunction with a side-viewing endoscope, Cook Endoscopy lithotripsy-compatible baskets and a Soehendra Lithotripter Handle.

Supplied non-sterile.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and mechanical lithotripsy.

Contraindications associated with biliary stone extraction include, but are not limited to: an ampullary opening inadequate to allow for the unimpeded passage of the stone and basket.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction to contrast or medication, respiratory depression or arrest, cardiac arrhythmia or arrest.

Potential complications associated with basket extraction include, but are not limited to: impaction of the object, localized inflammation, pressure necrosis.

PRECAUTIONS

This device is intended for use with a side-viewing endoscope only.

Assessment of the size of the stone and ampullary orifice must be made to determine the necessity of sphincterotomy prior to this procedure. The length of the intramural segment must be evaluated before determining the advisability of mechanical lithotripsy. In the event sphincterotomy is required, all necessary cautions, precautions and contraindications must be observed.

Cook Endoscopy wire guides are recommended for use with this product.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

Wire cutters must be available when using this device.

WARNINGS

Due to the mechanical pressure generated with this device, basket fragmentation in the common bile duct is a possibility that may require surgical intervention.

Only specified Cook Endoscopy Soft lithotripsy-compatible baskets are recommended for use with the Conquest TTC Lithotripter Cable. Refer to basket label.

Due to the varying compositions of biliary stones, stone fracture may not be possible. If the stone cannot be fractured, continued rotation of the handle may cause the basket wire to break, thus requiring surgical intervention.

INSTRUCTIONS FOR USE

1a. For baskets with removable handle.

Loosen the sliding handle from the basket wires by holding the thumb grip in one hand and turning the basket wire lock in the direction of the arrow. (See fig. 1) Hold the basket sheath firmly at the accessory channel port of the endoscope and slowly remove. **Note:** If desired, remove protective coil to expose basket wires. The basket wires can be used to facilitate basket control.

1b. For other lithotripsy compatible baskets.

Once the stone is completely captured in the basket, make a complete cut, using wire cutters, at the junction of the 10 French sheath and the 7 French sheath (See fig. 2) and dispose of basket handle per institutional guidelines for biohazardous medical waste.

2. Keeping the endoscope as straight as possible, remove the outer sheath of the basket slowly. As distal tip of the basket sheath clears endoscope elevator, gently lift the elevator to lock the basket wires in place. Then, continue removing the outer sheath. Dispose of outer sheath per institutional guidelines for biohazardous medical waste.

3. Introduce the lithotripsy cable over the basket wires until the wires exit the proximal end of the cable.

4. Grasp the basket wires exiting the cable. Release the elevator and continue to advance the cable until the distal tip reaches the basket within the duct.

5. If contrast injection is desired, attach adapter provided. **Note:** To attach adapter, insert the basket wires through the male Luer lock of the Tuohy-Borst adapter. Continue to slide the adapter over the basket wires until the male Luer lock engages the female Luer lock of the cable. Tighten the male Luer lock to prevent leakage during injection. The Tuohy-Borst adapter need not be attached if additional contrast is not required for the procedure.

6. Pass the basket wires through the male Luer lock on the lithotripsy handle; continue to pass the wires through the rotating-rod opening of the handle.

7. Connect the male Luer lock of the lithotripter handle to the female Luer lock of the adapter (if applicable). Tighten female Luer lock to prevent back flow during injection. Wrap the basket wires around the rotating rod to secure. (See fig. 3)

8. To fracture the stone, apply pressure by slowly rotating the handle. **Note:** Slow, constant pressure is often the most effective method of fracturing a stone. Sudden release of pressure may indicate the fracture of the stone, fracture of the basket, or loss of the stone from the basket. **Warning:** If basket fragmentation should occur, surgical intervention may be necessary. **Note:** The handle will only move in one direction unless the ratchet is released.

9. After the stone has been fractured, assess the need to remove the fractured stone from the duct.

10. Release the ratchet on the lithotripter handle to remove the basket. Examine the basket to confirm that it is not broken or damaged in any way, then dispose of the basket per institutional guidelines for biohazardous medical waste.

11. Dispose of the lithotripsy cable per institutional guidelines for biohazardous medical waste.

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Wilson-Cook Medical, Inc.
Cook Ireland Ltd.