EN

Reusable Hot Biopsy Forceps

CS

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ENGLISH

INTENDED USE

This device is used endoscopically in conjunction with monopolar electrosurgical current to obtain gastrointestinal mucosal tissue biopsies and for the removal of sessile polyps.

NOTES

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

We recommend this device be sterilized prior to first use following set points included in these instructions. Prior to each subsequent use, device should be cleaned and sterilized following cleaning and sterilization instructions provided in this booklet.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired biopsy or polypectomy site.

Contraindications to gastrointestinal mucosal biopsy and polypectomy include, but are not limited to: coagulopathy and insufficiently prepped bowel.

POTENTIAL COMPLICATIONS

Additional complications that can occur with endoscopic mucosal biopsy or polypectomy include, but are not limited to: transmural burns, thermal injury to the patient, explosion.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Endoscope must remain as straight as possible when inserting or withdrawing forceps.

Forceps cups must remain closed during introduction into, advancement through, and removal from endoscope. If cups are open, damage to forceps and endoscope may occur.

Exercising handle while forceps is coiled may result in damage to forceps.

Mucosal biopsy and polypectomy should be performed under direct endoscopic visualization.

Before using this device, follow recommendations provided by electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of patient return electrode. Ensure a proper path from patient return electrode to electrosurgical unit is maintained throughout the procedure.

Reusability of a device depends, in large part, upon care of device by user. Factors involved in prolonging life of this device include, but are not limited to: advancing through endoscope accessory channel in short increments, maintaining handle in alignment with endoscope accessory channel, gently withdrawing device from channel, avoiding loops and/or kinks in cable, avoiding use of excessive force, thorough cleaning following instructions included in this booklet.

This device should never be coiled in less than an 8 inch (20 cm) diameter. (See fig. 1)

WARNING

Failure to isolate tissue to be biopsied or polyp to be removed by pulling it away from mucosal wall may result in fulguration of normal mucosa and/or perforation.

PRODUCT INSPECTION

1. Uncoil forceps, open and close cups, verifying smooth handle operation and cup action.
2. With electrosurgical unit off, prepare equipment. Securely connect active cord to handle and electrosurgical unit.

INSTRUCTIONS FOR USE

1. Endoscopically visualize area to be biopsied or polyp to be removed.
2. With cups closed, insert forceps into accessory channel of endoscope.
3. Advance forceps in 1-2 cm increments until it is visualized exiting endoscope. **Note:** Keep end of forceps that is extending from accessory channel straight at all times. Allowing forceps to hang from accessory channel may cause damage to device.
4. Advance forceps to desired biopsy site or polyp, then open cups and advance into tissue to be biopsied or polyp to be removed.
5. Following electrosurgical unit manufacturer's instructions for settings, verify desired settings and activate electrosurgical unit. **Note:** Maximum rated input voltage for this device is 2.25 kVp-p.
6. Using slight pressure on handle, close forceps around tissue or polyp. Isolate tissue by gently pulling away from mucosal wall. **Caution:** When applying current, ensure metal tip of forceps does not come in contact with endoscope. Contact of forceps tip with endoscope may result in grounding, injury to patient and/or operator, as well as damage to endoscope and/or forceps.
7. Maintain gentle handle pressure to keep cups closed and gently withdraw forceps from site.
8. Continue to apply slight pressure on handle and withdraw forceps from channel. While withdrawing forceps from endoscope, wipe excess secretions from cable.
9. Prepare specimen for examination per institutional guidelines.
10. Upon completion of procedure, turn electrosurgical unit off and disconnect active cord from handle.
11. Wipe active cord with a damp cloth to remove all foreign matter. Store in a loose coil. **Note:** Wrapping tightly may damage device.

CLEANING INSTRUCTIONS

Please follow carefully to ensure proper reusability of this device.

Note: Cleaning must begin **immediately** after use of forceps. Allowing blood, feces or other body matter to dry on instrument will cause sticking of mechanism and eventual destruction of device. Device should be wiped clean of visible soil at point of use. If forceps cannot be cleaned immediately after use, they should be rinsed thoroughly with clean, running water, then **coiled loosely** and placed in a tray of warm water until they can be cleaned. Never allow dirty forceps to dry before cleaning. Wear rubber gloves to protect from infection or skin damage.

1. Rinse forceps thoroughly with clean, running water.
2. Fill an ultrasonic bath with warm water and enzymatic cleaner that has been mixed according to enzymatic cleaner manufacturer's instructions.
3. Draw detergent mixture into a clean plastic syringe and flush internal diameter of device through flush port. If internal diameter is blocked with debris, allow device to soak until debris is cleared.
4. Immerse forceps, cups first, into ultrasonic bath.
5. Ultrasonically clean device at 100 +/- 5°F (38 +/- 3°C) for not less than 10 minutes.
6. Draw detergent mixture into a clean plastic syringe and flush internal diameter of device through flush port again.
7. Test functionality of device by opening and closing forceps several times. If jaws respond smoothly and effortlessly and jaw mechanism is visually clean, proceed with step #8. If not, tip mechanism must be cleaned thoroughly using a nylon bristle brush until visible soils are removed and jaws are functional. (See fig. 2) Repeat ultrasonic cleaning, exercising jaws in cleaner until mechanism works efficiently. **Caution:** If forceps are not adequately cleaned, operation may be restrained and any subsequent disinfection or sterilization process may not achieve desired results.
8. Rinse both internal diameter of device and outer surface thoroughly with generous amounts of clean, running water, then dry. (See fig. 3) **Note:** After each procedure, forceps must be sterilized.
9. Flush internal diameter through flush port and rinse device in 70% ethyl alcohol to facilitate drying of inner drive cable.
10. Gently wipe external portions of forceps dry with a lint-free cloth. Air dry forceps by hanging for a minimum of three hours, or by using medical grade (oil free) forced air along cable and around handle and cups.
11. Lightly lubricate forceps hinges with medical grade silicone lubricant. (See fig. 4) **Note:** Proper lubrication of forceps after each use will help maintain performance characteristics of device.
12. Package device for sterilization following AAMI recommended practices.
13. Sterilize device using set points shown in table below.

AUTOCLAVE CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (24.8" Hg Vac)
Chamber Exposure Temperature	134°C (273.2°F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed usable chamber space.

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Wilson-Cook Medical, Inc.
Cook Ireland Ltd.