EN

Cunningham-Cotton Sleeve

CS

Návlek Cunningham-Cotton

DA

Cunningham-Cotton manchet

NL

Cunningham-Cotton huls

FR

Cathéter Cunningham-Cotton

DE

Cunningham-Cotton-Hülse

EL

Χιτώνιο Cunningham-Cotton

HU

Cunningham-Cotton hüvely

IT

Catetere Cunningham-Cotton

NO

Cunningham-bomullshylse

PL

Tuleja Cunningham-Cotton

PT

Cateter Cunningham-Cotton

ES

Manguito Cunningham-Cotton

SV

Cunningham-Cotton hölje

18966/1013

ENGLISH

INTENDED USE

This device is used for endoscopic cannulation of the ductal system.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Outer sleeve allows maintained access to biliary duct to facilitate papillotomy/sphincterotomy, biliary stent placement and/or nasal biliary drainage.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Cook Endoscopy wire guides are recommended for use with this device.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

INSTRUCTIONS FOR USE

1. Introduce and advance the device over a prepositioned wire guide.
2. Under fluoroscopic monitoring, advance the device to the desired location. **Note:** This device is equipped with two radiopaque bands on the dilator portion of the device and one radiopaque band on the sleeve portion of the device in order to fluoroscopically visualize the position of the catheter in the duct.
3. Remove the dilator portion of the device and the wire guide (if necessary), then fluoroscopically confirm proper position of the device.
4. Proceed with the exchange of other compatible Cook Endoscopy devices through the sleeve portion of the device.
5. Upon completion of the procedure, remove the device from the accessory channel of the endoscope and dispose of per institutional guidelines for biohazardous medical waste.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.