EN
3

Glo-Tip® ERCP Catheters

CS
3

Katetry Glo-Tip® pro ERCP

DA
4

**Glo-Tip® katetre til endoskopisk retrograd
kolangiopancreatikografi**

NL
5

Glo-Tip® ERCP-katheters

FR
6

Cathéters de CPRE Glo-Tip®

DE
7

Glo-Tip® ERCP-Katheter

EL
8

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
INTENDED USE

This device is used for endoscopic cannulation of the ductal system.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling pancreatic duct may cause pancreatitis.

INSTRUCTIONS FOR USE

1. Gently advance stylet into catheter. Lock stylet cap into catheter fitting. (*See fig. 1*)
Note: When using GT-5-4-3, remove and discard precurved wire from tip of catheter prior to introducing device into endoscope. (*See fig. 2*)
2. Introduce catheter into endoscope and advance in short increments until device is endoscopically visualized exiting scope.
3. Following cannulation, contrast may be injected to fluoroscopically confirm position of device by attaching a syringe to side port of catheter fitting.
4. A wire guide may be placed in duct to maintain position and facilitate exchange of compatible wire guided devices. Remove stylet from catheter and advance a compatible wire guide prior to removing catheter from endoscope. **Note:** Before each introduction of a wire guide into a device, wet wire guide in a sterile water or saline bath. A contrast filled catheter may make wire guide difficult to advance. Flush contrast from catheter with sterile water or saline before each introduction of the wire guide, observing necessary cautions if device is positioned in pancreatic duct.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.