EN

Biliary Stent

CS

Biliární stent

DA

Galdestent

NL

Biliaire stent

FR

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DE

Gallenstent

EL

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ENGLISH

INTENDED USE

This device is used to drain obstructed biliary ducts.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

This device must be stored in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

For 10 FR Sof-Flex® stents, lubricate the guiding catheter with water-soluble lubricant.

Retain the positioning sleeve (if any)* for use when introducing the stent flaps into the accessory channel.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Inability to pass wire guide or stent through obstructed area.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

This stent must be placed under fluoroscopic monitoring.

The positioning sleeve* or the pigtail straightener* is not intended for use in the accessory channel of the endoscope.

Care must be exercised when straightening the pigtail curls* in order to avoid kinking or breaking the stent.

The tapered tip end of the stent* or side holes* must be positioned in the common bile duct while the other end remains in the duodenum.

The longer flaps of the ST-2 Biliary Stent should be positioned in the duct while the shorter flaps remain in the duodenum. (See fig. 1)

This device should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Select the Cook stent introducer system (if not included) in the appropriate French size.

WARNING

For stents with radiopaque bands, MRI compatibility has not been established.

INSTRUCTIONS FOR USE

I. FOR NON-PIGTAIL STENTS:

1. Gently ensure full extension of all side flaps*.
2. Load positioning sleeve onto duodenal flap end of stent.

For 5 to 7 FR stents:

- 3.a. Introduce stent, tapered tip first, and positioning sleeve onto a pre-positioned wire guide.

For 8.5 FR and larger stents:

- 3.b. Remove Tuohy-Borst adapter from end of guiding catheter, then introduce guiding catheter into accessory channel over a pre-positioned wire guide.

- 3.c. Introduce stent, tapered tip first, and positioning sleeve onto guiding catheter and pre-positioned wire guide.

THEN REFER TO STEPS 4-7 BELOW

II. FOR PIGTAIL STENTS:

1. Advance pigtail straightener along shaft of stent to tapered tip end in order to straighten pigtail curl. (See fig. 2)
2. Introduce stent, tapered tip first, and pigtail straightener onto pre-positioned wire guide until straightener reaches second curl.
3. Advance pushing catheter over wire guide to advance pigtail stent into accessory channel.

THEN REFER TO STEPS 4-7 BELOW

4. Advance guiding catheter* and/or pushing catheter in 1-2 cm increments until stent is in desired position.
5. Fluoroscopically and endoscopically confirm desired stent position.

Note (for non-pigtail stents only): Replace Tuohy-Borst adapter onto end of guiding catheter*. Inject contrast, if desired, to fluoroscopically visualize stent position.

6. After confirming stent position, gently remove wire guide, then guiding catheter* from endoscope while maintaining position of stent with pushing catheter.
7. Gently remove pushing catheter from accessory channel.

Upon completion of procedure, dispose of introducer components per institutional guidelines for biohazardous medical waste.

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Wilson-Cook Medical, Inc.
Cook Ireland Ltd.