

EN 5	ERCP Catheter
CS 6	Katetr pro endoskopickou retrográdní cholangiopankreatografii (ERCP)
DA 8	Kateter til endoskopisk retrograd kolangiopancreatikografi
NL 10	ERCP-katheter
FR 12	Cathéter CPRE
DE 14	ERCP-Katheter
EL 16	Καθετήρας ERCP
HU 18	ERCP-katéter
IT 20	Catetere per ERCP
NO 22	ERCP-kateter
PL 24	Cewnik do endoskopowej cholangiopankreatografii wstecznej (ECPW)
PT 26	Cateter de CPRE
ES 28	Catéter de ERCP
SV 30	ERCP-kateter


- 1 Injection Port
 Injekční port
 Injektionsport
 Injectiepoort
 Orifice d'injection
 Injektionsport
 Θύρα έγχυσης
 Folyadékbefecskendezési
 nyílás
 Raccordo di iniezione
 Injeksjonsport
 Port do wstrzykiwań
 Orifício de injeção
 Orificio de inyección
 Injektionsport

- 2 Stylet
 Stilet
 Stilet
 Stilet
 Stylet
 Mandrin
 Στειλεός
 Mandrin
 Mandrino
 Stilett
 Mandryn
 Estilete
 Estilete
 Mandräng


2.

- | | |
|--|---|
| <p>1 Injection Port
 Injekční port
 Injektionsport
 Injectiepoort
 Orifice d'injection
 Injektionsport
 Θύρα έγχυσης
 Folyadékbefecskendezési nyílás
 Raccordo di iniezione
 Injeksjonsport
 Port do wstrzykiwań
 Orifício de injeção
 Orificio de inyección
 Injektionsport</p> | <p>2 Proximal Wire Port
 Proximální vodící port
 Proksimal kateterlederport
 Proximale voerdraadpoort
 Orifice pour guide proximal
 Proximaler Drahtport
 Εγγύς θύρα σύρματος
 Proximális drótnyílás
 Ingresso prossimale della guida
 Proksimal vaierport
 Proksymalny port dla prowadnika
 Orifício proximal do fio
 Acceso proximal de la guía
 Proximal ledarport</p> |
| <p>3 Wire Stop
 Zarážka drátu
 Kateterlederstop
 Draadstop
 Butée du guide
 Drahtstopper
 Στοπ σύρματος
 Dróthatároló
 Meccanismo di arresto della guida
 Vaiersperre
 Nakładka blokująca
 Batente do fio guia
 Tope de la guía
 Ledarspärr</p> | <p>4 IDE Port
 Port pro intraduktální výměnu (IDE)
 IDE-port
 IDE-poort
 Orifice IDE
 IDE-Port
 Θύρα IDE
 Intraductaliscser nyílás
 Accesso IDE
 IDE-port
 Port wymiany wewnątrzprzewodowej
 Orifício de IDE
 Acceso de IDE
 IDE-port</p> |
| <p>5 Radiopaque Band
 Rentgenkontrastní proužek
 Røntgenfast stribe
 Radiopake band
 Bande radio-opaque
 Röntgendichte Markierung
 Ακτινοσκοπική ταινία</p> | <p>Sugárfogó sáv
 Banda radiopaca
 Radioopakt bånd
 Opaska cieniodajna
 Banda radiopaca
 Banda radiopaca
 Röntgentätt band</p> |

3.

- | | | |
|---|---|--|
| <p>1 Injection Port
 Injekční port
 Injektionsport
 Injectiepoort
 Orifice d'injection
 Injektionsport
 Θύρα έγχυσης
 Folyadékbecskendezési nyílás
 Raccordo di iniezione
 Injeksjonsport
 Port do wstrzykiwań
 Orifício de injeção
 Orificio de inyección
 Injektionsport</p> | <p>2 Wire Control Port
 Kontrolní port pro drát
 Kateterlederkontrolport
 Draadcontrolepoort
 Orifice pour contrôle du guide
 Drahtkontroll-Port
 Θύρα ελέγχου σύρματος
 Drótvezérlő nyílás
 Apertura di controllo della guida
 Vaierkontrollport
 Port kontroli przewodnika
 Orificio de controlo do fio
 Acceso de control de la guía
 Ledarkontrollport</p> | <p>3 Slide Cover
 Posuvný kryt
 Glidehyster
 Schuifhoesje
 Protection coulissante
 Schiebeabdeckung
 Συρόμενο κάλυμμα
 Elcsúsztatható fedél
 Copertura scorrevole
 Glidehyster
 Ruchoma pokrywa
 Cobertura deslizante
 Cubierta deslizante
 Glidhölje</p> |
|---|---|--|

INTENDED USE

This device is used for endoscopic cannulation of the ductal system.

NOTES

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

The double lumen feature allows injection of contrast while providing wire guide access to desired duct.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

WARNINGS

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

SYSTEM PREPARATION

Attach the Wire Guide Locking Device to the endoscope accessory port (if applicable).

Flush catheter injection port with contrast to expel all air.

INSTRUCTIONS FOR USE

I. If using a non-Fusion device. Note: For best results, the wire guide should be kept wet.

1. Gently advance stylet into wire guide port. Lock stylet cap into catheter fitting. (See fig. 1)
2. Introduce catheter into endoscope and advance in short increments until device is endoscopically visualized exiting scope.
3. Following cannulation, contrast may be injected to fluoroscopically confirm position of device. **Note:** Catheter has a radiopaque tip to assist with fluoroscopic visualization.
4. A wire guide may be placed in wire guide port and into duct to maintain position and facilitate exchange of compatible wire guided devices. Remove stylet from catheter and advance a compatible wire guide prior to removing catheter from endoscope.

Fusion® Devices

II. If using Intra Ductal Exchange (IDE) port and a short wire. (See fig. 2)

1. Disengage Wire Stop from Proximal Wire Port (PWP) hub and ensure that distal tip of Wire Stop is proximal to IDE port.
2. Insert distal tip of wire guide into IDE port and advance until it is flush with distal tip of catheter.
3. Advance Wire Stop handle until it reaches PWP hub and connect Luer lock securely to PWP hub.
4. Advance tip of catheter through cap of Wire Guide Locking Device and continue advancing until it is endoscopically visible.

5. Disengage wire stop from PWP prior to cannulation and withdraw a minimum of 1 cm.
6. Cannulate common bile duct.
7. Following cannulation, contrast may be injected through injection port to fluoroscopically confirm position of device.
8. Prior to removing device, utilize reference marks on catheter to ensure that IDE port is within ductal system.
9. Fluoroscopically visualize radiopaque band at IDE port. Retract wire guide until radiopaque distal tip of wire guide passes band; a disengagement from wire guide lumen will occur.
10. Advance disengaged wire guide to maintain ductal access.
11. Lock wire guide into Wire Guide Locking Device and remove catheter from endoscope accessory channel.

III. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide. (See fig. 2)

Note: For best results, the wire guide should be kept wet.

1. Remove Wire Stop (if applicable).
2. Advance device over a pre-positioned wire guide, ensuring that wire exits catheter at PWP. Cannulate as desired.
3. Following cannulation, contrast may be injected through injection port to fluoroscopically confirm position of device.
4. Remove catheter using standard long wire exchange technique.

IV. If using an Omni device with Wire Control Port (WCP) and a short wire or long wire guide. (See fig. 3) **Note:** For best results, wire guide should be kept wet.

1. Advance tip of catheter through cap of Wire Guide Locking Device and continue advancing until it is endoscopically visible. Cannulate as desired.
2. Following cannulation, contrast may be injected through injection port to fluoroscopically confirm position of device.
3. Introduce a wire guide into WCP and advance into duct of choice, separating wire guide from catheter. (See figs. 4 & 5)
4. To withdraw catheter from the endoscope, ensure wire guide is locked in place, and then pull back on catheter, allowing wire guide lumen to separate from wire guide until metal band is visible at Wire Guide Locking Device and resistance is felt. (See fig. 6)
5. Unlock wire guide from Wire Guide Locking Device, completely remove catheter from wire guide, and re-lock wire guide.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.