COOK®

MEDICAL

CE

0123

EN
2

Fusion® Lithotripsy Basket

CS
4

Košík pro litotrypsi Fusion®

DA
5

Fusion® litotripsikurv

NL
7

Fusion®-lithotripsy basket

FR
8

Panier de lithotripsie Fusion®

DE
10

Fusion®-Lithotripsie-Körbchen

EL
12

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
ENGLISH

INTENDED USE

This device is used for endoscopic removal of biliary stones and foreign bodies.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to primary procedure to be performed in conjunction with removal of a stone or a foreign body.

If this device is to be used for biliary stone extraction, contraindications include an ampullary opening inadequate to allow for unimpeded passage of stone and basket.

Additional contraindications include, but are not limited to: biliary perforation.

POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy and ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction, respiratory depression or arrest, cardiac arrhythmia or arrest.

Other potential complications associated with use of an extraction basket include, but are not limited to: impaction of object, aspiration of foreign body, localized inflammation, pressure necrosis.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

If device is to be used for removal of biliary stones, assessment of stone size and ampullary orifice must be made to determine necessity of sphincterotomy. If sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

This device should never be coiled in less than an 8-inch (20 cm) diameter.

WARNINGS

If a stone cannot be removed endoscopically with this basket, a compatible inflation/lithotripter device and/or the Conquest TTC® Lithotripter Cable with Soehendra® Lithotripter Handle are needed to mechanically crush stone and aid in removal.

Basket wires may fragment and/or stone impaction may occur during lithotripsy, requiring surgical intervention.

SYSTEM PREPARATION

Flush injection port with contrast to expel all air.

INSTRUCTIONS FOR USE (See diagram at end of document.)

1. Unlock wire guide from Wire Guide Locking Device and advance tip of basket catheter over pre-positioned wire guide, ensuring wire guide exits catheter at Zip port.
2. Introduce device into endoscope accessory channel, relock wire guide into Wire Guide Locking Device. With elevator open, advance basket in short increments until endoscopically visualized exiting scope.
3. Endoscopically and fluoroscopically visualize extraction target and advance basket to appropriate position. **Note:** Fluoroscopic visualization may be enhanced by injecting contrast via the contrast injection port. Position basket sheath beyond target before extending basket. **Note:** Contrast injection may be performed with wire guide in place.
4. Confirm desired position of basket sheath relative to target.
5. Slowly advance basket out of sheath. Position basket around extraction target and slowly retract until entrapment is achieved.
6. Withdraw basket into sheath while maintaining entrapment. **Caution:** If difficulty is encountered when removing basket from duct, **do not use excessive force.** Moderate force may be applied by closing handle to manually fracture stone. If passage is still restricted, surgical intervention may be necessary.

Note: After removing all stones, previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices.

7. With elevator down and basket retracted into sheath, slowly withdraw device from endoscope accessory channel until catheter color change is visible and Zip port is exposed. Unlock wire guide from Wire Guide Locking device and perform wire guide exchange. **Caution:** Damage to wire guide lumen may result if Zip port is visible and device is pulled from accessory channel with wire guide locked in Wire Guide Locking Device. If wire guide lumen is damaged, discard device.
8. Relock wire guide into Wire Guide Locking Device and remove device from endoscope accessory channel.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.