EN

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19972/1013

ENGLISH**INTENDED USE**

This device is used to dilate strictures of the gastrointestinal tract, including strictures of the esophagus, pylorus, duodenum and colon.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Balloon can be inflated to three distinct diameters, as indicated on package and catheter tag.

This balloon is used in conjunction with an inflation device or manometer and may be filled with sterile water, sterile saline or up to a 1:1 contrast medium mixture consisting of contrast and sterile saline.

CONTRAINDICATIONS

Contraindications to dilation include, but are not limited to: uncooperative patient; asymptomatic rings, webs, or strictures; known or suspected perforation; acute corrosive injury of less than one week duration; failure to advance device through strictured area.

POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Do not pre-inflate the balloon.**WARNINGS**

During dilation, do not inflate balloon beyond the maximum indicated inflation pressure, as this could result in overextension or bursting of the balloon.

Do not advance the balloon dilator if resistance is encountered. Assess cause of resistance to determine if dilation should be reattempted.

INSTRUCTIONS FOR USE

1. Attach the balloon to a 60 ml (cc) inflation device with gauge to monitor balloon pressure.
2. To facilitate passage through the endoscope, apply negative pressure to the device.
3. Remove the protective sheath from the balloon. **Note:** Wire guide can be removed from balloon catheter prior to use and replaced with a standard .035" wire guide to aid in bridging difficult strictures or device can be back loaded over a pre-positioned wire guide. **Note:** Prior to advancing device over a wire guide, flush wire guide lumen.
4. Apply a lubricating agent to the balloon to facilitate passage through the endoscope accessory channel.
5. Maintain balloon deflation with negative pressure and introduce into endoscope accessory channel, advancing in short increments until balloon is completely visualized endoscopically. Positioning indicator (*fig. 1*) on catheter indicates balloon clearance from endoscope accessory channel. **Note:** If using wire guide, move slide to "on" position and advance wire guide to desired location. Approximate wire guide advancement can be determined using markers located at 5-cm increments on wire guide. **Note:** Fluoroscopy is recommended if advancing wire guide without direct visualization.
6. Monitoring endoscopically, advance device until balloon is in desired position within stricture. **Note:** When balloon and wire guide are in desired location, wire guide can be secured to maintain its position by moving slide to "off" position.
7. Inflate balloon to pressure corresponding to smallest balloon diameter and maintain until desired dilation is achieved. To achieve increasingly larger balloon diameters, increase pressure as indicated on catheter tag.
8. To deflate balloon, apply negative pressure and remove all fluid from balloon while observing balloon endoscopically.
9. Remove deflated balloon from accessory channel. **Caution:** A compromised balloon may prohibit removal from endoscope accessory channel. Removal of endoscope along with compromised balloon may be required.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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Wilson-Cook Medical, Inc.
Cook Ireland Ltd.