
EN 2	Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube Instructions for Use
CS 6	Balónková retenční gastrostomická výživová sonda Entuit™ Thrive Návod k použití
DA 11	Entuit™ Thrive ballonfastholdt gastrostomi-ernæringssonde Brugsanvisning
DE 16	Entuit™ Thrive Gastrostomie-Ernährungssonde mit Retentionsballon Gebrauchsanweisung
EL 22	Σωλήνας σίτισης γαστροστομίας με μπαλόνι συγκράτησης Entuit™ Thrive Οδηγίες χρήσης
ES 27	Sonda de alimentación por gastrostomía de retención con balón Entuit™ Thrive Instrucciones de uso
FR 33	Sonde d'alimentation pour gastrostomie à ballonnet de rétention Entuit™ Thrive Mode d'emploi
HU 38	Entuit™ Thrive ballonnal megtartott gastrostomiás táplálósonda Használati utasítás
IT 44	Sondino gastrostomico con sistema di ritenzione a palloncino Entuit™ Thrive Istruzioni per l'uso
NL 49	Entuit™ Thrive-gastrostomie-voedingssonde met fixatieballon Gebruiksaanwijzing
NO 55	Entuit™ Thrive ernæringssonde for gastrostomi med ballong og festemekanisme Bruksanvisning
PL 60	Gastrostomijny cewnik żywieniowy utrzymywany balonem Entuit™ Thrive Instrukcja użycia
PT 66	Tubo de alimentação por gastrostomia com balão de retenção Entuit™ Thrive Instruções de utilização
SV 71	Entuit™ Thrive gastrostomi med retentionsballong och matningsslang för tillförsel Bruksanvisning

ENGLISH

ENTUIT™ THRIVE BALLOON RETENTION GASTROSTOMY FEEDING TUBE

Instructions For Use:

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Contents of package:

1 - Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube

DEVICE DESCRIPTION

The Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube is a sterile device consisting of a silicone balloon, bi-lumen shaft, three-port funnel and bolster that allows for proper retention during enteral feeding, medication administration and decompression. The device contains one valve that allows for inflation and deflation of the silicone balloon.

INTENDED USE/INDICATIONS FOR USE

The Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube is intended to provide gastric access for enteral feeding, medication administration and decompression through an established gastrointestinal stoma tract.

The Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube is indicated for use in percutaneous placement of an enteral feeding tube in adult and pediatric patients that require enteral feeding, medication administration or decompression through an established gastrointestinal stoma tract.

CONTRAINDICATIONS

Device placement is contraindicated for patients with evidence of granulation tissue, infection, and/or irritation in the stoma site.

WARNING

- Inspect package integrity before use. Do not use device if package is damaged or if the sterile barrier has been compromised. Do not use if labeling is incomplete or illegible.
- Do not resterilize or reprocess this medical device as this may have an adverse effect on the known characteristics of the structural integrity, performance and biocompatibility of the device.
- For single patient use only. Do not reuse this medical device as this may increase the risk of contamination leading to transmission of infectious diseases which has the potential of resulting in patient injury, illness or death.
- **FOR ENTERAL USE ONLY.** The device is intended to connect to enteral giving sets for enteral feeding, and syringes for feeding, medication and inflation of balloon.
- This device has **KNOWN** misconnections with connectors found in the following medical devices/healthcare applications:
 - Intravascular devices;
 - Hypodermic applications;
 - Breathing systems and driving gas devices;
 - Urethral/urinary devices;
 - Limb cuff inflation devices;
 - Neuraxial devices
- Do **NOT** use this product in the vasculature.
- Once used, dispose of packaging and package contents in accordance with healthcare institution guidelines and/or local government policy.

PRECAUTIONS

- This device should only be used for its intended purpose by or under the supervision of trained healthcare professionals with a comprehensive understanding in clinical principles, procedures and risks associated with percutaneous placement of enteral feeding devices. It is recommended to adhere to the instructions for use provided with this device, the enteral feeding instructions for the healthcare facility and instructions recommended by physicians.

COMPLICATIONS/ADVERSE EVENTS

Complications and adverse events associated with placement and use of a balloon gastrostomy tube include, but are not limited to:

- Aspiration, reflux, sepsis, ascites, bleeding, peritonitis and perforation
- Granulation tissue, pressure necrosis and ulcers
- Irritation and infection such as redness, edema or purulent drainage
- Severe gastroesophageal reflux or diffuse inflammatory, infectious, or neoplastic disease involving the walls of the abdomen or anterior stomach
- Gastrointestinal obstruction and proximal small bowel fistulae
- Tube clogging, kinking, malposition, migration, leakage and unintended tube dislodgement.

HOW SUPPLIED

Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube is provided sterile with one (1) Instructions for Use for each device.

HANDLING AND STORAGE

Store in a cool, dry, and secure environment to prevent damage.

DEVICE PREPARATION

Warning: *An inappropriately sized tube may cause leakage of gastric contents, necrosis, granulation tissue, sepsis, associated sequelae and/or buried bumper syndrome.*

1. Select the appropriately sized Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube based on the diameter of patient's stoma.
2. Inspect package integrity before use. Do not use device if package is damaged or if the sterile barrier has been compromised.
3. Remove device from package and inspect contents for damage. Do not use contents if damaged.
4. Inflate the Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube balloon to recommended fill volume (*printed on the device*) with sterile or distilled water using a sterile syringe.
Caution: *Only use sterile or distilled water for balloon inflation. Do not use air, saline, feeding formula, medication, or radiopaque contrast for balloon inflation.*
5. Verify balloon integrity by visually inspecting inflated balloon for symmetry and by gently squeezing the balloon to determine if the balloon is leaking. Do not use device if any suspected characteristics that may affect the performance of the balloon are noted.
6. Deflate balloon by reinserting syringe and removing all water from balloon. Do not attempt to use device if balloon cannot be deflated.
7. If necessary, lubricate the tip of the device with a water soluble lubricant for ease of insertion. Sterile gauze pads may be used for ease of tip lubrication.

Note: *Do not use petroleum jelly or mineral oil for tip lubrication.*

DEVICE PLACEMENT

Caution: Prior to device placement, a Gastropexy must have been performed to affix the stomach wall to the abdominal wall, an incision site must be identified for creation of stoma tract, then the stoma must be dilated and finally measured, in order to ensure patient safety and comfort.

Caution: Confirm that the balloon of the Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube will not occlude the Pylorus when inflated, specifically, when device will be inserted into infants or children.

Warning: Do not attempt to affix the stomach to the abdominal wall using the Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube. Excessive tension, tissue necrosis and inadvertent tube removal may result in separation of the gastric and abdominal wall.

1. Select and prepare the appropriately sized Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube in accordance with the "DEVICE PREPARATION" instructions in the section above.
2. Carefully insert device through the established stoma and into the stomach lumen.
Note: Maintain insertion angle perpendicular to the surface of the skin.
3. Verify that device is inside the stomach lumen using fluoroscopic procedures to determine position of device using the radiopaque mark located on the tip of the device or by following the instructions provided in the "VERIFICATION OF TUBE PLACEMENT" section listed below.
4. Inflate the balloon to recommended fill volume (*printed on the device*) with sterile or distilled water using a sterile syringe.
Caution: Only use sterile or distilled water for balloon inflation. Do not use air, saline, feeding formula, medication, or radiopaque contrast for balloon inflation.
5. Gently pull device toward abdomen until balloon is up against the inner stomach wall.
6. Slide the external bolster down to the skin.
7. Verify that there is no leakage of gastric contents around the stoma site and that there is not excessive pressure between the balloon and external bolster.
Warning: Inflating balloon with less volume than the recommended volume may result in leakage of gastric contents and/or device migration. Inflating balloon with more volume than the recommended volume may decrease balloon performance.
8. Clean any remaining contents around the stoma site with mild soap and water.

VERIFICATION OF TUBE PLACEMENT

Warning: Verify proper tube placement before feeding administration.

1. Draw 10mL of water into a catheter tip syringe and insert into feeding port of the funnel.
2. Gently pull back on syringe plunger until stomach contents (*typically yellow or clear unless food is in patient's stomach*) are aspirated and visible in the feeding lumen.
Note: The presence of gastric contents in the feeding lumen confirms correct positioning and that feeding tube is inside the stomach lumen.
3. Flush contents with the 10mL of water in the syringe.
4. Remove syringe from device.

FEEDING ADMINISTRATION

Caution: Only begin feeding administration after confirmation of proper feeding tube placement and patency. Administration of feeding should be performed according to instructions provided by physician.

Caution: Do not attempt to administer solid nutrients through feeding lumen.

Warning: This device is intended **FOR ENTERAL USE ONLY**. Do **NOT** attempt to connect with the devices mentioned in the main "WARNING" section. Device funnel is printed "FEED ONLY" and "ENTERAL MEDS ONLY" to further clarify that this device is intended for enteral applications only.

1. Remove the plug from the feeding port.
2. If feeding with a **catheter tip syringe**:
 - 2.1 Fill catheter tip syringe with water and connect to the feeding port.
 - 2.2 Prime the feeding lumen with water using the water-filled syringe and inspect for leakage.
Note: The device may need to be replaced if leakage is present.
 - 2.3 Disconnect syringe. Remove remaining water, if any, and plunger from syringe.
 - 2.4 Reconnect the syringe without the plunger to the feeding port.
 - 2.5 Pour feeding formula into syringe.
 - 2.6 Raise or lower syringe to increase or decrease feeding flow rate.
 - 2.7 Flush the feeding tube with the prescribed amount of water when feeding is complete.
 - 2.8 Remove syringe.
3. If feeding **continuously**:
 - 3.1 Attach enteral giving set connector to the feeding port.
 - 3.2 Fill feeding container with feeding formula.
 - 3.3 Allow for filling of feeding lumen. Inspect for leakage.
Note: The device may need to be replaced if leakage is present.
 - 3.4 Connect enteral giving set to pump.
 - 3.5 Set desired flow rate and begin feeding.
 - 3.6 If feeding continuously for long periods of time, flush feeding tube with prescribed amount of water every 4-6 hours.
 - 3.7 Disconnect enteral giving set from pump and feeding port when feeding is complete.
 - 3.8 Flush the feeding tube with prescribed amount of water.
4. Reconnect plug to feeding port.

MEDICATION ADMINISTRATION

Caution: Only administer medication after confirmation of proper feeding tube placement and patency. Administration of medication should be performed according to instructions provided by physician.

Caution: Use medication in liquid form whenever possible. If only solid medication is available, seek advice from a physician to determine if it is safe to crush solid medication and mix with water.

Warning: Do not crush medication that has an enteric coating or attempt to mix medication with feeding formula.

Warning: This device is intended **FOR ENTERAL USE ONLY**. Do **NOT** attempt to connect with the devices mentioned in the main **"WARNING"** section. Device funnel is printed **"FEED ONLY"** and **"ENTERAL MEDS ONLY"** to further clarify that this device is intended for enteral applications only.

1. Before administering medication, flush feeding tube with prescribed amount of water using a catheter tip syringe.
2. Remove plug from the medication port adaptor.
3. Medication may be administered through the medication port or through the medication port adaptor using the appropriate syringe.
Note: If safe, solid medication may be pulverized into a fine powder and dissolved with water before being administered.
4. After medication administration is complete, flush the route used for medication with prescribed amount of water.

TUBE REMOVAL

Caution: Never attempt to remove feeding tube unless trained by physician or healthcare provider. Before removing feeding tube, verify that the feeding tube may be replaced at bedside.

1. Obtain equipment and supplies required for removal of device. Clean hands and use powder-free gloves before touching stoma site.
2. Gently rotate the feeding tube to ensure the tube moves freely and easily.
Warning: Do not apply excessive force when removing device. If resistance is felt, apply water soluble lubricant and manipulate tube by pushing and rotating until tube moves freely and easily. Contact physician if not able to remove feeding tube.
3. Remove any contents inside the balloon using a syringe.
4. Pull device out of stoma while simultaneously applying counter-pressure to the patient's abdomen.
5. Clean stoma site with mild soap and water. Allow area to dry.
6. Dispose of device following the guidelines for biohazardous waste provided by the healthcare institution.
7. If a replacement device is required, follow instructions provided in "DEVICE REPLACEMENT" instructions below.
Caution: Stoma tract may begin to close in the first hour following device removal.

DEVICE REPLACEMENT

1. Clean stoma site with mild soap and water. Allow area to dry.
2. Select the appropriate size for replacement device, inspect device prior to use and determine if water-soluble lubrication is needed following the instructions provided in the "DEVICE PREPARATION" section listed above.
3. Carefully insert device through the established stoma and into the stomach lumen.
Note: Maintain insertion angle perpendicular to the surface of the skin.
4. Verify that device is inside the stomach lumen using fluoroscopic procedures to determine position of device using the radiopaque mark located on the tip of the device or by following the instructions provided in the "VERIFICATION OF TUBE PLACEMENT" section listed above.
5. Inflate the balloon to recommended fill volume (*printed on the device*) with sterile or distilled water using a sterile syringe.
Caution: Only use sterile or distilled water for balloon inflation. Do not use air, saline, feeding formula, medication, or radiopaque contrast for balloon inflation.
6. Gently pull device toward abdomen until balloon is up against the inner stomach wall.
7. Slide the external bolster down to the skin.
8. Verify that there is no leakage of gastric contents around the stoma site and that there is not excessive pressure between the balloon and external bolster.
Warning: Inflating balloon with less volume than the recommended volume may result in leakage of gastric contents and/or device migration. Inflating balloon with more volume than the recommended volume may decrease balloon performance.
9. Clean any remaining contents around the stoma site with mild soap and water.

TUBE MAINTENANCE

1. The stoma site should be inspected and cleaned regularly. Cleaning should be performed with mild soap and water, remove moisture afterwards. Notify physician if signs of infection or irritation.
2. Clean residual contents from the external area of the feeding tube with mild soap and water as needed, remove moisture afterwards.
3. Rotate device daily to allow for adequate air circulation. Consult with physician if device does not rotate freely.
4. In order to prevent feeding tube clogging, always flush feeding lumen with prescribed amount of water after feeding administration, every 4-6 hours during continuous feeding, before administering medication and after medication administration is complete.

Note: The amount of water used for flushing will depend on patient's need, clinical condition and if patient is an infant or adult.

5. Assess the feeding tube daily for damage, clogging or abnormal discoloration.
6. Check balloon volume weekly, or as recommended by the physician, and compare to initial volume used to inflate balloon in order to determine balloon integrity.
7. Replace device as or if necessary following the "DEVICE REPLACEMENT" instructions in the section above.

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ACRA Regulatory Services Pty Ltd.