COOK

MEDICAL

EN
5

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









9

10

ARNDT AIRWAY EXCHANGE CATHETER SET WITH RAPI-FIT® ADAPTERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Arndt Airway Exchange Catheter Set (AAE) contains a wire guide and a 14 French, taper-tipped, radiopaque exchange catheter with centimeter marks to help facilitate the exchange of a supraglottic airway device (SAD), such as a laryngeal mask airway (LMA) or endotracheal tube (ETT). **Does not contain natural rubber latex.**

The device consists of the following components:

- Radiopaque exchange catheter
- Wire guide
- Double swivel connector
- Rapi-Fit adapter – 15 mm connector
- Rapi-Fit adapter – Luer lock connector

INTENDED USE

The Arndt Airway Exchange Catheter Set is intended for exchange of a supraglottic airway device (SAD) to an endotracheal tube (ETT) under bronchoscopic assistance, and for ETT exchange in adult and pediatric patients.

The product may be used for emergency, urgent, and elective airway management.

The 14 French catheter is recommended for placement of an endotracheal tube or a supraglottic airway device with an inner diameter of 5 mm or larger.

When used for high-pressure oxygenation with a Luer lock connector, the 14 French catheter is recommended for patients older than 12 years of age.

CONTRAINDICATIONS

None known

WARNINGS

- Advancing the catheter beyond the carina may introduce additional procedural risks and trauma.
- Attention should be paid to the insertion depth of the Arndt Airway Exchange catheter (AAE) into the patient's airway and correct tracheal position of replacement endotracheal tube (ETT). Markers on the AAE catheter refer to distance from the distal tip of catheter.
- Take care to avoid injuring the epiglottis, glottis, trachea, bronchi, or lung parenchyma, and to avoid perforating the sinus pyriformis while using this device.
- Ensure proper sizing of the AAE within the SAD or ETT.
- To avoid barotrauma, ensure that the tip of the exchange catheter is always above the carina, preferably 2-3 cm.
- Use of the Rapi-Fit adapter for oxygenation may be associated with a risk of barotrauma.
- Use of an oxygen source should be considered only if the patient has sufficient egression of the insufflated gas volume.
- If a high-pressure oxygen source is used for insufflation (e.g., jet ventilator), begin at a lower pressure (i.e., 5 psi) and work up gradually. Rising chest wall, pulse oximetry, and oral air flow should be carefully monitored.
- Ensure that the Rapi-Fit Adapter is securely connected to the AAE prior to oxygen delivery. Failure to properly secure the adapter to the AAE may result in hypoxia and serious adverse events.
- High-pressure oxygenation with a Luer lock connector should only be used in patients older than 12 years of age. If used in patients 12 years old or younger, the maximum airway pressure may be higher than 28 cm H₂O.

PRECAUTIONS

- This product is intended for use by clinicians trained and experienced in airway management techniques. Standard techniques for placement and exchange of SADs or ETTs should be employed.
- This product is not intended for intravascular use.
- This product should not be used in replacement of SADs or ETTs whose inner diameter is smaller than 5 mm.
- The use of capnography is highly recommended to confirm intubation.
- A pulse oximeter is recommended during the use of this device.
- It is recommended that a sterile lubricant or water be applied to the AAE, bronchoscope, and new SAD or ETT prior to their introduction into the patient.
- When using an SAD, ensure the cuff is properly inflated (if applicable) and breath sounds are present.
- Ensure that wire guide position is maintained during placement of the AAE.
- Patients with airway injuries may have significant spinal cord injury. When possible, the cervical spine should be immobilized before beginning the procedure.

- This device is not recommended for definitive, long-term airway management.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Barotrauma
- Pneumothorax
- Hypoxia
- Perforation of the trachea, bronchi, or lung parenchyma
- Injury to the epiglottis
- Airway bleeding
- Catheter dislodgment or migration
- Failed endotracheal tube placement

INSTRUCTIONS FOR USE

1. After ensuring that the SAD or ETT is correctly positioned within the trachea, insert the double swivel connector between the SAD or ETT and the ventilation source. **(Fig. 1) NOTE:** Ventilation should be continued during the exchange process, through the double swivel connector. Discontinue ventilation only during the actual exchange over the AAE.
2. Advance the bronchoscope into the SAD or ETT through the double swivel connector. If a SAD is in place, the bronchoscope should be advanced through the vocal cords and into the trachea. **(Fig. 2)** If an ETT is in place, the bronchoscope should be advanced directly into the trachea through the ETT.
3. Advance the soft, flexible end of the wire guide through the working port of the bronchoscope until it is seen to exit the bronchoscope into the patient's trachea. Advance the wire guide to the level of the carina. **(Fig. 3)**
4. Carefully remove the bronchoscope. Leave the wire guide in place in the trachea, exiting the bronchoscopic port. **(Fig. 4)**
5. Maintaining wire guide position, advance the AAE over the wire guide through the double swivel connector and through the SAD or ETT into the trachea. **(Fig. 5)**
6. Properly position the AAE within the SAD or ETT by aligning the appropriate centimeter mark on the AAE with the corresponding centimeter mark (if present) on the SAD or ETT.
7. Fully deflate the cuff of the SAD (if applicable) or ETT while maintaining position of the AAE. Remove the SAD or ETT, leaving the AAE and wire guide in place. **(Fig. 6)**
8. Advance the new ETT over the AAE and wire guide into the patient's airway. **(Fig. 7) NOTE:** It is recommended that a sterile lubricant be applied to the tip of the ETT prior to advancing the tube.
9. Remove the AAE including the wire guide and inflate the balloon cuff of the new ETT. **(Fig. 8)** Reestablish ventilatory device use and secure the new ETT in place. Verify the ETT's position using standard methods (e.g., capnography, breath sounds and chest x-ray).

Use of the Rapi-Fit® Adapter

Rapi-Fit adapters should only be used when oxygen requirements are high and intubation is unsuccessful. Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume. If an oxygen source is used for insufflation, begin at a lower pressure and work up gradually. Observe the chest for outward and inward movements to confirm oxygen insufflation and egression. Pulse oximetry and oral air flow should be carefully monitored as well. In cases of upper airway obstruction, gas discharge from the patient's lungs may require more time.

1. To attach the Rapi-Fit Adapter, position the adapter on the catheter, then push the white collar forward and lock into position. **(Fig. 9)**
2. To remove the adapter, pull the white collar back to release, and then remove from the catheter. **(Fig. 10)**

Catheter Oxygenation

The Arndt Airway Exchange Catheter Set is designed for both positive airway pressure ventilation (Rapi-Fit adapter – 15 mm connector) and jet ventilation (Rapi-Fit adapter – Luer lock connector) for both adult and pediatric patients. In the table below, the delivered minute volume and the measured airway pressure are given for jet ventilation in adult and pediatric patients with healthy lung tissue.

Real Part Number (RPN)	Patient Subgroup and Age Range	Delivered Minute Volume ¹ (L/min)	Measured Airway Pressure ¹ (cm H ₂ O)	
			Mean	Maximum
C-CAE-14.0-70-LMA-01	Adolescent >12 years through 21 years	6.0	16.6	17.3
	Adult >21 years	10.2	13.3	14.4

¹ The test conditions were used in an active model. See Testing Conditions table for additional details.

The presented numbers are mean values. The following testing conditions were used in an active model mode, tested with ASL 5000, IngMar Medical, Ltd., for each patient subgroup.

Testing Conditions

Input pressure was set to 50 psi for both patient type models.

Patient type	Model Body Weight (kg)	Inspiratory time (s)	Expiratory time (s)	Breaths per minute	Resistance (cm H ₂ O/L/s)	Lung Compliance (mL/cm H ₂ O)
Adolescent (> 12 years through 21 years)	50	0.9	2.1	20	5	40
Adult (> 21 years)	80	1.0	4.0	12	3	100

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from clinicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2018
Bloomington, IN 47404 U.S.A.