EN
3

Arndt Airway Exchange Catheter Set
Instructions for Use

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SV
20

Arndt-utbyteskateterset
Bruksanvisning

ARNDT AIRWAY EXCHANGE CATHETER SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Arndt Airway Exchange Catheter Set consists of a radiopaque exchange catheter, an extra stiff wire guide with positioning marks, a removable Rapi-Fit® Adapter, and a bronchoscopic port.

INTENDED USE

The Arndt Airway Exchange Catheter Set is intended for the exchange of laryngeal mask airways (LMA) and endotracheal tubes using a fiberoptic bronchoscope, bronchoscopic port, wire guide and exchange catheter. The product may be used for emergency, urgent and elective airway management. The removable Rapi-Fit Adapter permits use of ventilatory device, if necessary, during exchange procedure.

CONTRAINDICATIONS

None known

WARNINGS

- Attention should be paid to insertion depth of catheter into patient's airway and correct tracheal position of replacement endotracheal tube. Markers on the exchange catheter refer to distance from tip of catheter. Catheter and endotracheal tube should not be advanced beyond the carina.
- Ensure proper sizing of the exchange catheter within a **double lumen endotracheal tube**. Failure to do so may cause small fragments to be shaved off during removal of the exchange catheter.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in airway anatomy and in airway management and fiberoptic bronchoscopy techniques. Standard techniques for use of bronchoscopes, endotracheal tubes, and laryngeal mask airways should be employed.
- This product should not be used in replacement of endotracheal tubes whose inner diameter is smaller than 5 mm.
- The use of either electronic or chemical capnography is highly recommended to confirm intubation.
- A pulse oximeter is recommended during the use of this device.
- This device is not recommended for definitive long-term airway management. It is used for the exchange of an LMA or endotracheal tube.
- Care should be taken to maintain wire guide position during placement of exchange catheter.
- Patients with airway injuries may have significant spinal cord injury. When possible, the cervical spine should be immobilized before beginning the procedure.
- When using an LMA, ensure the cuff is properly inflated and breath sounds are present. If possible, confirm the position of the LMA by using chemical or electronic capnography.
- It is recommended that a sterile lubricant be applied to the exchange catheter, bronchoscope and new endotracheal tube prior to their introduction into patient.
- Before airway exam, ensure that fiberoptic bronchoscope is properly focused.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Airway bleeding
- Pneumothorax
- Infection
- Catheter dislodgement
- Hematoma

INSTRUCTIONS FOR USE

Use of Rapi-Fit® Adapter

A ventilatory device may be used at any time during the endotracheal tube exchange procedure by utilizing the attached Rapi-Fit Adapter.

1. To attach the Rapi-Fit Adapter (with either a 15 mm connector or a Luer lock connector), position the adapter on the catheter, then push the white collar forward and lock into position.

NOTE: Use of a high-pressure oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume.

2. To remove the adapter, pull the white collar back to release, then remove from the catheter.

LMA or Endotracheal Tube Exchange

NOTE: Ventilation should be continued during the exchange process, through the bronchoscopic port. Discontinue ventilation only during the actual exchange over the exchange catheter.

1. After ensuring that the LMA or endotracheal tube is correctly positioned within the trachea, insert the bronchoscopic port between the LMA or endotracheal tube and the ventilation source.
2. Advance the fiberoptic bronchoscope into the LMA or endotracheal tube through the bronchoscopic port. If an LMA is in place, the bronchoscope should be advanced through the vocal cords and into the trachea. If an endotracheal tube is in place, the bronchoscope should be advanced directly into the trachea through the endotracheal tube.
3. Advance the soft, flexible end of the wire guide through the working port of the bronchoscope until it is seen to exit the bronchoscope into the patient's trachea. Advance the wire guide to the level of the carina.
4. Carefully remove the bronchoscope, leaving the wire guide in place in the trachea, exiting the bronchoscopic port.
5. Maintaining wire guide position, advance the exchange catheter over the wire guide through the bronchoscopic port and through the LMA or endotracheal tube into the trachea.
6. Properly position the exchange catheter within the endotracheal tube by aligning the appropriate centimeter mark on the exchange catheter with the corresponding centimeter mark on the endotracheal tube. This placement is determined by visualizing the indicated centimeter length of the endotracheal tube, in place, as shown on its surface scale. (For example, an endotracheal tube that has been shortened to 24 cm should have the 24 cm marker of the exchange catheter aligned at the 24 cm mark of the endotracheal tube.)
7. Fully deflate the cuff of the endotracheal tube or LMA and remove the Rapi-Fit Adapter (if on the exchange catheter) while maintaining position of the exchange catheter. Remove the endotracheal tube or LMA, leaving the exchange catheter in place.
8. Advance the new endotracheal tube over the exchange catheter into the patient's airway. Passage of the tube may be aided by advancing the operator's fingers into the patient's airway and directing the tongue forward.

NOTE: It is recommended that a sterile lubricant be applied to the tip of the endotracheal tube prior to advancing the tube.

9. Remove the exchange catheter and inflate the balloon cuff of the new endotracheal tube. Re-establish ventilatory device use and secure the new endotracheal tube in place. Confirm its position using standard methods (e.g., capnography, breath sounds and chest x-ray).

NOTE: In cases where the endotracheal tube does not easily advance, the exchange catheter with a Rapi-Fit Adapter attached may be used as a source of jet ventilation to provide oxygenation. If a high-pressure oxygen source is used for insufflation (e.g., jet ventilation), begin at lower pressure and increase gradually. Rising chest wall, pulse oximetry and oral air flow should be carefully monitored.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd