## Arndt Endobronchial Blocker Set

Instructions for Use

EN  
3

## Arndt endobronkialt blokeringsanordningssæt

Brugsanvisning

DA  
8

## Arndt Endobronchialblocker-Set

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DE  
13

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C - T - A E B S - R E V - 3

## ARNDT ENDOBRONCHIAL BLOCKER SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Arndt Endobronchial Blocker Set is designed for use with a standard endotracheal tube and a small diameter (pediatric) fiberoptic bronchoscope to provide one-lung ventilation. The endobronchial blocker uses a guide loop assembly that fits through the lumen of the endobronchial blocker and exits from the endobronchial blocker's distal end to form a small, adjustable guide loop. This assembly may be advanced or retracted at the endobronchial blocker's proximal end to enlarge or decrease the size of the guide loop. A small diameter (pediatric) fiberoptic bronchoscope is advanced through the guide loop to couple the bronchoscope and the endobronchial blocker together. This allows the endobronchial blocker to follow the bronchoscope. After the area to be blocked is entered, the guide loop exits the end of the bronchoscope. The lung may then be blocked by inflating the balloon within either the right or left mainstem bronchus. **The guide loop may then be removed and can be replaced if necessary (9.0 French size only).** A suction adapter is provided for an alternative method of lung deflation, if desired.

### INTENDED USE

The Arndt Endobronchial Blocker Set is intended to differentially intubate a patient's bronchus in order to isolate the left or right lung for procedures that require one-lung ventilation. The product is intended for use by physicians trained and experienced in the use of fiberoptic bronchoscopes and airway anatomy. Standard techniques for use of fiberoptic bronchoscopes and endobronchial blockers should be employed.

### CONTRAINDICATIONS

- Airway diameter insufficient to allow passage of the Arndt Endobronchial Blocker
- Fiberoptic equipment unavailable

### WARNINGS

- **If any resistance is encountered in removing the blocker, the blocker tip should be inspected fiberoptically to ensure integrity.**
- The enclosed blocker balloon is a high-volume, low-pressure design. Excessive manipulation over a prolonged period may cause balloon rupture or deflation.
- To prevent airway damage, the balloon should never be overinflated.
- If ventilation should be acutely difficult during endobronchial blockade, the balloon should be immediately deflated and the balloon position inspected.
- The balloon should be retracted prior to bronchial clamping during pneumonectomy procedures.
- It may not be possible to adequately ventilate a patient during placement of the endobronchial blocker due to the presence of a bronchoscope and the endobronchial blocker in the endotracheal tube; in addition, higher positive end expiratory pressures may result.
- When using the 5.0 French Arndt Pediatric Endobronchial Blocker with smaller endotracheal tubes (i.e., 4.5 mm ET), higher airway pressure may result from the presence of the blocker within the lumen of the endotracheal tube as shown below:

| Control                 | Volume<br>4.5 mm ET<br>Tube | Ventilation<br>4.5 mm ET Tube<br>with 5.0 FR Blocker |
|-------------------------|-----------------------------|------------------------------------------------------|
| Tidal Volume            | 379 ± 1 cc                  | 379 ± 1 cc                                           |
| Peak Airway<br>Pressure | 21 ± 0 cm H2O               | 25 ± 0 cm H2O                                        |
| Mean Airway<br>Pressure | 8 ± 0 cm H2O                | 8 ± 0 cm H2O                                         |

- When using the 7.0 French Arndt Endobronchial Blocker with smaller endotracheal tubes (i.e., 6.0 mm ET), higher airway pressure may result from the presence of the blocker within the lumen of the endotracheal tube as shown below:

| Control                 | Volume<br>6.0 mm ET<br>Tube | Ventilation<br>6.0 mm ET Tube<br>with 7.0 FR Blocker |
|-------------------------|-----------------------------|------------------------------------------------------|
| Tidal Volume            | 724 ± 2 cc                  | 722 ± 2 cc                                           |
| Peak Airway<br>Pressure | 25 ± 0 cm H2O               | 28 ± 0 cm H2O                                        |
| Mean Airway<br>Pressure | 10 ± 0 cm H2O               | 11 ± 0 cm H2O                                        |


## PRECAUTIONS

- The use of pulse oximetry is recommended when using this device.
- Caution is recommended when working near the hilum. The balloon position should be verified to prevent inadvertent balloon damage.
- During the use of this device for one-lung ventilation, an increased FIO<sub>2</sub> for the ventilated lung is often needed to help maintain adequate arterial oxygen saturation. However, in some circumstances, an increased FIO<sub>2</sub> for the ventilated lung may not be sufficient to maintain adequate arterial oxygen saturation, and other techniques may also be necessary.
- Inflate balloon initially under direct vision to ensure correct position and placement.  
**NOTE:** Placement in the right mainstem bronchus may result in herniation of the balloon into the right upper lobe bronchus and thereby occlude it.
- Care should be taken to ensure the balloon remains fully inflated during longer procedures.
- During one-lung ventilation, the patient should be paralyzed to help prevent dislodgement of the balloon.
- The endobronchial blocker should be placed with the patient in a supine position.
- To avoid damage to the small diameter fiberoptic bronchoscope, the guide loop size should be adjusted with care.
- The position of the balloon should be reconfirmed after positioning the patient for the procedure.
- Assure the endotracheal tube is securely taped in place during use of the endobronchial blocker.
- The capnograph waveform may become distorted during use of the endobronchial blocker.
- The Arndt Endobronchial Blocker Set can be placed through a standard endotracheal tube. For optimal performance, the largest endotracheal tube appropriate for the patient's anatomy is recommended.
- The smallest recommended size endotracheal tube for use with the Arndt Endobronchial Blocker size 9.0 French is 7.5 mm; 7.0 French is 6.0 mm; and 5.0 French is 4.5 mm.
- The 5.0 French Arndt Pediatric Endobronchial Blocker, when placed in an endotracheal tube with a diameter smaller than 6.0 mm, will cause slightly higher Peak Airway Pressures when attempting to achieve constant tidal volume.
- The 7.0 French Arndt Endobronchial Blocker, when placed in an endotracheal tube with a diameter smaller than 7.0 mm, will cause slightly higher Peak Airway Pressures when attempting to achieve constant tidal volume.
- Following insertion of the blocker balloon through the multiport adapter, the balloon should be test inflated.
- If using supplied suction adapter during lung deflation procedure, ensure that you do not inadvertently attach it to the pilot balloon assembly. Suction adapter should only be attached to Luer lock connector.

## INSTRUCTIONS FOR USE


1. An Arndt Multiport Airway Adapter should be placed on the endotracheal tube and the patient ventilated with 100% oxygen. **NOTE:** The guide loop should be adjusted to loosely approximate the diameter of the bronchoscope.
2. Fully deflate the balloon on the endobronchial blocker.
3. Generously lubricate the bronchoscope, endobronchial blocker and inside of the endotracheal tube.

**NOTE:** Use of a silicone-based lubricant is suggested for airway instrumentation with bronchoscopes. Refer to the lubricant manufacturer's labeling for specific information concerning recommendations for use.


**Fig. 1**

4. Place the endobronchial blocker through the blocker port of the Arndt Multiport Airway Adapter, advancing the endobronchial blocker until the guide loop is visualized within the main body of the multiport airway adapter. **(Fig. 1)** **NOTE:** Air leakage can be controlled by screwing down the cap on the blocker port.


**Fig. 2**

5. Place the bronchoscope through the diaphragm of the bronchoscopy port of the Arndt Multiport Airway Adapter, advancing the bronchoscope until the guide loop comes into view. **(Fig. 2)**
6. Advance the bronchoscope into and through the guide loop, coupling the endobronchial blocker to the bronchoscope. This may be done within the multiport adapter, the endotracheal tube or the trachea proper. **(Fig. 2)** **NOTE:** The endobronchial blocker and bronchoscope may also be coupled together prior to attachment of the Arndt Multiport Airway Adapter to the endotracheal tube.


**Fig. 3**

7. Advance the bronchoscope into the section of the lung to be blocked. **(Fig. 3)** **NOTE:** The blocker guide loop may be cinched onto the bronchoscope by unscrewing the guide loop port and pulling back on the snare. The snare may be affixed in place by inserting syringe or CPAP adapter into the guide loop port.


**Fig. 4**

8. Keeping the bronchoscope position stable (loosen guide loop if snare is cinched onto bronchoscope), advance the endobronchial blocker until the guide loop is seen to exit the end of the bronchoscope. **(Fig. 4)**

**Fig. 5**

9. Retract the bronchoscope. The endobronchial blocker may be either advanced or retracted to place the balloon within either the right or left mainstem bronchus. **(Fig. 5)**

**Fig. 6**

10. Place patient in final position. Under bronchoscopic vision, inflate the balloon with air using the pilot balloon assembly. The balloon should fill the entire endobronchial lumen to be blocked and not herniate into the mainstem trachea. Do not overinflate balloon. **(Fig. 6)**

## AVERAGE INFLATION VOLUMES

| French Size | Balloon             | Inflation Volume |
|-------------|---------------------|------------------|
| 5.0         | Pediatric Spherical | .5 cc - 2.0 cc   |
| 7.0         | Adult Spherical     | 2.0 cc - 6.0 cc  |
| 9.0         | Adult Spherical     | 4.0 cc - 8.0 cc  |
| 9.0         | Adult Elliptical    | 6.0 cc - 12.0 cc |

**WARNING:** The lungs should be carefully auscultated following initial endobronchial blocker placement and balloon inflation to ensure proper functioning of the endobronchial blocker. If ventilation should be acutely difficult during endobronchial blockade, the balloon should be immediately deflated.

11. Following correct balloon placement, the balloon may be deflated until one-lung ventilation is required.
12. To achieve one-lung ventilation, the balloon is inflated with the same volume of air that was used during initial placement.

**WARNING:** To prevent airway damage, the balloon should never be overinflated.

13. Once the blocker is in final position, tighten the blocker port on the multiport adapter. This contains a Tuohy-Borst valve that will prevent the blocker from moving and will prevent a gas leak.
14. After the balloon is properly inflated, remove the guide loop assembly from the endobronchial blocker and inspect the position of the endobronchial blocker using the bronchoscope. Reconfirm the position of the balloon if the patient is repositioned. **NOTE:** When removing the guide loop assembly from the endobronchial blocker, excessive force should not be used as this may dislodge the endobronchial blocker.
15. **The guide loop assembly may be re-inserted if repositioning of the blocker is necessary (9.0 French size only).**
16. Upon completion of one-lung ventilation, deflate the balloon and remove the catheter from bronchus. **NOTE:** Assure complete deflation of the balloon before attempted removal of the endobronchial blocker. The Arndt Multiport Airway Adapter should be detached upon removal of the endobronchial blocker.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.