**COOK**

**MEDICAL**

**CE**  
0088

**EN**  
**3**

## **Cohen Tip Deflecting Endobronchial Blocker**

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Bruksanvisning

C \_ T \_ A E B T \_ R E V 5

## COHEN TIP DEFLECTING ENDOBRONCHIAL BLOCKER

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Cohen Tip Deflecting Endobronchial Blocker consists of a blocker catheter and a multi-port airway adapter. A suction adapter is provided for an alternative method of lung deflation, if desired. The endobronchial blocker uses a tip-deflecting mechanism to allow the tip to be directed into the right or left mainstem bronchus. A small-diameter (pediatric) fiberoptic bronchoscope is used to visualize and aid in the placement of the endobronchial blocker. The lung may then be blocked by inflating the balloon within either the right or left mainstem bronchus.

The endobronchial blocker does not obstruct the working channel of the bronchoscope during placement. Refer to label for further information.

### INTENDED USE

The Cohen Tip Deflecting Endobronchial Blocker is intended for use to differentially intubate a patient's bronchus in order to isolate the left or right lung for procedures that require one-lung ventilation. The product is intended for use by physicians trained and experienced in one-lung ventilation. Standard techniques for placement of endobronchial blockers should be employed.

### CONTRAINDICATIONS

- Airway or endotracheal tube diameter insufficient to allow passage of the Cohen Tip Deflecting Endobronchial Blocker.
- Fiberoptic equipment is unavailable. It is not possible to properly position the blocker without a fiberoptic bronchoscope.

### WARNINGS

- Always advance the endobronchial blocker under direct fiberoptic vision.
- Do not advance the endobronchial blocker if resistance is encountered.
- If ventilation should become acutely difficult during endobronchial blockade, the balloon should be immediately deflated and the balloon position inspected.
- Caution is recommended when working near the hilum. The balloon position should be verified to prevent inadvertent balloon damage or movement. This is particularly recommended during right-sided procedures in which the balloon is positioned near the carina.
- As per all cases of one-lung ventilation, the use of pulse oximetry is mandatory when using this device.
- It may not be possible to adequately ventilate a patient during placement of the endobronchial blocker due to the presence of a bronchoscope and the endobronchial blocker in the endotracheal tube.
- An increased FiO<sub>2</sub> for the ventilated lung is often needed to help maintain adequate arterial oxygen saturation.
- In some circumstances, an increased FiO<sub>2</sub> for the ventilated lung may not be sufficient to maintain adequate arterial oxygen saturation, and other standard techniques routinely employed during one-lung ventilation may be necessary.
- Placement in the right mainstem bronchus may result in herniation of the balloon into the right upper lobe bronchus and thereby occlude it.
- During one-lung ventilation, the patient should be paralyzed to help prevent dislodgement of the balloon.
- If any resistance is met in removing the blocker after surgery, the blocker tip should be inspected fiberoptically to ensure integrity.

### PRECAUTIONS

- The endobronchial blocker balloon is a high-volume, low-pressure design. Excessive manipulation over a prolonged period may cause balloon rupture or deflation.
- This device should only be used by those familiar with the use of small diameter (pediatric) fiberoptic bronchoscopes and airway anatomy.
- Retract balloon prior to bronchial clamping during pneumonectomy procedures.
- Inflate balloon initially under direct vision to ensure correct position and placement.
- Ensure the balloon remains fully inflated during longer (>1 hour) procedures.
- For best results, initially, the blocker should be placed with the patient in a supine position, although the blocker can be easily and safely directed to the desired bronchus when the patient is in the lateral decubitus.
- To avoid damage to the small-diameter fiberoptic bronchoscope, standard lubricant may be used for the bronchoscope and the blocker.
- Reconfirm balloon position after positioning the patient for the procedure.

- Securely tape in place the endotracheal tube during use of the endobronchial blocker.
- During use of the endobronchial blocker, the capnograph waveform may become distorted.
- The blocker can be placed through a standard endotracheal tube. For optimal performance, the largest endotracheal tube is recommended.
- The smallest recommended size endotracheal tube for use with the 9.0 French blocker is 7.0 mm.
- Following insertion through the multi-port adapter, test inflate the balloon.
- If using supplied suction adapter during lung deflation procedure, ensure that you do not inadvertently attach it to the pilot balloon assembly. Suction adapter should only be attached to Luer lock connector.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

## POTENTIAL ADVERSE EVENTS

The adverse events are those of lung isolation and one-lung ventilation:

- Hypoxia
- Hypoxemia
- Endobronchial irritation
- Tracheal injury

## INSTRUCTIONS FOR USE

Only those physicians familiar with the use of this device, the use of small diameter (pediatric) fiberoptic bronchoscopes and airway anatomy should use the Cohen Tip Deflecting Endobronchial Blocker.

**NOTE:** The use of pulse oximetry is mandatory when using this device.

**NOTE:** Ventilation parameters are affected by inserting a 9.0 French endobronchial blocker or a blocker and a scope into a 7.0 mm ETT. However, the increases in airway pressures and decreases in compliance that result from a 9.0 French endobronchial blocker only are not clinically significant.

**NOTE:** At physician's discretion, CPAP may be administered using standard techniques through attachment of the 15 mm CPAP adapter provided with the endobronchial blocker. The 15 mm CPAP adapter is compatible with commercially available CPAP circuits. Exercise care to monitor administration of oxygen and to not over-inflate the isolated lung.

1. A multi-port airway adapter should be placed on the endotracheal tube and the patient ventilated with 100% oxygen.
2. Fully deflate balloon on the endobronchial blocker.
3. Generously lubricate the bronchoscope, the endobronchial blocker and the inside of the endotracheal tube.

**NOTE:** Use of a silicone-based lubricant is suggested for airway instrumentation with bronchoscopes. Refer to the lubricant manufacturer's labeling for specific information concerning recommendations for use.

4. The bronchoscope should first be placed through the diaphragm of the multi-port airway adapter and advanced until the carina is well visualized.
5. Place Cohen Tip Deflecting Endobronchial Blocker through the blocker port of the multi-port airway adapter.
6. To facilitate guidance into the desired main bronchus, the blocker should be inserted with the rotating wheel facing the side that is to be blocked.
7. Advance the blocker until it is out of the endotracheal tube and is visualized by the fiberoptic bronchoscope.

**NOTE:** When using the blocker with a 7.0 mm or 7.5 mm endotracheal tube, the bronchoscope must be removed prior to introduction of the blocker. The blocker should then be placed through the blocker port of the multi-port airway adapter and advanced until the end of the blocker is 3 cm beyond the tip of the endotracheal tube (depth marks on the blocker and the endotracheal tube should be used to determine correct placement). Then the bronchoscope should be reintroduced to visualize the blocker.

8. The black mark at the proximal end of the blocker will indicate the direction of the tip deflection.
9. If necessary, the blocker may be rotated by twisting the proximal end of the blocker until the black mark faces the main bronchus to be blocked.
10. Keeping the fiberoptic bronchoscope stable, deflection of the blocker tip should be done with the wheel under direct vision.
11. In most cases, a 10–45 degree deflection will be sufficient to direct the blocker into the main bronchus to be blocked.
12. Advance the blocker into the main bronchus to be blocked.
13. Place patient in final position. Under bronchoscopic vision, inflate balloon with air using the pilot balloon assembly. The balloon should fill the entire endobronchial lumen to be blocked. Do not over-inflate balloon.

**NOTE:** In most patients, the balloon will require between 6–9 ml of air to seal the lumen of the bronchus.

**WARNING: The lungs should be carefully auscultated following initial endobronchial blocker placement and balloon inflation to ensure proper functioning of the endobronchial blocker. If ventilation should become acutely difficult during endobronchial blockade, the balloon should be immediately deflated.**

14. Following correct balloon placement, the balloon may be deflated until one-lung ventilation is required.

15. To achieve one-lung ventilation, the balloon is inflated with the same volume of air that was used during initial placement.

**WARNING: To prevent airway damage, the balloon should never be over-inflated.**

16. After balloon is properly inflated, inspect the position of the endobronchial blocker using the bronchoscope. Reconfirm the position of the balloon if the patient is repositioned.

17. Upon completion of one-lung ventilation, deflate the balloon and remove the catheter from bronchus.

**NOTE:** Assure complete deflation of balloon before attempted removal of the endobronchial blocker.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.