COOK

MEDICAL

**EN
5**

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









9

10

AINTREE INTUBATION CATHETER WITH RAPI-FIT® ADAPTERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Aintree Intubation Catheter (AIC) is a 19 French, blunt-tipped, radiopaque catheter with centimeter markings. The catheter has a 4.7 mm inner diameter (ID), appropriate for a bronchoscope with a maximum outer diameter of 4.2 mm. **Does not contain natural rubber latex.**

The device consists of the following components:

- Radiopaque catheter
- Rapi-Fit adapter - 15 mm connector
- Rapi-Fit adapter - Luer lock connector
- Double swivel connector

INTENDED USE

The Aintree Intubation Catheter is intended for exchange of a supraglottic airway device (SAD) to an endotracheal tube (ETT) under bronchoscopic assistance, and for ETT exchange in adult and pediatric patients.

The 19 French catheter is recommended for use with a supraglottic airway device and placement of a single-lumen endotracheal tube with an inner diameter of 7 mm or larger.

When used for high-pressure oxygenation with a Luer lock connector, the 19 French catheter is recommended for patients older than 12 years of age.

CONTRAINDICATIONS

None known

WARNINGS

- Advancing the catheter beyond the carina may introduce additional procedural risks and trauma.
- Attention should be paid to the insertion depth of the AIC into the patient's airway and correct tracheal position of replacement ETT. Markers on the AIC refer to distance from the distal tip of catheter.
- Take care to avoid injuring the epiglottis, glottis, trachea, bronchi, and lung parenchyma, and to avoid perforating the sinus pyriformis, while using this device.
- Ensure proper sizing of the AIC within the ETT or SAD.
- To avoid barotrauma, ensure that the tip of the AIC is always above the carina, preferably 2-3 cm.
- Use of the Rapi-Fit adapter for oxygenation may be associated with a risk of barotrauma.
- Use of an oxygen source should be considered only if the patient has sufficient egression of the insufflated gas volume.
- If a high-pressure oxygen source is used for insufflation (e.g., jet ventilator), begin at a lower pressure (i.e., 5 psi) and work up gradually. Rising chest wall, pulse oximetry, and oral air flow should be carefully monitored.
- Ensure that the Rapi-Fit Adapter is securely connected to the AIC prior to oxygen delivery. Failure to properly secure the adapter to the AIC may result in hypoxia and serious adverse events.
- High-pressure oxygenation with Luer lock connector should only be used in patients older than 12 years of age. If used in patients 12 years old or younger, the maximum airway pressure may be higher than 28 cm H₂O.

PRECAUTIONS

- This product is intended for use by clinicians trained and experienced in airway management techniques. Standard techniques for placement and exchange of ETTs or SADs should be employed.
- This product is not intended for intravascular use.
- The AIC should slide easily onto the bronchoscope. If it does not slide easily, there is a risk that the catheter may become lodged on the bronchoscope and cause damage.
- This product should not be used in replacement of ETTs or SADs whose inner diameter is smaller than 7 mm.
- The use of capnography is highly recommended to confirm intubation.
- A pulse oximeter is recommended during the use of this device.
- It is recommended that a sterile lubricant or sterile water be applied to the AIC, bronchoscope, and new ETT or SAD prior to their introduction into the patient.
- Do not attempt to place a bronchoscope through the Rapi-Fit adapter. The adapter must be removed prior to scope insertion.
- When using an SAD, ensure the cuff is properly inflated (if applicable) and breath sounds are present.
- This device is not recommended for definitive, long-term airway management.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Barotrauma
- Pneumothorax
- Hypoxia
- Perforation of the trachea, bronchi, or lung parenchyma
- Injury to the epiglottis
- Catheter dislodgment or migration
- Failed endotracheal tube placement

INSTRUCTIONS FOR USE

Bronchoscopic Assisted Exchange/Intubation

1. After a supraglottic airway device (SAD) has been placed, if the patient requires a single-lumen ETT, proceed to add the double swivel connector into the anesthetic circuit, if not already present. (**Fig. 1**) This allows the patient to be adequately ventilated during the bronchoscopic intubation procedure. **NOTE:** Ventilation should be continued during the exchange process, through the double swivel connector. Discontinue ventilation only during the actual exchange over the AIC.
2. Lubricate the bronchoscope using sterile medical-grade lubricant approved for use in the airway, according to hospital protocol.
3. Place the appropriately sized bronchoscope through the AIC and tape the catheter securely to the scope. **NOTE:** The catheter should slide easily onto the scope. If it does not slide easily, there is a risk that the catheter may become lodged on the bronchoscope and cause damage. The distal 3–4 cm of the bronchoscope should extend beyond the distal tip of the AIC. This allows adequate tip deflection of the bronchoscope.
4. While maintaining adequate ventilation, pass the bronchoscope and AIC assembly through the bronchoscopic port of the double swivel connector. (**Fig. 2**)
5. Direct the bronchoscope and AIC assembly through the vocal cords and down toward the carina.
6. While maintaining constant position of the AIC, remove the bronchoscope. Leave the AIC and the SAD *in situ*.
7. Deflate the cuff of the SAD (if applicable). While maintaining constant position of the AIC, remove the SAD, leaving just the AIC in the trachea. (**Fig. 3**)
8. Using sterile lubricant, lubricate the proximal end of the AIC, and load the ETT onto the AIC's proximal end.
9. While maintaining the position of the AIC, use the outer margin of the patient's mouth as a landmark and advance the ETT over the AIC into appropriate position within the trachea. (**Fig. 4**) **NOTE:** It may be necessary to rotate the ETT 90 degrees counter-clockwise while advancing it. This will position the bevel of the tube toward the anterior aspect of the vocal cords.
10. Hold the endotracheal tube in position and remove the AIC.
11. Inflate the cuff of the endotracheal tube. Reestablish ventilation and secure the new ETT in place.
12. Reinsert the bronchoscope to confirm the correct position of the ETT.
13. Remove the bronchoscope.

Endotracheal Tube Exchange

1. Before advancing the Aintree Intubation Catheter (AIC) into the endotracheal tube (ETT) to be replaced, confirm correct ETT position.
2. Using the outer margin of the patient's mouth or nasal orifice as a landmark, note the marking on the ETT. A piece of tape or other marker may be placed on the AIC at the corresponding distance from the tip to aid in correct placement within the ETT.
3. Advance the AIC, sideported end first, into the ETT to be replaced. (**Fig. 5**) **NOTE:** It is recommended that a sterile lubricant or sterile water be applied to the orifice of the ETT prior to introduction of the AIC.
4. Properly position the AIC within the ETT by aligning the appropriate centimeter mark on the AIC with the corresponding centimeter mark on the ETT. This placement is determined by visualizing the indicated centimeter length of the ETT, in place, as shown on its surface scale. (For example, an ETT that has been shortened to 24 cm should have the 24 cm marker of the AIC aligned at the 24 cm mark of the ETT.)
5. Fully deflate the cuff of the ETT. While maintaining the position of the AIC, remove the ETT, leaving the AIC in place. (**Fig. 6**)
6. While maintaining the position of the AIC, use the patient's mouth or nasal orifice (depending on approach) as a landmark and advance the new ETT over the AIC into appropriate position. (**Fig. 7**) **NOTE:** It is recommended that a sterile lubricant be applied to the tip of the ETT prior to advancing the ETT.
7. Remove the AIC and inflate the cuff of the new ETT. Reestablish ventilation and secure the new ETT in place. Confirm the new ETT's position using standard methods (e.g., capnography, breath sounds, and chest x-ray). (**Fig. 8**)

Use of the Rapi-Fit® Adapter

Rapi-Fit adapters should only be used when oxygen requirements are high and intubation is unsuccessful. Use of an oxygen source should only

be considered if the patient has sufficient egression of the insufflated gas volume. If an oxygen source is used for insufflation, begin at a lower pressure and work up gradually. Observe the chest for outward and inward movements to confirm oxygen insufflation and egression. Pulse oximetry and oral air flow should be carefully monitored as well. In cases of upper airway obstruction, gas discharge from the patient's lungs may require more time.

1. To attach the Rapi-Fit Adapter, position the adapter on the catheter, then push the white collar forward and lock into position. (**Fig. 9**)
2. To remove the adapter, pull the white collar back to release, and then remove from the catheter. (**Fig. 10**)

Catheter Oxygenation

The Aintree Intubation Catheter is designed for both positive airway pressure ventilation (Rapi-Fit Adapter - 15 mm connector) and jet ventilation (Rapi-Fit Adapter - Luer lock connector) for both adult and pediatric patients. In the table below, the delivered minute volume and the measured airway pressure are given for jet ventilation in adult and pediatric patients with healthy lung tissue.

Real Part Number (RPN)	Patient Subgroup and Age Range	Delivered Minute Volume ¹ (L/min)	Measured Airway Pressure ¹ (cm H ₂ O)	
			Mean	Maximum
C-CAE-19.0-56-AIC-01	Adolescent > 12 years through 21 years	7.5	20.2	21.4
	Adult > 21 years	11.4	14.6	15.1

¹ The test conditions were used in an active model. See Testing Conditions table for additional details.

The presented numbers are mean values. The following testing conditions were used in an active model mode, tested with ASL 5000, Ingmar Medical, Ltd., for each patient subgroup.

Testing Conditions

Input pressure was set to 50 psi for both patient type models.

Patient type	Model Body Weight (kg)	Inspiratory Time (s)	Expiratory Time (s)	Breaths per minute	Resistance (cm H ₂ O/L/s)	Lung Compliance (mL/cm H ₂ O)
Adolescent (> 12 years through 21 years)	50	0.9	2.1	20	5	40
Adult (> 21 years)	80	1.0	4.0	12	3	100

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from clinicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2018
Bloomington, IN 47404 U.S.A.