COOK

MEDICAL

CE
0088

EN
3

Aintree Intubation Catheter

Instructions for Use

DA
4

Aintree intubationskateter

Brugsanvisning

DE
6

Aintree Intubationskatheter

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EL
8

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C _ T _ A I C _ R E V 2

AINTREE INTUBATION CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Aintree Intubation Catheter (AIC) is a blunt-tipped, radiopaque catheter with centimeter marks to facilitate placement. The catheter is supplied with Rapi-Fit™ adapters (with 15 mm connector and Luer lock connector), permitting use of ventilatory devices during exchange procedure. The lumen design with distal sideports enhances air flow. This set also includes a bronchoscopic port.

The Aintree Intubation Catheter is a 19 French catheter for exchange of endotracheal tubes with an inner diameter (ID) of 7 mm or larger. The catheter has a 4.7 mm ID, allowing it to be pre-loaded onto a fiberoptic scope with a maximum outer diameter (OD) of 4.2 mm.

INTENDED USE

The Aintree Intubation Catheter has been designed for assisted fiberoptic intubation and for uncomplicated, atraumatic endotracheal tube exchange.

CONTRAINDICATIONS

None known

WARNINGS

- Possible allergic reactions should be considered.
- Do not advance the catheter beyond the carina. Attention should be paid to the depth of insertion of the catheter into the patient's airway, and to the correct tracheal position of the endotracheal tube. Markers on the catheter refer to distance from the tip of the catheter.
- Take care not to injure the epiglottis and glottis, or to perforate the sinus pyramidalis, trachea or bronchus.
- Use of the Rapi-Fit adapter for oxygenation may be associated with a risk of barotrauma.
- Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume.
- If a high-pressure oxygen source is used for insufflation (e.g., jet ventilator), begin at a lower pressure and work up gradually. Rising chest wall, pulse oximetry and oral air flow should be carefully monitored.
- Ensure proper sizing of the catheter within an endotracheal tube. Failure to do so may cause small fragments to be shaved off during removal of the catheter.
- To avoid barotrauma, ensure that the tip of the AIC catheter is always above the carina, preferably 2-3 cm.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in airway management techniques. Standard techniques for placement and exchange of endotracheal tubes should be employed.
- This product is not intended for intravascular use.
- Lubricate the catheter and endotracheal tube before use.
- Do not attempt to place a fiberoptic scope through the Rapi-Fit adapter. The adapter must be removed prior to scope insertion.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Barotrauma
- Perforation of the bronchi or lung parenchyma
- Pneumothorax

INSTRUCTIONS FOR USE

Use of Rapi-Fit™ Adapter

Rapi-Fit adapters should only be used when oxygen requirements are high and intubation is unsuccessful. Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume. If an oxygen source is used for insufflation, begin at a lower pressure and work up gradually. Observe the chest for outward and inward movements to confirm oxygen insufflation and egression. Pulse oximetry and oral air flow should be carefully monitored as well. In cases of upper airway obstruction, gas discharge from the patient's lungs may require more time.

1. To attach the Rapi-Fit Adapter, position the adapter on the catheter, then push the white collar forward and lock into position.
2. To remove the adapter, pull the white collar back to release, then remove from the catheter.

Supraglottic Airway Device/Endotracheal Tube Exchange

1. Insert a laryngeal mask airway (LMA) of size 2 or larger, or a cuffed oropharyngeal airway (COPA), into the patient's airway. This allows the patient to be adequately ventilated during the fiberoptic intubation

procedure. Add a bronchoscopic swivel adapter into the anesthetic circuit, if not already present.

2. Lubricate the fiberoptic bronchoscope using either a sterile lubricant or sterile water.
3. Mount the Aintree Intubation Catheter onto the bronchoscope and tape securely. The catheter should slide easily onto the scope. If it does not slide easily, there is a risk that the catheter may become lodged on the bronchoscope and cause damage. The distal 3-4 cm of the bronchoscope should extend beyond the distal tip of the catheter. This allows adequate tip deflection of the bronchoscope.
4. While maintaining adequate ventilation, pass the bronchoscope and AIC assembly through the bronchoscopic port on the swivel adapter and advance through the LMA/COPA through the vocal cords.
5. Direct the bronchoscope/AIC assembly through the vocal cords and down to the carina.
6. While maintaining constant position of the AIC, remove the bronchoscope, leaving the AIC and the LMA/COPA *in situ*.
7. Deflate the cuff of the LMA/COPA and, while maintaining constant position of the AIC, remove the LMA/COPA, leaving just the AIC in the trachea.
8. Lubricate the AIC, and load the endotracheal tube onto its proximal end.
9. While maintaining position of the AIC, using patient's upper incisors as a landmark, advance the endotracheal tube over the AIC into the trachea, and position appropriately. **NOTE:** It may be necessary to rotate endotracheal tube 90 degrees counter-clockwise while advancing, which will position the bevel of the tube toward the anterior aspect of the vocal cords.
10. Hold the endotracheal tube in position and remove the AIC.
11. Inflate the cuff of the tube. Re-establish ventilation and secure the new tube in place.
12. Reinsert the bronchoscope to confirm the correct position of the endotracheal tube.
13. Remove the bronchoscope.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2013
Bloomington, IN 47404 U.S.A.