EN
4

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
Fig. 6


Fig. 7

ARNDT EMERGENCY CRICOTHYROTOMY CATHETER SET

Caution: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Arndt Emergency Cricothyrotomy Catheter is a 9.0 French, 6 cm flexible sheath with a built-in 15 mm Luer lock connector. The set includes various components used for Seldinger placement of a cricothyrotomy catheter.

INTENDED USE

The Arndt Emergency Cricothyrotomy Catheter Set is used for emergency airway access when endotracheal intubation and/or ventilation cannot be performed. Airway access is achieved utilizing percutaneous entry (Seldinger) technique via the cricothyroid membrane. Subsequent dilation of the tract and tracheal entrance site permits passage of the emergency airway catheter.

CONTRAINDICATIONS

None known.

WARNINGS

- The Arndt Emergency Cricothyrotomy Catheter Set is not specifically designed for pediatric applications. Use with pediatric patients should be determined by attending physician. Use in children below the age of 12 is not recommended.
- Consideration should be given to the following medical and anatomic conditions:
 - Distorted anatomy
 - Subcutaneous abscess
 - Hematoma
 - Post-operative scarring/radiation
 - Coagulopathies or systemic thrombolytic therapy

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in proper emergency airway techniques. Standard techniques for emergency performance of a Seldinger cricothyrotomy should be employed.
- The use of either chemical or electronic capnography is recommended to confirm correct airway placement.
- Pulse oximetry is recommended during placement and use of this device.
- This device is not recommended for definitive, long-term airway management. Intubation, tracheostomy or the return of spontaneous ventilation is recommended before removing the airway catheter.
- Patients in need of cricothyrotomy may have significant spinal injury. In patients who have sustained significant trauma, the cervical spine should be immobilized throughout the procedure if possible.
- Whenever possible and appropriate, utilize aseptic technique and local anesthetic for the procedure.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Bleeding
- Pneumothorax
- Pneumomediastinum
- Subcutaneous emphysema
- Tracheoesophageal fistula
- Cutaneous fistula
- Infection
- Hematoma
- Catheter dislodgement

INSTRUCTIONS FOR USE

1. Identify the cricothyroid membrane between the cricoid and thyroid cartilages. (**Fig. 1**)
2. Carefully palpate the cricothyroid membrane and, while stabilizing the cartilage, make a vertical incision in the midline using the short handled scalpel blade. A sufficient incision to allow introduction of the dilator and catheter is recommended. (**Fig. 2**)

3. Attach the supplied syringe to the introducer needle, and advance it through the incision into the airway at a 45 degree angle to the frontal plane in the midline in a caudad direction. **When advancing the needle forward, verification of correct entrance into the airway can be confirmed by aspiration on the syringe, resulting in free air return. (Fig. 3)**
4. Remove the syringe from introducer needle, leaving the introducer needle in place. Advance the floppy end of the wire guide through the needle and into the airway several centimeters. **(Fig. 4)**
5. Remove the introducer needle, leaving the wire guide in place. **(Fig. 5)**
6. Advance the dilator, tapered end first, into the connector end of the airway catheter and secure at the hub by the female Luer lock connector. **(Fig. 6)**
7. Advance the emergency airway access assembly over the wire guide until the proximal stiff end of the wire guide is completely through and visible at the handle end of the dilator. **It is important to always visualize the proximal end of the wire during the airway insertion procedure to prevent its inadvertent loss into the trachea.**
8. Maintaining the wire guide position, advance the emergency airway access assembly over the wire guide with a reciprocating motion, and completely into the trachea. **Care should be taken not to advance the tip of the dilator beyond the tip of the wire within the trachea.**
9. Remove the wire guide and dilator simultaneously. **(Fig. 7)**
10. Fix the emergency airway catheter in place with the cloth tracheostomy tape strip, in a standard fashion.
11. Connect the airway catheter, using its standard 15-22 mm adapter, to an appropriate ventilatory device.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2011
Bloomington, IN 47404 U.S.A.