
COOK

MEDICAL


- EN**
4 Cook Airway Exchange Catheters with Rapi-Fit® Adapters
Instructions for Use
- CS**
7 Výměnné katetry Cook pro dýchací cesty s adaptéry Rapi-Fit®
Návod k použití
- DA**
10 Cook luftvejsudskiftningskateter med Rapi-Fit® adaptere
Brugsanvisning
- DE**
13 Cook Umintubationskatheter mit Rapi-Fit® Adaptern
Gebrauchsanweisung
- EL**
16 Καθετήρες αλλαγής αεραγωγών της Cook με προσαρμογείς Rapi-Fit®
Οδηγίες χρήσης
- ES**
19 Catéteres de intercambio para vía aérea Cook con adaptadores Rapi-Fit®
Instrucciones de uso
- FR**
23 Cathéters d'échange de sonde endotrachéale Cook avec adaptateurs Rapi-Fit®
Mode d'emploi
- IT**
26 Cateteri Cook per la sostituzione del tubo endotracheale con adattatori Rapi-Fit®
Istruzioni per l'uso
- NL**
29 Cook wisselkatheters voor tubes in de luchtwegen met Rapi-Fit® adapters
Gebruiksaanwijzing
- NO**
32 Cook luftveisutskiftingskatetre med Rapi-Fit®-adaptere
Bruksanvisning
- PL**
35 Cewniki do wymiany AEC Cook ze złączkami Rapi-Fit®
Instrukcja użycia
- PT**
39 Cateteres de substituição para vias aéreas Cook com adaptadores Rapi-Fit®
Instruções de utilização
- SV**
42 Cook-utbyteskatetrar med Rapi-Fit®-adaptrar
Bruksanvisning


C - T - C A E 2 - R E V O


1

2

3

4

COOK AIRWAY EXCHANGE CATHETERS WITH RAPI-FIT® ADAPTERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician or a properly licensed practitioner.

DEVICE DESCRIPTION

The Cook Airway Exchange Catheter (CAE) is a blunt-tipped radiopaque catheter with centimeter markings to help facilitate the removal and exchange of an endotracheal tube (ETT). **Does not contain natural rubber latex.**

The device consists of the following components:

- Radiopaque catheter
- Rapi-Fit adapter – 15 mm connector
- Rapi-Fit adapter – Luer lock connector

INTENDED USE

The Cook Airway Exchange Catheter is intended for endotracheal tube exchange in adult and pediatric patients.

The 8 French catheter is recommended for placement of an endotracheal tube with an inner diameter of 3 mm or larger.

When used for high-pressure oxygenation with a Luer Lock connector, the 8 French catheter is recommended for patients older than 1 month of age.

The 11 French catheter is recommended for placement of an endotracheal tube with an inner diameter of 4 mm or larger.

When used for high-pressure oxygenation with a Luer Lock connector, the 11 French catheter is recommended for patients older than 2 years of age.

The 14 French catheter is recommended for placement of an endotracheal tube with an inner diameter of 5 mm or larger.

When used for high-pressure oxygenation with a Luer Lock connector, the 14 French catheter is recommended for patients older than 2 years of age.

The 19 French catheter is recommended for placement of an endotracheal tube with an inner diameter of 7 mm or larger.

When used for high-pressure oxygenation with a Luer Lock connector, the 19 French catheter is recommended for patients older than 12 years of age.

CONTRAINDICATIONS

None known

WARNINGS

- Possible allergic reactions (e.g., to butyl rubber) should be considered.
- Advancing the catheter beyond the carina may introduce additional procedural risks and trauma.
- Attention should be paid to the insertion depth of a Cook Airway Exchange Catheter (CAE) into the patient's airway and correct tracheal position of replacement endotracheal tube (ETT). Markers on the CAE refer to distance from the distal tip of catheter.
- Take care to avoid injuring the epiglottis, glottis, trachea, bronchi, or lung parenchyma, and to avoid perforation the sinus pyriformis while using this device.
- Ensure proper sizing of the CAE within an ETT (**including double-lumen ETT**).
- To avoid barotrauma, ensure that the tip of the CAE is always above the carina, preferably 2-3 cm.
- Use of the Rapi-Fit adapter for oxygenation may be associated with a risk of barotrauma.
- Use of an oxygen source should be considered only if the patient has sufficient egression of the insufflated gas volume.
- Oxygen insufflation may not be appropriate for all patient subgroups and with all products. Please refer to the Catheter Oxygenation Table for product and patient subgroup details.
- If a high-pressure oxygen source is used for insufflation (e.g., jet ventilator), begin at a lower pressure (i.e., 5 psi) and work up gradually. Rising chest wall, pulse oximetry, and oral air flow should be carefully monitored.
- Ensure that the Rapi-Fit Adapter is securely connected to the CAE prior to oxygen delivery. Failure to properly secure the adapter to the CAE may result in hypoxia and serious adverse events.
- Using high-pressure oxygenation with a Luer lock connector in the infant patient subpopulation (> 1 month to 2 years) has only been evaluated with the 8 French catheter; using larger catheters with this subpopulation is not recommended.
- For the 19 French catheter, high-pressure oxygenation with a Luer lock connector should only be used in patients older than 12 years of age. If used in patients 12 years old or younger, the maximum airway pressure may be higher than 28 cm H₂O.

PRECAUTIONS

- The product is intended for use by clinicians trained and experienced in airway management techniques. Standard techniques for placement and exchange of ETTs should be employed.
- This product is not intended for intravascular use.
- Lubricate the catheter and ETT (with a sterile lubricant) before use.
- This device is not recommended for definitive, long-term airway management.

POTENTIAL ADVERSE EVENTS

- Barotrauma
- Pneumothorax
- Hypoxia
- Perforation of the trachea, bronchi, or lung parenchyma
- Catheter migration
- Failed endotracheal tube placement

INSTRUCTIONS FOR USE

1. Before advancing the Cook Airway Exchange Catheter (CAE) into the ETT to be replaced, verify correct ETT position.
2. Using the outer margin of the patient's mouth or nasal orifice as a landmark, note the marking on the ETT. A piece of tape or other marker may be placed on the CAE at the corresponding distance from the tip to aid in correct placement within the ETT.
3. Advance the CAE, sideported end first, into the ETT to be replaced. **(Fig. 1) NOTE:** It is recommended that a sterile medical-grade lubricant, approved for use in the airway, according to hospital protocol, be applied to the orifice of the ETT prior to introduction of the CAE.
4. Properly position the CAE within the ETT by aligning the appropriate centimeter mark on the CAE with the corresponding centimeter mark on the ETT. This placement is determined by visualizing the indicated centimeter length of the ETT, in place, as shown on its surface scale. (For example, an ETT that has been shortened to 24 cm should have the 24 cm marker of the CAE catheter aligned at the 24 cm mark of the ETT.)
5. Fully deflate the cuff of the ETT. While maintaining the position of the CAE, remove the ETT, leaving the CAE in place. **(Fig. 2)**
6. While maintaining the position of the CAE, use the patient's mouth or nasal orifice (depending on approach) as a landmark and advance the new tube over the CAE into appropriate position. **(Fig. 3) NOTE:** It is recommended that a sterile lubricant be applied to the tip of the ETT prior to advancing the ETT.
7. Remove the CAE and inflate the cuff of the new ETT. Verify the new ETT's position using standard methods (e.g., capnography, breath sounds, and chest x-ray). Resume ventilation and secure the new ETT in place. **(Fig. 4)**

Use of Rapi-Fit® Adapter

Rapi-Fit adapters should only be used when oxygen requirements are high and intubation is unsuccessful. Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume. If an oxygen source is used for insufflation, begin at a lower pressure and work up gradually. Observe the chest for outward and inward movements to confirm oxygen insufflation and egression. Pulse oximetry and oral air flow should be carefully monitored as well. In cases of upper airway obstruction, gas discharge from the patient's lungs may require more time.

1. To attach the Rapi-Fit Adapter, position the adapter on the catheter, then push the white collar forward and lock into position. **(Fig. 5)**
2. To remove the adapter, pull the white collar back to release, and then remove from the catheter. **(Fig. 6)**

Catheter Oxygenation

The Cook Airway Exchange Catheters are designed for both positive airway pressure ventilation (Rapi-Fit Adapter – 15 mm connector) and jet ventilation (Rapi-Fit Adapter – Luer lock connector) for both adult and pediatric patients. In the table below, the delivered minute volume and the measured airway pressure are given for jet ventilation in adult and pediatric patients with healthy lung tissue.

Real Part Number (RPN)	Patient Subgroup and Age Range		Delivered Minute Volume ¹ (L/min)	Measured Airway Pressure ¹ (cm H ₂ O)	
				Mean	Maximum
C-CAE-8.0-45-01	Infant	> 1 month to 2 years	0.6	8.9	9.5
	Child	> 2 years to 12 years	0.1	2.1	2.5
	Adolescent	> 12 years through 21 years	0.1	2.2	2.6
	Adult	> 21 years	0.7	3.0	3.2

Real Part Number (RPN)	Patient Subgroup and Age Range		Delivered Minute Volume ¹ (L/min)	Measured Airway Pressure ¹ (cm H ₂ O)	
				Mean	Maximum
C-CAE-11.0-83-01	Infant	> 1 month to 2 years	*	*	*
	Child	> 2 years to 12 years	4.2	8.7	9.4
	Adolescent	> 12 years through 21 years	0.9	4.6	5.2
	Adult	> 21 years	1.3	4.5	4.9
C-CAE-11.0-100-DLT-EF-ST-01	Infant	> 1 month to 2 years	*	*	*
	Child	> 2 years to 12 years	1.0	7.3	7.8
	Adolescent	> 12 years through 21 years	0.6	3.9	4.4
	Adult	> 21 years	0.7	3.9	4.1
C-CAE-14.0-83-01	Infant	> 1 month to 2 years	*	*	*
	Child	> 2 years to 12 years	3.3	23.2	23.8
	Adolescent	> 12 years through 21 years	3.7	11.0	11.3
	Adult	> 21 years	6.3	9.5	10.2
C-CAE-14.0-100-DLT-EF-ST-01	Infant	> 1 month to 2 years	*	*	*
	Child	> 2 years to 12 years	3.0	20.9	21.8
	Adolescent	> 12 years through 21 years	3.3	10.1	10.5
	Adult	> 21 years	5.5	8.6	8.9
C-CAE-19.0-83-01	Infant	> 1 month to 2 years	*	*	*
	Child	> 2 years to 12 years	*	*	*
	Adolescent	> 12 years through 21 years	7.2	19.7	20.3
	Adult	> 21 years	11.2	14.2	14.5

¹ The test conditions were used in an active model. See Testing Conditions table for additional details.

* Testing was not performed within this subgroup with this product.

The presented numbers are mean values. The following testing conditions were used in an active model mode, tested with ASL 5000, Ingmar Medical, Ltd., for each patient subgroup.

Testing Conditions

Input pressure was set to 50 psi for all patient type models.

Patient type	Model Body Weight (kg)	Inspiratory time (s)	Expiratory time (s)	Breaths per minute	Resistance (cm H ₂ O/L/s)	Lung Compliance (mL/cm H ₂ O)
Infant (> 1 month to 2 years)	10	0.46	1.54	30	25	10
Child (> 2 years to 12 years)	20	0.55	1.85	25	20	20
Adolescent (> 12 years through 21 years)	50	0.9	2.1	20	5	40
Adult (> 21 years)	80	1.0	4.0	12	3	100

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from clinicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

This product is protected by U.S. patent 5,052,386.

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© COOK 2018
Bloomington, IN 47404 U.S.A.