COOK

MEDICAL

CE
0088

EN
2

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Bruksanvisning

C _ T _ C A S P _ R E V 4

CATHETER ASPIRATION SET FOR SIMPLE PNEUMOTHORAX

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner.)

DEVICE DESCRIPTION

The Catheter Aspiration Set consists of an introducer needle, wire guide, radiopaque catheter and three-way stopcock.

INTENDED USE

The Catheter Aspiration Set is intended for relief of simple spontaneous or iatrogenic pneumothorax and tension pneumothorax. The product is intended for use by physicians trained and experienced in treatment of pneumothorax. Standard techniques for placement of wire guides and pneumothorax catheters should be employed.

CONTRAINDICATIONS

The product is not recommended for large fluid accumulations or hemothorax.

WARNINGS

Use of the three-way stopcock with this set requires extreme care. Inadvertent opening of stopcock sidearm to air may permit air aspiration.

PRECAUTIONS

None known

INSTRUCTIONS FOR USE

1. Prep the access site with antiseptic solution and drape in standard fashion.
2. Introduce local anesthesia through the skin, subcutaneous tissue, muscle and fascia down to the parietal pleura. Make a 2-3 mm incision through the skin and subcutaneous tissue, but not to the parietal pleura.
3. With a syringe attached to the introducer needle, insert the needle into the pleural cavity. When the needle tip enters the pleural space, remove the syringe and introduce the flexible end of the wire guide into the needle 10-15 cm.
4. Remove the needle, leaving the wire guide in place.
5. Advance the pneumothorax catheter over the wire guide into the pleural cavity. The catheter should be advanced far enough that all catheter sideports are within the pleural space.
6. Remove the wire guide and attach the catheter to suction or a Heimlich valve, using the three-way stopcock if needed.
NOTE: Monitor the catheter tip position fluoroscopically. The distal tip of the catheter should be positioned so that it points toward the extreme apex of the pleural space.
NOTE: All connections must be secure and airtight. Perform inspections of the catheter and connections regularly.
7. Attach the catheter to the chest with adhesive tape and skin suture.
8. Obtain a chest roentgenogram within 24-48 hours before catheter removal.
9. To remove the catheter, remove any catheter skin attachments. Have patient hold his/her breath in full inspiration. Withdraw the catheter and apply an occlusive dressing, such as petrolatum-impregnated gauze.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.