COOK

MEDICAL

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
Fig. 7


Fig. 8

NATHANSON TRANSCYSTIC BILE DUCT STONE EXPLORATION PACK

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Nathanson Transcystic Bile Duct Stone Exploration Pack is packaged preloaded and sterilized. The preloaded assembly includes: **(Fig. 1)**

- 5.5 French, 70 cm long radiopaque nylon catheter with a soft tip and five visual markers spaced 2 cm apart, starting at the tip of the catheter.
- 75 cm long, four-wire Flatwire Stone Extractor and plastic pin vise.
- "Y" Tuohy-Borst Adapter

This set is designed for use through a laparoscopic cholangiography instrument clamp (e.g., Reddick-Olsen Clamp) with a minimum 6 French working channel.

INTENDED USE / INDICATIONS FOR USE

The Nathanson Transcystic Bile Duct Stone Exploration Pack is intended for extraction of bile duct stones via laparoscopic transcystic approach.

CONTRAINDICATIONS

Extraction of bile duct stones greater in diameter than the cystic duct.

WARNINGS

- Rigidity of the catheter assembly may cause cystic duct and bile duct perforation during biliary cannulation, especially when acute inflammation is present or deployment past an impacted stone is encountered.
- Wire deployment within the ampulla of Vater and duodenum may entrap mucosa anchoring the catheter assembly *in situ*, potentially injuring mucosa and triggering pancreatitis.
- This device is designed for the removal of stones and debris ranging from 0.3 cm to 1.0 cm in diameter.
- Common bile duct exploration and alternative surgical intervention may need to be performed if basket becomes impacted after capturing a distal common bile duct stone.
- This device contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in the anatomy of the biliary tree and stone extraction techniques.
- It is important that the surgeon and operating room personnel become completely familiar with the operation of this device prior to use.
- Adequate lead shielding for all operating room staff is required during fluoroscopic X-ray screening.
- Patient lead shielding is necessary if any suspicion of early pregnancy exists.
- Kinking of the catheter's soft tip may obstruct advancement of the extractor.
- Stone extraction via the cystic duct may lead to stone fragmentation.
- Withdrawal of the extractor back into the catheter, in preparation for biliary re-cannulation, may lead to entanglement in the laparoscopic cholangiography instrument clamp's tip mechanism, damaging the wires.
- Avoid deploying the extractor in the ampulla of the distal common duct. **(Fig 6)**

POTENTIAL ADVERSE EVENTS

- Fragmentation of stones in the bile duct
- Ampullary edema
- Bile duct perforation
- Pancreatitis
- Retained bile duct stones
- Retained stones in the peritoneal cavity
- Splitting of cystic duct
- Retained choledocholithiasis
- Biliary leak
- Late biliary strictures

CLINICAL STUDY

A 2007 retrospective study evaluated 160 consecutive patients with bile duct stones (Taylor et al, 2007). The Nathanson Transcystic Bile Duct Stone Exploration Pack (previously known as Nathanson Basket Kit) was used for transcystic exploration when 1 cm diameter stones were observed in the common bile duct. 131 cases (82%) achieved successful laparoscopic common bile duct exploration (LCBDE). 29 cases required additional intervention (i.e., ERCP) to achieve clearance at the time of surgery. Major morbidity occurred in 13.8%, including biliary leak in 7.5% and one late biliary stricture (0.6%).

INSTRUCTIONS FOR USE

Preparation for Use

1. Remove contents from the package in a sterile manner and place on a sterile field.
2. Examine the contents to ensure that all components are intact and assembled correctly.
3. Withdraw the extractor back into the catheter approximately 5 cm by grasping the pin vise and stiff end of the exposed wire and gently pulling backwards. Note the position of the pin vise in relation to the valve. Maintain this position to ensure that the tip of the extractor is not exposed at the distal end of catheter when inserted into the patient.
4. Gently tighten the Tuohy-Borst adapter to ensure that flush does not escape and to stabilize the wire extractor for insertion. Do not over-tighten on the wire.
5. Flush the catheter system with sterile saline through the side-arm adapter to remove any air bubbles.
6. Load the assembly into a pre-positioned, open-tipped laparoscopic cholangiography clamp, ready for use.

Access and Cholangiography

1. Initiate standard laparoscopic exposure with umbilical and three subcostal ports.
2. During dissection of the cystic duct, note the duct size and dissect down closer to the common hepatic duct to avoid spiral valves. **(Fig. 2)**
3. Using an appropriately sized open-ended cholangiography catheter under fluoroscopic imaging, note the cystic duct course into the common bile duct (CBD) (direct, into the right hepatic duct, or very low down near the ampulla) and the CBD diameter, coupled with the size, shape and number of bile duct stones. **NOTE:** Identification of stones in the biliary tree proximal to the cystic duct entry may be reason to choose either choledochotomy or ERCP clearance.
4. If there is uncertainty in interpreting images, it may be helpful to wash contrast out with saline and re-run the cholangiogram. Imaging during the washout phase can be helpful, as during the contrast run, to demonstrate details of the stones.

Transcystic Stone Extraction

1. Once bile duct stones have been confirmed and are considered suitable for transcystic extraction, exchange the *in situ* cholangiography catheter for the pre-prepared Nathanson 5.5 French catheter with pre-loaded extractor, loaded into an appropriate laparoscopic cholangiography instrument clamp. Please refer to the "Preparation for Use" section of these Instructions for Use for information on preparing the catheter.
2. Only after the catheter has passed the stone, loosen the Tuohy-Borst adapter and advance extractor for stone capture.
3. Under fluoroscopic monitoring, advance the catheter tip just beyond the first stone, then advance the extractor through the catheter until its tip emerges. **(Fig. 3)**
4. Hold the extractor's position steady and withdraw the catheter and Tuohy-Borst adapter until the adapter abuts the extractor handle. This will allow the extractor to deploy around the stone. **(Fig 4)**
5. Engage the extractor wires around the stone by manipulating the extractor wires with the extractor handle. Upon capturing the stone, slowly withdraw the extractor and catheter as a unit from the cystic duct. **(Fig. 5)**
6. Place the stone in a suitable area of the abdominal cavity for its later removal, along with the gallbladder. **NOTE:** Calculi may be fragmented upon removal. Further visual inspection with a small-caliber, flexible choledochoscope may be helpful to ensure complete clearance of the stone.
7. After stone removal, patency of CBD can be confirmed by injecting contrast media through the side-arm of the adapter through the catheter and into the CBD.
8. For multiple stones, repeat the process in sequence, working down toward the ampulla. **CAUTION:** Avoid deploying the extractor in the ampulla of the distal common duct, as its mucosa may catch in the extractor wires, potentially inducing acute pancreatitis or resulting in long-term ampullary scarring. Early resistance upon extractor withdrawal and fluoroscopic inversion of the ampulla and the distal bile duct are indicators of entanglement. **(Fig. 7)** Should this occur, discontinue catheter withdrawal and disentangle the ampulla mucosa from the extractor wires by pushing the extractor and the tip of the catheter into the duodenum. **(Fig. 8)** Withdraw the extractor into the catheter, and withdraw the catheter into the bile duct.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

- These instructions for use are based on experience from physicians and (or) their published literature.
- Taylor CJ, Kong J, Ghush M, et al. Laparoscopic bile duct exploration: results of 160 consecutive cases with 2-year follow up. *ANZ J Surg.* 2007;77(6):440-445.
- Nathanson LK, O'Rourke NA, Martin IJ, et al. Postoperative ERCP versus laparoscopic choledochotomy for clearance of selected bile duct calculi: a randomized trial. *Ann Surg.* 2005;242(2):188-192.
- Sgourakis G, Karaliotas K. Laparoscopic common bile duct exploration and cholecystectomy versus endoscopic stone extraction and laparoscopic cholecystectomy for choledocholithiasis: a prospective randomized study. *Minerva Chir.* 2002;57(4):467-474.
- Urbach DR, Khajanchee YS, Jobe BA. Cost-effective management of common bile duct stones: a decision analysis of the use of endoscopic retrograde cholangiopancreatography (ERCP), intraoperative cholangiography, and laparoscopic bile duct exploration. *Surg Endosc.* 2001;15(1):4-13.
- Liu TH, Consorti ET, Kawashima A, et al. Patient evaluation and management with selective use of magnetic resonance cholangiography and endoscopic retrograde cholangiopancreatography before laparoscopic cholecystectomy. *Ann Surg.* 2001;234(1):33-40.
- Traverso LW. A cost analysis of the treatment of common bile duct stones discovered during cholecystectomy. *Semin Laparosc Surg.* 2000;7(4):302-307.
- Cuschieri A, Lezoche E, Morino M, et al. E.A.E.S. multicenter prospective randomized trial comparing two-stage vs single-stage management of patients with gallstone disease and ductal calculi. *Surg Endosc.* 1999;13(10):952-957.
- Rhodes M, Sussman L, Cohen L, et al. Randomised trial of laparoscopic exploration of common bile duct versus postoperative endoscopic retrograde cholangiography for common bile duct stones. *Lancet.* 1998;351(9097):159-161.

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Taylor CJ, Kong J, Ghush M, et al. Laparoscopic bile duct exploration: results of 160 consecutive cases with 2-year follow up. *ANZ J Surg.* 2007;77(6):440-445.

Nathanson LK, O'Rourke NA, Martin IJ, et al. Postoperative ERCP versus laparoscopic choledochotomy for clearance of selected bile duct calculi: a randomized trial. *Ann Surg.* 2005;242(2):188-192.

Sgourakis G, Karaliotas K. Laparoscopic common bile duct exploration and cholecystectomy versus endoscopic stone extraction and laparoscopic cholecystectomy for choledocholithiasis: a prospective randomized study. *Minerva Chir.* 2002;57(4):467-474.

Urbach DR, Khajanchee YS, Jobe BA. Cost-effective management of common bile duct stones: a decision analysis of the use of endoscopic retrograde cholangiopancreatography (ERCP), intraoperative cholangiography, and laparoscopic bile duct exploration. *Surg Endosc.* 2001;15(1):4-13.

Liu TH, Consorti ET, Kawashima A, et al. Patient evaluation and management with selective use of magnetic resonance cholangiography and endoscopic retrograde cholangiopancreatography before laparoscopic cholecystectomy. *Ann Surg.* 2001;234(1):33-40.

Traverso LW. A cost analysis of the treatment of common bile duct stones discovered during cholecystectomy. *Semin Laparosc Surg.* 2000;7(4):302-307.

Cuschieri A, Lezoche E, Morino M, et al. E.A.E.S. multicenter prospective randomized trial comparing two-stage vs single-stage management of patients with gallstone disease and ductal calculi. *Surg Endosc.* 1999;13(10):952-957.

Rhodes M, Sussman L, Cohen L, et al. Randomised trial of laparoscopic exploration of common bile duct versus postoperative endoscopic retrograde cholangiography for common bile duct stones. *Lancet.* 1998;351(9097):159-161.

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Scourakis G, Karaliotas K. Laparoscopic common bile duct exploration and cholecystectomy versus endoscopic stone extraction and laparoscopic

cholecystectomy for choledocholithiasis: a prospective randomized study. *Minerva Chir.* 2002;57(4):467–474.

Urbach DR, Khajanchee YS, Jobe BA. Cost-effective management of common bile duct stones: a decision analysis of the use of endoscopic retrograde cholangiopancreatography (ERCP), intraoperative cholangiography, and laparoscopic bile duct exploration. *Surg Endosc.* 2001;15(1):4–13.

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Cuschieri A, Lezoche E, Morino M, et al. E.A.E.S. multicenter prospective randomized trial comparing two-stage vs single-stage management

of patients with gallstone disease and ductal calculi. *Surg Endosc.* 1999;13(10):952–957.

Rhodes M, Sussman L, Cohen L, et al. Randomised trial of laparoscopic exploration of common bile duct versus postoperative endoscopic retrograde cholangiography for common bile duct stones. *Lancet.* 1998;351(9097):159–161.

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Taylor CJ, Kong J, Ghusn M, et al. Laparoscopic bile duct exploration: results of 160 consecutive cases with 2-year follow up. *ANZ J Surg.* 2007;77(6):440-445.

Nathanson LK, O'Rourke NA, Martin IJ, et al. Postoperative ERCP versus laparoscopic choledochotomy for clearance of selected bile duct calculi: a randomized trial. *Ann Surg.* 2005;242(2):188-192.

Scourakis G, Karaliotas K. Laparoscopic common bile duct exploration and cholecystectomy versus endoscopic stone extraction and laparoscopic cholecystectomy for choledocholithiasis: a prospective randomized study. *Minerva Chir.* 2002;57(4):467-474.

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Urbach DR, Khajanchee YS, Jobe BA. Cost-effective management of common bile duct stones: a decision analysis of the use of endoscopic retrograde

cholangiopancreatography (ERCP), intraoperative cholangiography, and laparoscopic bile duct exploration. *Surg Endosc.* 2001;15(1):4–13.

Liu TH, Consorti ET, Kawashima A, et al. Patient evaluation and management with selective use of magnetic resonance cholangiography and endoscopic retrograde cholangiopancreatography before laparoscopic cholecystectomy. *Ann Surg.* 2001;234(1):33–40.

Traverso LW. A cost analysis of the treatment of common bile duct stones discovered during cholecystectomy. *Semin Laparosc Surg.* 2000;7(4):302–307.

Cuschieri A, Lezoche E, Morino M, et al. E.A.E.S. multicenter prospective randomized trial comparing two-stage vs single-stage management of patients with gallstone disease and ductal calculi. *Surg Endosc.* 1999;13(10):952–957.

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Bloomington, IN 47404 U.S.A.