- EN**
3 Cook Unimpregnated Central Venous Catheter with Power Injection
Instructions for Use
- DA**
10 Cook uimprægneret centralvenekateter med maskininjektion
Brugsanvisning
- DE**
17 Cook nicht imprägnierter zentraler Venenkatheter für Druckinjektionen
Gebrauchsanweisung
- EL**
24 Μη διαποτισμένος κεντρικός φλεβικός καθετήρας της Cook με
δυνατότητα διενέργειας έγχυσης με πίεση
Οδηγίες χρήσης
- ES**
31 Catéter venoso central no impregnado con inyección mecánica Cook
Instrucciones de uso
- FR**
38 Cathéter veineux central non imprégné Cook avec injection automatique
Mode d'emploi
- HU**
45 Cook nyomás alatt injektálható, impregnálatlan centrális vénás katéter
Használati utasítás
- IT**
52 Catetere venoso centrale non impregnato Cook con iniezione automatica
Istruzioni per l'uso
- NL**
59 Cook niet-geïmpregneerde centraalveneuze katheter voor gebruik met
injectiepomp
Gebruiksaanwijzing
- PL**
66 Nieimpregnowany centralny cewnik żylny Cook nadający się do
automatycznego wstrzykiwania
Instrukcja użycia
- PT**
73 Cateter venoso central Cook não impregnado com injeção elétrica
Instruções de utilização
- SV**
80 Cook icke-impregnerad central venkateter med tryckinjektion
Bruksanvisning

COOK UNIMPREGNATED CENTRAL VENOUS CATHETER WITH POWER INJECTION

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Cook Unimpregnated Central Venous Catheters with power injection are available in various French sizes, lengths, and lumen configurations. Multi-lumen catheters incorporate separate, non-communicating vascular access lumens within a single catheter body. **The device and its components do not contain natural rubber latex.**

- 4.0 French catheters are recommended for patients from birth and older (or weighing at least 2.4 kilograms).
- 5.0 French catheters are recommended for patients aged 2 years and older (or weighing at least 10 kilograms).
- 7.0 French and 8.0 French catheters are recommended for patients aged 10 years and older (or weighing at least 50 kilograms).
- 9.0 French and 10.0 French catheters are recommended for patients aged 21 years and older (adults only).

INTENDED USE

The Cook Unimpregnated Central Venous Catheter with power injection is used for:

1. Continuous or intermittent drug infusions
2. Central venous blood pressure monitoring (CVP)
3. Acute hyperalimentation
4. Blood sampling
5. Delivery of whole blood or blood products
6. Simultaneous, separate infusion of drugs for multi-lumen catheters only
7. Power injection of contrast media*

*The flow rate may not exceed 3 mL/sec for 4.0 and 5.0 French catheters and 10 mL/sec for 7.0, 8.0, 9.0, and 10.0 French catheters. Verify prior to use that the maximum safety cut-off pressure limit is set at or below 250 psi for 4.0 and 5.0 French catheters and 325 psi for 7.0, 8.0, 9.0 and 10.0 French catheters.

The device is a short-term use catheter, intended for less than 30 days.

The 9.0 and 10.0 French catheters include an inner catheter to facilitate insertion of the main catheter.

The Cook dilator is used for dilating puncture sites or catheter tracts.

The Cook Unimpregnated Central Venous Catheter with power injection is intended for adult and pediatric populations.

CONTRAINDICATIONS

The use of this device is contraindicated when:

- The presence of device-related infection, bacteremia, or septicemia is known or suspected.
- The patient's body size is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device.
- Severe chronic obstructive lung disease exists (percutaneous subclavian placement only).
- There has been past irradiation of prospective insertion site.
- There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.
- There are local tissue factors that may prevent proper device stabilization and/or access.

WARNINGS

- Complications arising from the use of this device can result in serious injury or death; the catheter tip can erode or perforate vascular walls.
- Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of central venous system. Tip position should be verified by x-ray and monitored on a routine basis. Periodic lateral view x-ray is suggested to assess tip location in relation to vessel wall. Tip position should appear to be parallel to vessel wall.
- To avoid vascular injury, a short needle length should be used for initial access, especially in pediatric populations.
- To avoid vascular injury, do not use excessive force when advancing dilators. Use the smallest size dilator catheter placement will allow, especially in pediatric populations. The wire guide must always lead the dilator by several centimeters. Do not advance the dilator more than a few centimeters into the vessel.
- The safe and effective use of central venous catheters with power injector pressures (safety cut-off) set above 250 psi or 325 psi have not been established for 4.0 and 5.0 French or 7.0, 8.0, 9.0, and 10.0 French, respectively.
- **Do not power inject if maximum injection rate cannot be verified to be limit of 3 mL/sec for 4.0 and 5.0 French, and 10 mL/sec for 7.0, 8.0, 9.0, and 10.0 French.**
- Use only the distal lumen for power injection. The safe and effective use of power injectors with the proximal or mid lumen has not been established.
- To safely use catheters with a power injector, the technician/health care professional must verify prior to use that the maximum safety cut-off pressure limit is set at or below 250 psi and maximum flow rate of 3 mL/sec for 4.0 and 5.0 French and 325 psi and maximum flow rate of 10 mL/sec for 7.0, 8.0, 9.0, and 10.0 French.
- Do not exceed maximum flow rate of 3 mL/sec for 4.0 and 5.0 French and 10 mL/sec for 7.0, 8.0, 9.0, and 10.0 French. Exceeding the maximum flow rate may result in catheter failure and/or catheter tip displacement.

- Power injector machine pressure limiting (safety cut-off) settings may not prevent over-pressurization of an occluded catheter.
- This catheter's indication for power injection of contrast media implies the catheter's ability to withstand the procedure, but does not imply appropriateness of the procedure for a particular patient. A suitably trained clinician is responsible for evaluating the health status of a patient as it pertains to a power injection procedure.
- To distend great vessels and to prevent inadvertent air aspiration during catheter insertion, patient should be placed in Trendelenburg position.
- Do not place the CVC into or allow it to remain in the right atrium or right ventricle. Verify the catheter tip position using radiography or appropriate technology based on hospital protocol.
- In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.
- When placed via the femoral vein, the catheter tip should be in the inferior vena cava (IVC) above the level of the diaphragm. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.
- Avoid entanglement or contact of wire guide with previously implanted devices. Review and check patient history and/or anatomy prior to procedure.
- Be sure to maintain wire guide integrity by avoiding kinking and aggressive usage.
- Discard and manage sharps devices (needles, scalpels, etc.) safely during and after procedure.
- Be sure to avoid inadvertent arterial puncture, especially in pediatric populations.

PRECAUTIONS

- The catheter is intended for use by healthcare professionals trained and experienced in the placement of central venous catheters using percutaneous entry (Seldinger) technique. Standard Seldinger technique for placement of percutaneous vascular access sheaths, catheters, and wire guides should be employed during the placement of a central venous catheter.
- Do not re-sterilize the catheter.
- Do not cut, trim or modify the catheter or components prior to placement or intraoperatively.
- Patient movement can cause catheter tip displacement, especially in pediatric populations. Use should be limited to controlled hospital situations. Catheters placed via an antecubital vein have shown forward tip movement of up to 10 cm with motion of the extremity. Catheters placed from either a jugular or subclavian vein have demonstrated forward tip movement of 1-3 cm with neck and shoulder motion.
- The catheter should not be used for long-term indwelling applications.
- **If lumen flow is impeded, do not force injection or withdrawal of fluids. Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure. Notify attending physician immediately.**
- The catheter should not be used for chronic hyperalimentation.
- Select puncture site and length of catheter needed by assessing patient anatomy and condition.
- The catheter should not be used for patients of low birth weight to extremely low birth weight. **NOTE:** The catheter is intended for patients weighing at least 2.4 kilograms.
- The recommended minimum age/weight for each catheter French size is established based on published literature. Physicians should ultimately decide which catheter size and configuration is to be used in patients.
- The left subclavian and left jugular veins should be used only when other sites are not available.
- Use of ECG, ultrasound and/or fluoroscopy is suggested for accurate catheter placement, especially in pediatric populations.
- Radiographic confirmation of catheter placement prior to each power injection procedure is recommended.
- If a tray contains a lidocaine label, the label should be used to label the syringe for lidocaine injection. Any unused lidocaine should be disposed of at the conclusion of the procedure.
- Do not cut the wire guide to alter the length.
- **Do not withdraw the wire guide against needle bevel;** doing so could result in severing or damaging the wire guide.
- Do not use excessive force when introducing wire guide or tissue dilator as this can lead to vessel perforation and bleeding, especially in pediatric populations.
- Passage of wire guide into the right heart can cause dysrhythmias, right bundle branch block, and a perforation of vessel, atrial, or ventricular wall.
- Do not apply excessive force during removal of the wire guide. When withdrawal cannot be easily accomplished; a visual image should be obtained and further consultation requested.
- When making a skin incision, position the blade of the scalpel away from the wire guide.

POTENTIAL ADVERSE EVENTS

- Catheter perforation / displacement / leakage / occlusion
- Bleeding at catheter insertion site
- Arterial puncture
- Infiltration
- Pleural effusion
- Migration
- Infection
- Sepsis
- Hematoma
- Supraventricular tachycardia

- Pneumothorax
- Hemothorax
- Vein laceration
- Thrombosis
- Phlebitis

PRODUCT RECOMMENDATIONS

Catheter Size*	Maximum Power Injection Flow Rate (mL/sec)	Injector Safety Cut-off Pressure (psi)	Average Maximum Internal Catheter Pressure During Maximum Flow Rate (psi)**	Average Maximum Static Burst Pressure (psi) ***
4 French Double-Lumen	3	250	98.1	332.3
5 French Double-Lumen	3	250	217.8	382.9
5 French Triple-Lumen	3	250	197.6	342.8
7 French Triple-Lumen	10	325	143.3	181.6
8 French Double-Lumen	10	325	55.1	208.1
9 French Triple-Lumen	10	325	35.8	162.8
10 French Five-Lumen	10	325	32.1	228.4

*All power injection testing was performed on the distal lumen of each catheter.

**Pressures determined using room temperature Omnipaque 300™ or Visipaque 270™ contrast and verified using a Medrad Stellant™ CT or Crescent injector system. Omnipaque 300™ has a viscosity of 11.8 centipoise at 20 degrees C. Visipaque 270™ has a viscosity of 8.8 centipoise at 20 degrees C. A change in temperature or viscosity of the contrast medium used will result in a change in achievable flow rates. Omnipaque 300™ and Visipaque 270™ are registered trademarks of GE Healthcare.

*** Maximum Static Burst Pressure is the static burst pressure failure point of the catheter. When the catheter was occluded, failure occurred at the specified pressure.

Catheter Size and Puncture Site

Preliminary reports indicate that catheter size can influence clotting; larger diameter catheters have more tendency to promote clots. As reported by Amplatz and others, clot formation has less relation to type of catheter material than to size of catheter.¹ The angle of the catheter tip to the vessel wall should be checked carefully. Gravenstein and Blackshear in an *in vitro* study confirmed that perforation was significantly more likely when the angle of incidence between catheter and insertion surface was greater than 40 degrees than when it was 40 degrees.²

Another critical factor that can cause a catastrophic event is the choice of puncture site. Findings by Tocino and Watanabe indicate that the left subclavian and left jugular veins should be avoided when practical.³ Eighty percent of the perforations or erosions were found when these vessels were used. In addition, they have observed that the tip curve of a wedged catheter can be detected with lateral view x-ray. Findings by de Jonge indicate the femoral, jugular, and subclavian veins are the most frequently used sites for CVC insertion in children. The femoral vein is the preferential insertion site for quick central venous access in non-sedated children if a relatively short (<5 days) indwelling time is expected. The jugular vein is the preferred site for CVC access in sedated children with a controlled airway. For the subclavian vein, a left-sided approach is recommended, providing there are no contraindications.⁴

The above discussion is meant to be a guide for catheter size and puncture site. As more data become available, other causal factors may become evident, but present information suggests that:

- The catheter size should be as small as the use will allow.
- Left subclavian and left jugular veins should be used only when other sites are not available.

The following variables should be considered when selecting appropriate catheter and length for both adult and pediatric patients:

- Patient history
- Size and age of patient
- Access site available
- Unusual anatomical variables
- Proposed use and duration of treatment plan

For pediatric patients, especially in an emergent event, refer to standard weight-based tape (Broselow Tape) for infusion of fluids and medications.

CATHETERS AVAILABLE FOR USE

Depth Marker Information

Catheters smaller than 6.0 French have markers indicated by dashes and dots on the shaft to indicate depth. The 5 cm long catheter has dash marks starting at 1 cm from the distal tip and continuing at 1 cm increments. The 6 cm long catheter has dash marks starting at 1 cm from the distal tip and continuing at 1 cm increments, until the 5 cm mark. The 5 cm mark is designated with a single dot, denoting 5 cm from the distal tip. For catheters 8 cm and longer, the 5 cm dot mark is the first depth marker on the catheter. Dash marks continue past the 5 cm mark, continuing at 1 cm increments. Each dot configuration represents a unique distance from the distal tip.


Dot Number	Distance from Distal Tip
1	5 cm
2	10 cm
3	15 cm

Catheters 7.0 French or larger have markers indicated by dashes or the respective numeric markers.

Double-Lumen Information (Polyurethane Catheter)


The 4.0 and 5.0 French Double-Lumen Central Venous Catheters have circular lumen configurations.

The 8.0 French Double-Lumen Central Venous Catheter has a "double-D" lumen configuration.

Catheter Fr	Cross-section	Lumen No. / Port	Equivalent Gage	Lumen Volume
4.0		#1 Distal #2 Proximal	20 23	0.3 mL 0.3 mL
5.0		#1 Distal #2 Proximal	20 20	0.2 mL 0.2 mL
8.0		#1 Distal #2 Proximal	14 14	0.9 mL 1.0 mL

Triple-Lumen Information (Polyurethane Catheter)

The 5.0, 7.0, and 9.0 French Triple-Lumen Central Venous Catheters have circular lumen configurations.

Catheter Fr	Cross-section	Lumen No. / Port	Equivalent Gage	Lumen Volume
5.0		#1 Distal #2 Mid #3 Proximal	18 23 23	0.3 mL 0.2 mL 0.2 mL
7.0		#1 Distal #2 Mid #3 Proximal	16 18 18	0.6 mL 0.5 mL 0.6 mL
9.0		#1 Distal #2 Mid #3 Proximal	14 18 18	0.9 mL 0.4 mL 0.5 mL

Five-Lumen Information (Polyurethane Catheter)

The 10.0 French Five-Lumen Central Venous Catheters have circular lumen configurations.

Catheter Fr	Cross-section	Lumen No. / Port	Equivalent Gage	Lumen Volume
10.0		#1 Distal #2 Mid #3 Mid #4 Proximal #5 Proximal	14 17 17 19 19	0.9 mL 0.4 mL 0.4 mL 0.2 mL 0.2 mL

Suggested Lumen Utilization: Double-Lumen

- **#1 Distal exit port (endhole)** — whole blood or blood product delivery and sampling; any situation requiring more flow rate; CVP monitoring; medication delivery; power injection studies. **It is strongly recommended that this lumen be used for all blood sampling.** "CT" labeled on the distal #1 lumen hub indicates that this is the lumen that should be used for power injection.
- **#2 Proximal exit port** — medication delivery; acute hyperalimentation.

Suggested Lumen Utilization: Triple-Lumen

- **#1 Distal exit port (endhole)** — whole blood or blood product delivery and sampling; any situation requiring more flow rate; CVP monitoring; medication delivery; power injection studies. **It is strongly recommended that this lumen be used for all blood sampling.** "CT" labeled on the distal #1 lumen hub indicates that this is the lumen that should be used for power injection.

- #2 Middle exit port — medication delivery; acute hyperalimentation.
- #3 Proximal exit port — medication delivery.

Suggested Lumen Utilization: Five-Lumen

- #1 Distal exit port (endhole) — whole blood or blood product delivery and sampling; any situation requiring more flow rate; CVP monitoring; medication delivery; power injection studies. **It is strongly recommended that this lumen be used for all blood sampling.** "CT" labeled on the distal #1 lumen hub indicates that this is the lumen that should be used for power injection.
- #2 Middle exit port — medication delivery; acute hyperalimentation.
- #3 Middle exit port — medication delivery; acute hyperalimentation.
- #4 Proximal exit port — medication delivery.
- #5 Proximal exit port — medication delivery.



Suggested Catheter Maintenance

- The catheter entry site must be prepared and maintained in a manner consistent with standard procedure for central venous catheterization.
- To prevent clotting or possibility of air embolus, the double-lumen's #2 lumen, the triple-lumen's #2 and #3 lumens, and the five-lumen's #2, #3, #4, and #5 lumens should be filled with saline solution or heparinized saline solution (100 units of heparin per mL of saline is usually adequate), depending on institutional protocol, prior to catheter introduction.
- After the catheter is placed and prior to use, tip position and lumen patency should be confirmed by free aspiration of venous blood. **If blood is not freely aspirated, the physician should immediately reevaluate catheter tip position. Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure.**
- Any unused lumens should be maintained with continuous saline or heparinized saline drip or locked with heparinized saline solution. Heparin-locked lumens should be reestablished at least every 8 hours.
- Before using any lumen already locked with heparin, the lumen should be flushed with twice the indicated lumen volume using normal saline. Lumens should be flushed with normal saline between administrations of different infusates. After use, the lumen should again be flushed with twice the indicated lumen volume using normal saline before reestablishing heparin lock.
- Strict aseptic technique must be adhered to while using and maintaining the catheter.

INSTRUCTIONS FOR USE

Access Sites of Choice

1. High Internal Jugular
2. External Jugular
3. Low Internal Jugular
4. Supraclavicular
5. Infraclavicular
6. Femoral

- The instructions for use are applicable for both pediatric and adult patient populations. Catheter size with intended length and patient age, or weight is provided.

Patient		Catheter	
Age	Weight (kg)	French	Length (cm)
Birth and older	≥2.4kg	4.0	5-15
≥2 years	≥10kg	5.0	5-15
≥10 years	≥50kg	7.0-8.0	15-26.5
≥21 years	Adult	9.0-10.0	15-26.5

- If applicable, remove the Luer-lock end cap from each extension.
- Prepare the catheter for insertion by flushing each of the lumens and clamping or attaching the injection caps to the appropriate extensions. Leave the distal extension uncapped for wire guide passage.
- Introduce the thinwall percutaneous entry needle into the vessel. Venous blood should be easily aspirated to confirm the position of the needle tip within the vessel. (Fig. 1)
- Straighten the wire guide using the J-tip straightener. Pass the straightened wire guide through the needle; advance the wire guide 5-10 cm into the vessel. If a straight wire is used, always advance the soft, flexible end through the needle hub and into the vessel. (Fig. 2) **If you encounter resistance during the wire guide insertion, do not force the wire guide. Withdrawal of the wire guide through the needle should be avoided; breakage may result.**
- While maintaining wire guide position, withdraw the needle and the wire guide straightener. (Fig. 3)
- Enlarge puncture site with number 11 scalpel blade. (Fig. 4)
CAUTION: To avoid vascular injury, do not use excessive force when advancing dilators. Use the smallest size dilator catheter placement will allow. The wire guide must always lead the dilator by several centimeters. Do not advance the dilator more than a few centimeters into the vessel.
- The dilator can be advanced over the wire guide and removed prior to insertion of the central venous catheter.
- Measure the catheter to be used against the patient to determine the approximate length of the catheter needed from puncture site to the central venous tip position. For pediatric populations, weight and length of patient can also be used to confirm length of the catheter needed.⁵
- Introduce the central venous catheter over the wire guide. While maintaining the wire guide position, advance the catheter into the vessel with a gentle twisting motion. (Fig. 5)
NOTE: Do not advance the catheter tip beyond the distal tip of the wire guide. Always have the wire guide leading during catheter placement. Verify the catheter tip position using radiography or appropriate technology.
- After the catheter is in position, remove the wire guide and inner catheter, as applicable. The 9 French and 10 French central venous catheters include an inner catheter in the main distal lumen that should be removed with the wire guide. (Fig. 6) **Venous blood should be easily aspirated.** Lumens should now be flushed with 5-10 mL normal saline prior to use or establishment of heparin lock. The winged hub can now be sutured into place. If the catheter is not introduced to its full length, additional suture should be carefully placed around the catheter and affixed to the skin at the entry site (if a movable suture wing is not included). This will help prevent backward or forward catheter movement.
NOTE: A wire guide that is at least twice as long as the catheter is recommended for catheter exchange procedure.

POWER INJECTION PROCEDURE

- Confirm proper catheter tip position radiographically prior to injection.
- Remove any injection/needleless caps from the catheter.
- Attach a 10 mL (or larger) syringe filled with sterile normal saline to the hub of the distal (#1) extension tube to be used for power injection.
- Ensure adequate blood return and flush the catheter vigorously with the entire 10 mL of sterile normal saline to ensure lumen patency.
WARNING: Failure to ensure patency of the catheter lumen prior to injection may result in catheter failure.
- Remove the syringe and attach the power injection device to the catheter using the manufacturer's recommendations.
NOTE: Refer to the power injector manufacturer's recommendations regarding connecting tubes.
- Conduct a study using the power injector, making sure not to exceed the maximum flow rate of 3 mL/sec or pressure limit safety cut-off of 250 psi for 4.0 and 5.0 French catheters or the maximum flow rate of 10 mL/sec or pressure limit safety cut-off of 325 psi for 7.0, 8.0, 9.0, and 10.0 French catheters.
- Disconnect the power injection device and flush the catheter again with 10 mL of sterile normal saline.
- Replace the injection/needleless cap on the catheter.
- Confirm proper catheter tip position radiographically following power injection.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

1. Amplatz K. A simple non-thrombogenic coating. *Invest Radiol.* 1971;6(4):280–289.
2. Gravenstein N, Blakeshear RH. In vitro evaluation of relative perforating potential of central venous catheters: comparison of materials, selected models, number of lumens, and angles of incidence to simulated membrane. *J Clin Monit.* 1991;7(1):1–6.
3. Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenol.* 1986;146(3):487–490.
4. de Jonge RCJ, et al. Central venous catheter use in the pediatric patient: mechanical and infectious complications. *Ped Crit Care Med.* 2005;6(3):329–339.
5. Andropoulos DB, et al. Recommended length of central venous catheter (CVC) insertion in pediatric patients based on weight. *Anesth Analg.* 2001;93:883–886.

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1. Amplatz K. A simple non-thrombogenic coating. *Invest Radiol.* 1971;6(4):280–289.
2. Gravenstein N, Blackshear RH. In vitro evaluation of relative perforating potential of central venous catheters: comparison of materials, selected models, number of lumens, and angles of incidence to simulated membrane. *J Clin Monit.* 1991;7(1):1–6.
3. Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenol.* 1986;146(3):487–490.
4. de Jonge RCJ, et al. Central venous catheter use in the pediatric patient: mechanical and infectious complications. *Ped Crit Care Med.* 2005;6(3):329–339.
5. Andropoulos DB, et al. Recommended length of central venous catheter (CVC) insertion in pediatric patients based on weight. *Anesth Analg.* 2001;93:883–886.

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- Amplatz K. A simple non-thrombogenic coating. *Invest Radiol.* 1971;6(4):280-289.
- Gravenstein N, Blackshear RH. In vitro evaluation of relative perforating potential of central venous catheters: comparison of materials, selected models, number of lumens, and angles of incidence to simulated membrane. *J Clin Monit.* 1991;7(1):1-6.
- Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenol.* 1986;146(3):487-490.
- de Jonge RCJ, et al. Central venous catheter use in the pediatric patient: mechanical and infectious complications. *Ped Crit Care Med.* 2005;6(3):329-339.
- Andropoulos DB, et al. Recommended length of central venous catheter (CVC) insertion in pediatric patients based on weight. *Anesth Analg.* 2001;93:883-886.

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- Amplatz K. A simple non-thrombogenic coating. *Invest Radiol.* 1971;6(4):280–289.
- Gravenstein N, Blakeshear RH. In vitro evaluation of relative perforating potential of central venous catheters: comparison of materials, selected models, number of lumens, and angles of incidence to simulated membrane. *J Clin Monit.* 1991;7(1):1–6.
- Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenol.* 1986;146(3):487–490.
- de Jonge RCJ, et al. Central venous catheter use in the pediatric patient: mechanical and infectious complications. *Ped Crit Care Med.* 2005;6(3):329–339.
- Andropoulos DB, et al. Recommended length of central venous catheter (CVC) insertion in pediatric patients based on weight. *Anesth Analg.* 2001;93:883–886.

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Bloomington, IN 47404 U.S.A.