EN 3	VersaTube™ Tapered Tracheostomy Tube with Disposable Inner Cannula Instructions for Use
DA 6	VersaTube™ tilspidset trakeostomitube med engangsinderkanyle Brugsanvisning
DE 10	VersaTube™ konisch zulaufende Trachealkanüle mit Einweg-Innenkanüle Gebrauchsanweisung
EL 15	Κωνικός σωλήνας τραχειοστομίας VersaTube™ με αναλώσιμη εσωτερική κάνουλα Οδηγίες χρήσης
ES 19	Tubo de traqueostomía cónico VersaTube™ con cánula interior desechable Instrucciones de uso
FR 24	Canule de trachéotomie effilée VersaTube™ avec canule interne jetable Mode d'emploi
IT 28	Cannula tracheostomica rastremata VersaTube™ con cannula interna monouso Istruzioni per l'uso
NL 33	VersaTube™ tapse tracheostomietube met disposable binnencanule Gebruiksaanwijzing
PT 37	Tubo de traqueostomia cónico VersaTube™ com cânula interior descartável Instruções de utilização
SV 41	VersaTube™ avfasad trakeostomikanyl med innerkanyl för engångsbruk Bruksanvisning
ZH 45	带一次性内套管的 VersaTube™ 锥头气管造口导管 使用说明


Fig. 1



Fig. 2

VERSATUBE™ TAPERED TRACHEOSTOMY TUBE WITH DISPOSABLE INNER CANNULA

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The VersaTube Tapered Tracheostomy Tube consists of an outer cannula, inner cannula, obturator, and tracheostomy tube strap.

Tracheostomy Tube ID/OD (mm)	Tracheostomy Tube Length "C" (mm)	Angle	Inner Cannula ID (mm) / Tracheostomy Tube Color Code	Cuff Resting Diameter (mm)
7*/10	78.0	96°	6/Green	25
8*/11	86.0	96°	7/White	28
9*/12	98.0	97°	8/Blue	30

* 7 mm Tracheostomy Tube color code is green, 8 mm color code is white, and 9 mm color code is blue.

NOTE: Tracheostomy Tube ID is the inside diameter of the outer cannula. OD is the outside diameter of the outer cannula (not including cuff features). Length is the distance from the back of the neck plate to the distal tip of the tube. Cuff OD is the outer diameter of the cuff at the resting diameter.

INTENDED USE

The VersaTube Tapered Tracheostomy Tube is intended to provide an artificial airway to establish airway patency and to provide maintenance of the airway.

CONTRAINDICATIONS

None known

WARNINGS

- When using a laser or electrosurgery electrodes in conjunction with this or other tracheostomy tubes, ensure that the laser beam or electrodes do not

- come into contact with the tube. Such contact, especially in the presence of oxygen-enriched mixtures, could result in rapid combustion of the tube.
- The inner cannula included in this package should only be replaced with a VersaTube disposable inner cannula of the same size. Use of an incorrectly sized inner cannula may result in restricted gas flow, tracheal damage or occlusion, or accumulation of secretions at the tip of the tube which can lead to infection, blockage, or airway obstruction.
 - Do not use an inner cannula when using a tracheostomy tube in conjunction with a right-angled connector that incorporates a coaxial fresh gas tube which protrudes into the connector lumen (e.g., Norman elbow). The protruding fresh gas tube could occlude the inner cannula, restricting expiratory flow of gas and resulting in barotrauma/hypoxia.
 - Do not plug this tracheostomy tube with a decannulation plug; doing so may result in insufficient airway for the patient, even when the cuff is completely deflated.
 - Cuff pressures should be monitored and adjusted routinely, and should never exceed 25mm Hg. Over-inflation of the cuff may inhibit ventilation and may also result in permanent damage to the trachea. Underinflation of the cuff may result in aspiration of subglottic secretions leading to lung infections.
 - Use of this tracheostomy tube in patients with abnormal upper airway anatomy or pathology of the upper airway is not recommended; compression of the outer cannula may result in partial or total airway obstruction.
 - Verify tracheostomy tube position via bronchoscopic visualization and/or chest X-ray to ensure correct placement. Incorrect placement may result in trauma to the trachea or respiratory obstruction.
 - The products are considered sterile unless the protective package has been opened, damaged or broken. Do not attempt to resterilize products.

PRECAUTIONS

- This product is designed to be used in Percutaneous Dilational Tracheostomy (PDT) procedures and Balloon Assisted Tracheostomy procedures. Standard PDT or Balloon Assisted Tracheostomy techniques should be employed. Sterile technique should be adhered to for the handling and insertion of the tracheostomy tube.
- The VersaTube Tracheostomy Tube is classified as a single use device. The manufacturer recommends that the tracheostomy tube usage does not exceed twenty-nine days.
- For ventilator-dependent patients, a replacement tracheostomy tube should be kept at bedside and cuff inflation should be monitored regularly.
- Insertion of this tracheostomy tube in morbidly obese patients may result in tube obstruction or false passage of the tube.
- If lubricant is used prior to tracheostomy tube insertion, ensure that lubricant does not occlude the tube's lumen; lumen occlusion may prevent adequate ventilation.
- Patients breathing gas should be adequately humidified to minimize encrustation of the tracheostomy tube and/or inner cannula lumen and prevent mucosal damage.

- The inner cannula should be checked and replaced at regular intervals to avoid blockage which may reduce or obstruct airway lumen.
- Do not reposition the *in situ* tracheostomy tube while the cuff is inflated; doing so may result in trauma to the stoma and/or trachea.
- During anesthesia, nitrous oxide may diffuse into the cuff, causing an initial increase, with subsequent decrease, in cuff pressure. This can result in pressure trauma to the trachea and loss of airway seal, resulting in aspiration of subglottic secretions into the lungs.
- Prior to removal of cuffed tracheostomy tubes, ensure that all air is completely removed from the cuff; failing to do so may result in pain to the patient and damage to the trachea and/or stoma.
- During and after attachment of the breathing system to the tracheostomy tube connector, avoid excessive rotational or linear forces on the tube; application of such forces may result in tube kinking, disconnection, or occlusion.

PRODUCT RECOMMENDATIONS

The VersaTube Tapered Tracheostomy Tube is compatible with Ciaglia, Ciaglia Blue Rhino and the Ciaglia Blue Rhino G2 Percutaneous Tracheostomy Introducer Sets and Trays. The VersaTube Tapered Tracheostomy Tube is also compatible with the Dolphin Balloon Assisted Tracheostomy Introducer. Refer to the table below for sizing information.

Tracheostomy Tube Size (ID)	Compatible Loading Dilator
7.0 mm	21 French
8.0 mm	24 French
9.0 mm	27 French

INSTRUCTIONS FOR USE

Tube Preparation

1. Select the appropriately sized tracheostomy tube.
2. Test the cuff and inflation system for leakage.
3. With the cuff inflated, taper the cuff back by gently moving away from the distal tip of the outer cannula toward the swivel neck plate as the residual air is removed by deflation. This will ease insertion and guard against cuff perforation by sharp cartilage edges. **NOTE:** Do not use sharp instruments such as forceps or hemostats to taper the cuff; doing so may result in cuff damage.
4. Ensure that the cuff is fully deflated.
5. Remove inner cannula. (**Fig. 1**)

Insertion

6. With the inner cannula removed from the tube, insert an appropriately sized loading dilator into the tracheostomy tube. The tapered section of the loading dilator should extend from the distal tip of the tracheostomy tube by approximately 2 cm. **NOTE:** Water-soluble lubricant can be applied to the outer cannula, cuff and protruding portion of loading dilator to facilitate insertion.

7. Perform the dilational tracheostomy procedure per the introducer system manufacturer's instructions.
8. Using aseptic technique, insert the inner cannula into position. **(Fig. 2)**
NOTE: Two sterile inner cannulae are provided. The outer cannula can be used independently (without the inner cannula in place).
9. Inflate the low-pressure cuff by injecting air into the Luer valve of the inflation line using a syringe. Selection of a cuff inflation and deflation method is at the discretion of the physician/hospital.
10. Switch the breathing apparatus from the endotracheal tube to the tracheostomy tube and verify adequate ventilation.
11. Secure the tracheostomy tube per hospital protocol. **NOTE:** A tracheostomy tube strap is included with this device.

CAUTION: This product is composed of soft materials to conform to tracheal tissue for performance and patient comfort. Simple precautions in handling of this tube during insertion and while in place will facilitate proper function and minimize tears and breaks in the inflation system. Avoid pulling or manipulation of the inflation line, as it is designed to conduct and hold air as part of the cuff inflation system. It is recommended that the inflation line be maintained in a position allowing for patient mobility without placing tension on the line-to-cannula junction.

Cuff Deflation and Tube Removal

12. Before removing the tube, the cuff should be completely deflated. This will ensure that the cuff passes through the stoma with minimal resistance. **NOTE:** Accumulated secretions above the cuff may need to be suctioned before deflating the cuff (with a syringe), unless suctioning is contraindicated.
13. To deflate the cuff, withdraw air slowly from the Luer valve of the inflation line using a syringe. After the cuff is fully deflated, slowly remove the tracheostomy tube.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.