

|                        |                                                                                                                 |
|------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>EN</b><br><b>4</b>  | <b>Intraosseous Infusion Needles</b><br><b>For Pediatric Use (&lt; 24 months)</b><br>Instructions For Use       |
| <b>DA</b><br><b>6</b>  | <b>Intraossøse infusionskanyler til pædiatrisk brug</b><br><b>(&lt; 24 måneder)</b><br>Brugsanvisning           |
| <b>DE</b><br><b>8</b>  | <b>Intraossäre Infusionskanülen für die Pädiatrie</b><br><b>(&lt; 24 Monate)</b><br>Gebrauchsanweisung          |
| <b>EL</b><br><b>10</b> | <b>Βελόνες ενδοστικής έγχυσης για παιδιατρική χρήση</b><br><b>(&lt; 24 μηνών)</b><br>Οδηγίες χρήσης             |
| <b>ES</b><br><b>12</b> | <b>Agujas de infusión intraósea para uso pediátrico</b><br><b>(menores de 24 meses)</b><br>Instrucciones de uso |
| <b>FR</b><br><b>15</b> | <b>Aiguilles de perfusion intraosseuse à usage</b><br><b>pédiatrique (&lt; 24 mois)</b><br>Mode d'emploi        |
| <b>IT</b><br><b>17</b> | <b>Aghi per infusione intraossea per uso pediatrico</b><br><b>(&lt; 24 mesi)</b><br>Istruzioni per l'uso        |
| <b>NL</b><br><b>19</b> | <b>Intraossale infusienaalden voor pediatrisch gebruik</b><br><b>(&lt; 24 maanden)</b><br>Gebruiksaanwijzing    |
| <b>NO</b><br><b>21</b> | <b>Intraossøse infusjonsnåler til pædiatrisk bruk</b><br><b>(&lt; 24 måneder)</b><br>Bruksanvisning             |
| <b>PT</b><br><b>23</b> | <b>Agulhas de infusão intra-óssea para uso pediátrico</b><br><b>(&lt; 24 meses)</b><br>Instruções de utilização |
| <b>SV</b><br><b>26</b> | <b>Intraosseösa infusionsnålar för pædiatrisk</b><br><b>användning (&lt; 24 månader)</b><br>Bruksanvisning      |


**Fig. 1**

- |                                                                      |                                                                                 |                                                                                |
|----------------------------------------------------------------------|---------------------------------------------------------------------------------|--------------------------------------------------------------------------------|
| 1. Femur                                                             | 1. Fémur                                                                        | 1. Femur                                                                       |
| 2. Tibia                                                             | 2. Tibia                                                                        | 2. Tibia                                                                       |
| 3. Access site is 1-2 cm distal to midpoint                          | 3. El lugar de acceso está a 1-2 cm en dirección distal respecto al punto medio | 3. Tilgangsstedet er 1-2 cm distalt for midtpunktet                            |
| 1. Femur                                                             | 1. Fémur                                                                        | 1. Fémur                                                                       |
| 2. Tibia                                                             | 2. Tibia                                                                        | 2. Tibia                                                                       |
| 3. Indførsingssted er 1-2 cm distalt for midtpunktet                 | 3. Le site d'accès se trouve de 1 à 2 cm en aval du point central               | 3. O local de acesso situa-se a 1 a 2 cm distalmente em relação ao ponto médio |
| 1. Femur                                                             | 1. Femore                                                                       | 1. Lårben                                                                      |
| 2. Tibia                                                             | 2. Tibia                                                                        | 2. Skenben                                                                     |
| 3. Die Zugangsstelle befindet sich 1-2 cm distal vom Mittelpunkt     | 3. Sito di accesso di 1-2 cm distale rispetto al punto centrale                 | 3. Punktionsstället ligger 1-2 cm distalt om mittpunkten                       |
| 1. Μηριαίο οστό                                                      | 1. Femur                                                                        |                                                                                |
| 2. Κνημιαίο οστό                                                     | 2. Tibia                                                                        |                                                                                |
| 3. Το σημείο προσπέλασης είναι 1-2 cm περιφερικά προς το μέσο σημείο | 3. Introductieplaats 1 tot 2 cm distaal van het middelpunt                      |                                                                                |


**Fig. 2**

- |                       |                      |
|-----------------------|----------------------|
| 1. Access Site        | 1. Sito di accesso   |
| 1. Indførsingssted    | 1. Introductieplaats |
| 1. Zugangsstelle      | 1. Tilgangssted      |
| 1. Σημείο προσπέλασης | 1. Local de acesso   |
| 1. Lugar de acceso    | 1. Punktionsställe   |
| 1. Site d'accès       |                      |

**Fig. 3**

**Fig. 4**

## INTRAOSSUEOUS INFUSION NEEDLES FOR PEDIATRIC USE (< 24 MONTHS)

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Intraosseous Infusion Needle is a needle/handle assembly. The standard design (product prefix DIN-) has a positioning mark 1 cm proximal to the needle tip and is available with a trocar, lancet, or pencil point tip. The Sussmane-Raszynski design (product prefix SDIN-) has a fine-screw needle cannula and a trocar tip.

### INTENDED USE

The Intraosseous Infusion Needle is intended for use as an alternative to intravenous access during pediatric emergencies, permitting infusion of drugs and fluids.

### CONTRAINDICATIONS

- Skin infection at insertion site
- Trauma or comminuted fracture at insertion site or compromised extremity

### WARNINGS

- Intraosseous infusion is an emergency alternative to be utilized until standard or conventional venous access can be obtained. It is **not** recommended for use longer than a 24-hour period.
- Recommended for use in children under 24 months of age. If used in pediatric patients over 24 months, the high-density hub design (C-DIN-15.5-3.0-T45-DKM-HD) should be used to provide greater stability in bones that are more dense.
- Maintain needle orientation in line with the long axis of the bone and slightly away from the joint space and epiphyseal plate.

### PRECAUTIONS

- The product is intended for use by physicians trained and experienced in emergency vascular access techniques. Standard techniques for placement of intraosseous needles should be employed.
- Possible allergic reactions should be considered.
- To avoid extravasation, advance the needle far enough into the medullary cavity for ease of infusion.

### PRODUCT RECOMMENDATIONS

#### Needle Placement

- For placement in the **proximal tibia**, palpate the tibial tuberosity and the medial border of the tibia. Locate the mid-point of the tuberosity and orient the needle 1-2 cm distal to this midpoint in the broad, flat portion of the bone. **(Fig. 1) The needle should be pointed away from the joint space and epiphyseal plate, directed distally and in line with the long axis of the bone.**
- For placement in the **distal tibia**, palpate the medial malleolus or medial surface of the distal tibia. Needle placement should be performed on the medial surface of the distal tibia, at the broad, flat junction of the shaft of

the tibia with the malleolus. (**Fig. 2**) **The needle should be pointed away from the joint space and epiphyseal plate, directed proximally and in line with the long axis of the bone.**

## INSTRUCTIONS FOR USE

### Needle Placement

1. Prepare the intended insertion site with antiseptic solution. **NOTE:** If the patient is awake or alert, administer local anesthesia.
2. Grasp the handle with the thumb and middle finger. Stabilize the needle by placing the index finger on the skin surface next to the needle tip.
3. Stabilize the extremity and begin inserting the needle with firm, downward pressure into the bone, with the needle tip directed away from the joint space and epiphyseal plate. **Always maintain needle orientation in line with the long axis of the bone. (Fig. 3)**
4. As the needle tip enters the bone, continue firm downward pressure with a steady clockwise rotation of the needle assembly. (**Fig. 4**) Continue advancing the needle in this manner until slight give is felt. The needle should stand upright, unsupported in the bone.
5. Remove the needle trocar stylet by stabilizing the baseplate of the needle cannula and turning the handle counterclockwise to disengage.
6. Confirm intramedullary needle tip position. This can be done in several ways:
  - a. Aspiration of bone marrow (presents as thick serosanguinous fluid);
  - b. Test injection of sterile saline and aspiration of bone marrow; or
  - c. Free-flow infusion through the needle with no evidence of subcutaneous extravasation. **Re-check with aspiration of bone marrow.**
7. Stabilize the needle at skin level in standard fashion, and initiate infusion.

### In the Event of Obstruction or Failure

1. Reassess the insertion site landmarks.
2. Flush the needle cannula, reinsert the stylet and repeat the insertion procedure.
3. If unsuccessful, consider alternative insertion site on the opposite leg.

### Needle Removal

- **To remove the standard intraosseous infusion needle (product prefix DIN-),** apply gentle traction and rotate the needle cannula back and forth to disengage from the bone cortex. Dress the insertion site in standard fashion.
- **To remove the Sussman-Raszynski intraosseous infusion needle (product prefix SDIN-),** apply gentle traction and turn counterclockwise to disengage from the bone cortex. Dress the insertion site in standard fashion.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2015
Bloomington, IN 47404 U.S.A.