COOK
MEDICAL

CE
0088

EN 2	Emergency Transtracheal Airway Catheter Instructions for Use
DA 3	Transtrakealt nødluftvejskateter Brugsanvisning
DE 5	Transtrachealer Katheter Gebrauchsanweisung
EL 7	Καθετήρας επείγουσας διατραχειακής αεροφόρου οδού Οδηγίες χρήσης
ES 9	Catéter transtraqueal de urgencia para la vía aérea Instrucciones de uso
FR 10	Cathéter d'intubation transtrachéale d'urgence Mode d'emploi
IT 12	Catetere transtracheale per la ventilazione d'emergenza Istruzioni per l'uso
NL 14	Acute transtracheale luchtwegkatheter Gebruiksaanwijzing
PT 16	Cateter transtraqueal de emergência para vias aéreas Instruções de utilização
SV 18	Akut transtrakeal luftvägskateter Bruksanvisning

C _ T _ D T J V B T T _ R E V 4

EMERGENCY TRANSTRACHEAL AIRWAY CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Emergency Transtracheal Airway Catheter is a 6.0 French reinforced FEP catheter, introduced over a 15 gage needle.

INTENDED USE

The Emergency Transtracheal Airway Catheter is intended for emergency airway access when conventional endotracheal intubation cannot be performed.

CONTRAINDICATIONS

None known

WARNINGS

- Firm control must be applied to the catheter to maintain position fully within the trachea at all times during use. High oxygen pressure may cause an unsecured catheter to become dislodged from the airway, resulting in subcutaneous emphysema.
- The Emergency Transtracheal Airway Catheter is not specifically designed for pediatric applications. Use with pediatric patients should be determined by attending physician.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in emergency critical care procedures. Standard techniques for placement of transtracheal airway catheters should be employed.
- Patients with airway injuries may have significant spinal injuries. Whenever possible, the cervical spine should be immobilized before beginning the procedure. Care should always be exercised to avoid additional spinal injuries.
- Whenever possible and appropriate, aseptic technique and local anesthetic should be used for the procedure.
- Do not reinsert the needle into the catheter once the needle is withdrawn. The needle may exit catheter wall.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Barotrauma
- Bleeding
- Catheter dislodgement

- Hematoma
- Pneumomediastinum
- Subcutaneous emphysema
- Tracheoesophageal fistula

INSTRUCTIONS FOR USE

1. Prep the patient's anterior neck.
2. Identify the cricothyroid membrane between the cricoid and thyroid cartilages, and stabilize the trachea with thumb and index finger.
3. Cannulate the trachea with the catheter tip facing caudad. A loss of resistance is felt when the trachea is entered.
4. Aspirate air into a water-filled 5 cc syringe to confirm catheter position within the tracheal lumen.
5. Remove the needle and advance the catheter caudad. Aspirate again to confirm position.
6. Connect the catheter to an oxygen source.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.