

- EN**  
5 **Cook Staged Extubation Set**  
Instructions for Use
- DA**  
6 **Cook trindelt ekstubationsset**  
Brugsanvisning
- DE**  
8 **Cook Set für die stufenweise Extubation**  
Gebrauchsanweisung
- EL**  
11 **ΣΕΤ σταδιακής αποδιασώληνωσης της Cook**  
Οδηγίες χρήσης
- ES**  
13 **Equipo de extubación por etapas Cook**  
Instrucciones de uso
- FR**  
15 **Set d'extubation par étapes Cook**  
Mode d'emploi
- IT**  
17 **Set di estubazione a stadi Cook**  
Istruzioni per l'uso
- NL**  
19 **Cook gefaseerde extubatieset**  
Gebruiksaanwijzing
- PT**  
21 **Conjunto de extubação faseada Cook**  
Instruções de utilização
- SV**  
23 **Cook set för etappvis extubation**  
Bruksanvisning


**Fig. 1**

**Fig. 2**

**Fig. 3**


Fig. 7


Fig. 8


Fig. 9

## COOK STAGED EXTUBATION SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Cook Staged Extubation Set is a tapered Staged Re-Intubation Catheter, with centimeter markings and a Nitinol Staged Extubation Wire with depth markings. The set is supplied with two Rapi-Fit Adapters (15 mm connector or a Luer lock connector), an extra-firm, soft-tipped, Staged Re-Intubation Catheter, a floppy-tipped Staged Extubation Wire, securing tape, airway wire holder, and securement dressing.

### INTENDED USE

The Cook Staged Extubation Set is intended for staged extubation and subsequent re-intubation.

### CONTRAINDICATIONS

None known

### WARNINGS

- Extubation should be executed only according to standard hospital protocol.
- Strict attention should be paid to the insertion depth of the Staged Extubation Wire by using the featured depth markings.
- Ensure the Staged Extubation Wire is secured and placement is maintained to prevent distal migration or proximal displacement from the glottis.
- Depth markings on the Staged Extubation Wire refer to the distance from the distal tip.
- Attention should be paid to the insertion depth of the Staged Re-Intubation Catheter in the patient's airway and correct tracheal position of replacement endotracheal tube.
- Centimeter markings on the Staged Re-Intubation Catheter refer to the distance from the distal tip.
- The Staged Extubation Wire and Staged Re-Intubation Catheter should be advanced to the optimal position according to physician's preference.
- Use of the Rapi-Fit adapter for oxygenation may be associated with a risk of barotrauma.
- Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume.
- To avoid barotrauma and/or pneumothorax, examine the patient's anatomy to help determine the optimal placement for the Staged Re-Intubation Catheter. Ensure the catheter is in the physician-preferred location relative to the carina by referencing its centimeter markings.

### PRECAUTIONS

- This product is intended for use by physicians trained and experienced in airway anatomy and airway management. Standard techniques for staged extubation, re-intubation and oxygenation should be employed.
- This product should not be used with endotracheal tubes with an inner diameter smaller than 5 mm.
- This device is not recommended for definitive, long-term airway management. Intubation, tracheostomy or the return of spontaneous ventilation is recommended before removing this device.
- Care should be taken to maintain the Staged Extubation Wire's position during placement of the Staged Re-Intubation Catheter.
- This product is not intended for intravascular use.

### POTENTIAL ADVERSE EVENTS

- Barotrauma
- Pneumothorax
- Perforation of the bronchi or lung parenchyma
- Staged Extubation Wire migration or displacement
- Staged Re-Intubation Catheter migration

### INSTRUCTIONS FOR USE

#### Staged Extubation

1. Confirm correct endotracheal tube position before advancing the Staged Extubation Wire into the endotracheal tube.
2. Using the patient's teeth as a landmark, note the depth of the endotracheal tube. **NOTE:** This is the recommended depth to which the Staged Re-Intubation Catheter will need to be inserted during re-intubation.
3. Administer 100% oxygen.
4. Suction the oropharynx.
5. Advance the Staged Extubation Wire, floppy end first, into the endotracheal tube to the predetermined depth. (**Figs. 1 and 2**) The optimal insertion depth should be determined by the physician, taking into consideration the patient's anatomy. **NOTE:** The Staged Extubation Wire has 5 depth markings, indicating distance from the distal tip. The

double band is positioned 20 cm from the distal tip, a single band at 25 cm, a triple band at 30 cm, a single band at 35 cm and a quadruple band at 40 cm. (Fig. 3)

6. Fully deflate the cuff of the endotracheal tube.
7. While maintaining position of the Staged Extubation Wire, extubate the patient. (Figs. 4 and 5)
8. Introduce the proximal end of the Staged Extubation Wire into the smaller wire holder provided in the set. Secure the wire holder to the patient with the securement device provided. **NOTE:** Proper securement will minimize the potential for accidental dislodgement.
9. Remove and discard the Staged Extubation Wire according to hospital protocol when it is clinically determined the device is no longer necessary for airway access. **NOTE:** Usage of the Staged Extubation Wire should not exceed 120 hours. This is a one-time use device.

## Re-Intubation

**Re-Intubation should be executed only according to standard hospital protocol.**

1. Administer 100% Oxygen (if possible).
2. Carefully remove the Staged Extubation Wire from the wire holder while maintaining its position within the patient's airway. Re-confirm position relative to the patient's teeth by referencing depth marks on the Staged Extubation Wire; ensure the wire is at the pre-determined depth. (Fig. 6)
3. Advance the Staged Re-Intubation Catheter over the Staged Extubation Wire until the desired depth is reached. (Fig. 7) **NOTE:** Consideration should be given to the depth to which the original endotracheal tube was inserted. **NOTE:** If oxygen insufflation is needed, the Staged Extubation Wire must be removed before the Rapi-Fit Adapter can be attached to the catheter (see below for instructions). Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume. If an oxygen source is used for insufflation, begin at a lower pressure and work up gradually. Observe the chest for outward and inward movements to confirm oxygen insufflation and egression. Pulse oximetry and oral air flow should be carefully monitored as well. In cases of upper airway obstruction, gas discharge from the patient's lungs may require more time.
4. With the catheter in position, advance the endotracheal tube over the catheter and into position within the patient's trachea. (Fig. 8) **NOTE:** Always maintain control and position of the Staged Extubation Wire and Re-Intubation Catheter during advancement of the endotracheal tube.
5. Advance the endotracheal tube to the optimal position and inflate the balloon cuff (if provided).
6. Remove the Staged Extubation Wire and Staged Re-Intubation Catheter from the endotracheal tube. (Fig. 9)
7. Verify position and secure the endotracheal tube in standard fashion.

## Use of Rapi-Fit Adapter for Oxygenation

Rapi-Fit adapters should only be used when oxygen requirements are high and intubation is unsuccessful. Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume. If an oxygen source is used for insufflation, begin at a lower pressure and work up gradually. Observe the chest for outward and inward movements to confirm oxygen insufflation and egression. Pulse oximetry and oral air flow should be carefully monitored as well. In cases of upper airway obstruction, gas discharge from the patient's lungs may require more time. **NOTE:** If oxygen insufflation is needed, the Staged Extubation Wire must be removed before the Rapi-Fit Adapter can be attached to the Staged Re-Intubation Catheter.

1. To attach the Rapi-Fit Adapter, position the adapter on the catheter, then push the white collar forward and lock into position.
2. To remove the adapter, pull the white collar back to release, and then remove from the catheter.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd