Frova Intubating Introducers	Instructions for Use	EN 3
Intubační zavaděče Frova	Návod k použití	CS 4
Frova intubationsindførere	Brugsanvisning	DA 5
Frova-Intubationshilfen	Gebrauchsanweisung	DE 6
Εισαγωγείς διασωλήνωσης Frova	Οδηγίες χρήσης	EL 7
Introducciones para intubación Frova	Instrucciones de uso	ES 8
Introduceurs Frova pour intubation	Mode d'emploi	FR 10
Introduttori per intubazione Frova	Istruzioni per l'uso	IT 11
Frova intubatie-introducers	Gebruiksaanwijzing	NL 12
Introdutores de intubação Frova	Instruções de utilização	PT 13
Frova intubationsledare	Bruksanvisning	SV 14

CE
0123

COOK[®]
MEDICAL

C _ T _ F I I _ R E V 6

FROVA INTUBATING INTRODUCERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Frova Intubating Introducer is a radiopaque catheter introducer with centimeter markings and a blunt, curved tip that can be passed blindly into the trachea. Products with the suffix "-FII" in their order numbers also include a stiffening cannula and two Rapi-Fit® adapters for connection to a ventilatory device. Products with the suffix "-FI" include two Rapi-Fit® adapters for connection to a ventilatory device, but no stiffening cannula.

INTENDED USE

The Frova Intubating Introducer is intended to facilitate endotracheal intubation in patients where visualization of the glottis is inadequate. The product has been designed for use by physicians trained and experienced in airway management techniques. Standard techniques for placement of intubation catheters should be employed.

CONTRAINDICATIONS

None known

WARNINGS

- Attention should be paid to insertion depth of catheter into the patient's airway and correct tracheal position of endotracheal tube. Markers on the Frova Intubating Introducer refer to distance from tip of catheter. Catheter and endotracheal tube should not be advanced beyond the carina.
- Ensure that the Rapi-Fit Adapter is securely connected to the Frova Intubating Introducer prior to oxygen delivery. Failure to properly secure the adapter to the introducer may result in hypoxia, hypoxemia, and serious adverse events.
- Oxygen delivery through the Frova Intubating Introducer is not possible without a Rapi-Fit Adapter. Rapi-Fit Adapters are only included in products with an "-FI" or "-FII" suffix in their product numbers.

PRECAUTIONS

Sizing of the Rapi-Fit Adapter and the Frova Intubating Introducer must match in order to provide oxygenation through the catheter.

POTENTIAL ADVERSE EVENTS

- Injury to the epiglottis
- Perforation of the piriform fossa, trachea or bronchus
- Barotrauma

PRODUCT RECOMMENDATIONS

- The 14.0 French Intubating Introducer is designed for placement of an endotracheal tube with an inner diameter of 6 mm or larger.
- The 8.0 French Intubating Introducer is designed for placement of an endotracheal tube with an inner diameter of 3 mm or larger.

INSTRUCTIONS FOR USE

1. Preload an appropriately sized, lubricated endotracheal tube onto the proximal part of the introducer.
2. Using laryngoscopy, introduce the tip of the introducer beyond the epiglottis and advance it in a straight line toward the glottis.
3. Advance the introducer into the trachea. If resistance is encountered, **do not force the introducer**; gently rotate and advance. After 2-3 cm of introduction, remove the stiffening cannula (if included) and gently advance the introducer into the trachea. **NOTE:** During introduction, you may feel tracheal rings.
4. Confirm introducer position using standard methods (e.g., capnography, breath sounds, X-ray).
5. Maintaining position of the introducer, advance the preloaded endotracheal tube into the trachea to the appropriate distance.
6. Maintaining position of the endotracheal tube, remove the introducer and laryngoscope.
7. Confirm endotracheal tube position using standard methods (e.g., capnography, breath sounds, X-ray).

NOTE: Products with the suffix "-FII" or "-FI" in their order numbers include two Rapi-Fit adapters for use with a ventilatory device.

- To attach the adapter, position it on the catheter, push the white collar forward and lock into position.
- To remove the adapter, pull the white collar back to release, and remove it from the catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2020
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd