COOK

MEDICAL

CE
0088

EN

3

Silicone Central Venous Catheters

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Bruksanvisning

C _ T _ H A S T _ R E V 4

SILICONE CENTRAL VENOUS CATHETERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Silicone Central Venous Catheters are available in single-, double-, and triple-lumen configurations.

INTENDED USE

Silicone Central Venous Catheters are intended for vascular access infusion and withdrawal of blood, blood products and fluids; central venous blood pressure monitoring (CVP); acute hyperalimentation; and continuous or intermittent drug infusion.

CONTRAINDICATIONS

None known

WARNINGS

- Complications arising from the use of this device can result in serious injury or death; the catheter tip can erode or perforate vascular walls. Only medical practitioners licensed by law, trained, and experienced in proper positioning of catheters in the central venous system using percutaneous entry (Seldinger) technique should place this catheter. **Extreme caution must be used in placement and monitoring.**
- Placement of this catheter into the subclavian vein may result in compression of the catheter by the clavicle and first rib. Excessive compression may result in catheter damage, including rupture or catheter embolus. Compression of the catheter and consequent risks may be minimized by placement lateral to the junction of the clavicle and first rib or use of alternate venous sites.
- Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system. Tip position should be verified by X-ray and monitored on a routine basis. Periodic lateral view X-ray is suggested to assess tip location in relation to vessel wall. Tip position should appear to be parallel to vessel wall.
- To avoid vascular injury, do not use excessive force when advancing dilators. Use the smallest size dilator catheter placement will allow. The wire guide must always lead the dilator by several centimeters. Do not advance the dilator more than a few centimeters into the vessel.
- **Do not power inject contrast medium through the catheter.** Catheter rupture may result. Use of a 10 ml or larger syringe will reduce the risk of catheter rupture.
- To distend great vessels and to prevent inadvertent air aspiration during catheter insertion, the patient should be placed in Trendelenburg position.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in the proper positioning of catheters in the central venous system using

percutaneous entry (Seldinger) technique. Standard techniques for placement of central venous catheters should be employed.

- Select the puncture site and the length of catheter needed by assessing patient anatomy and condition.
- Patient movement can cause catheter tip displacement. Catheters placed via an antecubital vein have shown tip movement of up to 10 cm with motion of the extremity. Catheters placed via either a jugular or subclavian vein have demonstrated tip movement of 1-3 cm with neck and shoulder motion.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify the attending physician immediately.
- Use of ECG, ultrasound, and/or fluoroscopy is suggested for accurate catheter placement.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.


PRODUCT RECOMMENDATIONS

Suggested Catheter Maintenance

- Catheter entry site must be prepared and maintained in a manner consistent with standard procedure for central venous catheterization.
- To prevent clotting or possibility of air embolus, the double-lumen's #2 lumen, and the triple-lumen's #2 and #3 lumens should be filled with heparinized saline solution (100 units of heparin per ml of saline is usually adequate) prior to catheter introduction.
- After catheter is placed and prior to use, tip position and lumen patency should be confirmed by free aspiration of venous blood. **If blood is not freely aspirated, physician should immediately reevaluate catheter tip position.**
- Any unused lumens should be maintained with continuous saline or heparinized saline drip or locked with heparinized saline solution.
- Before using any lumen already locked with heparin, lumen should be flushed with twice the indicated lumen volume using normal saline. Lumens should be flushed with normal saline between administrations of different infusates. After use, lumen should again be flushed with twice the indicated lumen volume using normal saline before reestablishing heparin lock.
- Strict aseptic technique must be adhered to while using and maintaining catheter.

Access Sites of Choice

1. High Internal Jugular
2. External Jugular
3. Low Internal Jugular
4. Supraclavicular
5. Infraclavicular
6. Femoral (not shown)

INSTRUCTIONS FOR USE

1. Measure the catheter to be used against the patient to determine the approximate length of catheter needed from the puncture site to central venous tip position.
2. Introduce the access needle into the vessel. Venous blood should be easily aspirated to confirm position of the needle tip within the vessel.
3. Slide the Safe-T-J® wire guide straightener (positioned on the distal tip of the wire guide) over the "J" portion of the wire guide.
4. Pass the straightened wire guide through the needle and advance it 5-10 cm into the vessel.

NOTE: If a straight wire guide is used, always advance the soft, flexible end through the needle hub and into the vessel. **If resistance is encountered during wire guide insertion, do not force the wire guide. Withdrawal of the wire guide through the needle should be avoided; breakage may result.**

5. While maintaining wire guide position, withdraw the needle and Safe-T-J wire guide straightener.

6. Enlarge the puncture site with a scalpel if necessary. If dilation is required, the dilator can be advanced over the wire guide and removed prior to inserting the central venous catheter.

CAUTION: To avoid vascular injury, do not use excessive force when advancing dilators. Use the smallest size dilator catheter placement will allow. The wire guide must always lead the dilator by several centimeters. Do not advance the dilator more than a few centimeters into the vessel.

7. If a subcutaneous tunnel is to be used, the catheter can be passed through the tunnel at this time.
8. Introduce the sheath/dilator assembly over the wire guide and advance the assembly into the vessel with a twisting motion.
9. Leaving the sheath in place, remove the dilator and wire guide. To prevent inadvertent air aspiration, place the thumb or finger over the cuffed, proximal end of the sheath after removing the dilator and wire guide.
10. Introduce the catheter into the sheath and advance the catheter into position. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates. **NOTE:** After the catheter is in position, venous blood should be easily aspirated.
11. Grasp the two knobs of the sheath and pull outward and upward simultaneously to peel away and remove the sheath from the catheter.
12. Secure the catheter in place with suture or other method. If the catheter is not introduced to its full length, additional suture should be carefully placed around the catheter and affixed to the skin at the entry site to help prevent catheter movement.
13. The catheter lumens should be flushed with 5-10 ml of normal saline prior to use or establishment of heparin lock.

NOTE: A wire guide that is at least twice as long as the catheter is recommended for catheter exchange procedures.

NOTE: The 12 French triple-lumen catheter is designed for high-volume infusion and withdrawal through the two main lumens. The smaller third lumen is ideally suited for administration of therapeutic drugs or pressure monitoring.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.