ENGLISH

Cook Pressure Monitoring Catheter

Intended Use

The Cook Pressure Monitoring Catheter is used for arterial pressure monitoring and blood sampling. It is provided sterile in peel-open packages and is intended for one-time use.

Warnings

- Complications arising from the use of this device can result in serious injury or death; catheter tip can erode or perforate vascular walls.
- Only medical practitioners licensed by law, trained and experienced in proper positioning of catheters in the vascular system using percutaneous entry (Seldinger) technique should place this catheter. **EXTREME CAUTION** must be used in placement and monitoring.
- Select puncture site and length of catheter needed by assessing patient anatomy and condition.
- Catheter size should be as small as the use will allow.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify attending physician immediately.
- Catheter should not be used for long-term indwelling applications.
- Patient movement can cause catheter tip displacement. Use should be limited to controlled hospital situation.
- Do not cut, trim or modify catheter or components, prior to placement or intraoperatively.
- To avoid vascular injury, do not use excessive force when advancing dilators. Use the smallest size dilator catheter placement will allow. Wire guide must always lead dilator by several centimeters. Do not advance dilator more than a few centimeters into the vessel.

Suggested Instructions for use of the Cook Pressure Monitoring Catheter

- Introduce thinwall percutaneous entry needle into vessel. Arterial blood should be easily aspirated to confirm position of needle tip within vessel. (Figure 1)
- Slide Safe T-J[®] wire guide straightener (positioned on distal tip of wire guide) over “J” position of wire guide. Pass straightened wire guide through needle; advance wire guide 5-10 cm into vessel. If straight wire is used, always advance soft, flexible end through needle hub and into vessel. If resistance is encountered during wire guide insertion, do not force wire guide. Withdrawal of wire guide through needle should be avoided; breakage may result. (Figure 2)
- While maintaining wire guide position, withdraw needle and Safe T-J[®] wire guide straightener. (Figure 3)
- Enlarge puncture site with number 11 scalpel blade, if required. If dilation is required, dilator can be advanced over wire guide and removed prior to insertion of arterial catheter. (Figure 4)

Caution: To avoid vascular injury, do not use excessive force when advancing dilator. Use the smallest size dilator catheter placement will allow. Wire guide must always lead dilator by several centimeters. Do not advance dilator more than a few centimeters into the vessel.
- Introduce Single Lumen Arterial Catheter over wire guide. While maintaining wire guide position, advance catheter into vessel with a gentle twisting motion. Do not advance catheter tip beyond distal tip of wire guide. Always have wire guide leading during catheter placement. (Figure 5)
- After catheter is in position, remove wire guide along with inner dilator. Arterial blood should be easily aspirated. Winged hub can now be sutured into place. If catheter is not introduced to its full length, additional suture should be carefully placed around catheter and affixed to the skin at entry site. This will help prevent backward or forward catheter movement. Lumen should now be flushed with 5-10 cc normal saline prior to use or establish a heparin lock. (Figure 6)

Note: A wire guide that is at least twice as long as the catheter is recommended for catheter exchange procedure.

Packaging/Sterilization

Supplied in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Store in cool, dry place.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

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