COOK

MEDICAL

CE
0088

EN 3	Left Atrial Pressure Monitoring Set Instructions for Use
DA 5	Sæt til monitorering af tryk i venstre atrium Brugsanvisning
DE 7	Set zur Überwachung des linken Vorhofdrucks Gebrauchsanweisung
EL 9	Σετ παρακολούθησης πίεσης αριστερού κόλπου Οδηγίες χρήσης
ES 11	Equipo de monitorización de la presión auricular izquierda Instrucciones de uso
FR 13	Set de contrôle de la pression auriculaire gauche Mode d'emploi
IT 15	Set per il monitoraggio della pressione atriale sinistra Istruzioni per l'uso
NL 17	Set voor drukbewaking van het linker atrium Gebruiksaanwijzing
PT 19	Conjunto de monitorização da pressão auricular esquerda Instruções de utilização
SV 21	Set för tryckövervakning i vänster förmak Bruksanvisning

C _ T _ L A P _ R E V 3

LEFT ATRIAL PRESSURE MONITORING SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Left Atrial Pressure Monitoring Set consists of a polyethylene or polyurethane pigtail catheter, a wire guide introducer needle, a curved catheter introducer needle, and a wire guide.

INTENDED USE

The Left Atrial Pressure Monitoring Set is intended for postoperative monitoring of left atrial pressure.

CONTRAINDICATIONS

None known

WARNINGS

- Complications arising from the use of this device can result in serious injury or death.
- Every effort must be made to properly place and monitor catheters in order to prevent erosion or perforation of the vascular walls.
- Tip position should be verified by X-ray and monitored on a routine basis.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in atrial pressure monitoring techniques. Standard techniques should be employed.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify attending physician immediately.
- The catheter should not be used for long-term indwelling applications.
- Patient movement can cause catheter tip displacement. Use should be limited to controlled hospital situations.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Air embolus
- Thrombosis
- Bleeding
- Migration out of the left atrium
- Dislodgement

INSTRUCTIONS FOR USE

1. Following sternotomy, place a circumferential (purse string) suture with a pledget at the desired site of access to the left atrium. Typically, this is the right superior pulmonary vein.

2. With ventilation paused, insert the straight needle tip a short way (usually less than 1cm) into the insertion site.
 3. Pass the wire guide through the needle hub and advance the wire guide into the atrium.
 4. Withdraw the needle, leaving the wire guide in place.
 5. Using the curved catheter introducer needle, make a skin puncture at the desired exit site for the atrial catheter. Typically, this is below the right subcostal margin.
 6. Carefully advance the needle into the mediastinum, taking care to avoid any vital structures.
 7. Insert the proximal end of the wire guide into the tip of the curved needle and advance the wire guide until it exits the needle hub.
 8. Withdraw the needle, leaving the wire guide in place.
 9. Introduce the catheter over the wire guide and into the left atrium.
 10. Remove wire guide.
 11. Attach a transducer pressure measurement system. Ensure that all air is removed from the catheter and the remainder of the pressure transducing fluid path. Aspirate from the catheter, checking for free flow of blood.
 12. Transduce the pressure and examine the pressure waveform to ensure that the catheter tip is in the left atrium.
- NOTE:** Frequently during catheter placement, the initial catheter tip position is in the left ventricle. It may be necessary to reposition the catheter tip to the left atrium.
13. Verify that the catheter can move without undue friction.
 14. After the chest is closed, secure the catheter and ensure that no air is in the fluid path.

NOTE: The catheter should be easily removed with a gentle tug. If the catheter cannot be removed in this manner, it may be excessively restrained by the circumferential suture. If this is the case, mediastinal exploration may be necessary for catheter removal. When the catheter is removed, examine the tip to ensure it is intact.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2011
Bloomington, IN 47404 U.S.A.