EN 3	Pneumothorax Sets Instructions for Use
DA 4	Pneumothorax-sæt Brugsanvisning
DE 5	Pneumothorax-Sets Gebrauchsanweisung
EL 7	ΣΕΤ ΠΝΕΥΜΟΘΩΡΑΚΑ Οδηγίες χρήσης
ES 8	Equipos para neumotórax Instrucciones de uso
FR 10	Sets pour pneumothorax Mode d'emploi
IT 12	Set per pneumotorace Istruzioni per l'uso
NL 13	Pneumothoraxsets Gebruiksaanwijzing
PT 15	Conjuntos de pneumotórax Instruções de utilização
SV 16	Pneumotoraxset Bruksanvisning


Fig. 1

Fig. 3

Fig. 2

Fig. 4

PNEUMOTHORAX SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Pneumothorax Sets consist of:

- 1 stainless steel needle
- 1 radiopaque catheter
- 1 Chest Drain Valve
- 1 polyvinylchloride connecting tube
- 1 plastic one-way stopcock

INTENDED USE

Pneumothorax Sets are intended for relief of simple spontaneous or iatrogenic pneumothorax and tension pneumothorax.

CONTRAINDICATIONS

Pneumothorax Sets are not recommended for large fluid accumulation or hemothorax.

WARNINGS

None known

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in emergency treatment of pneumothorax. Standard techniques for placement of pneumothorax catheters in the pleural cavity should be employed.
- Lung puncture may result in an air embolus, which could lead to ischemia or infarction of major organs including the brain or cardiac system.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Preassemble catheter, connecting tube with stopcock and Chest Drain Valve. **NOTE:** Chest Drain Valve may be obviated if catheter is to be connected to a water seal suction apparatus or similar mechanical suction device. Do not connect catheter directly to wall suction. Attach Chest Drain Valve in direction indicated by arrow on valve.
2. Prep the access site with antiseptic and drape in a standard fashion. **NOTE:** The optimal insertion site of the pneumothorax catheter is the second anterior intercostal space either to the right or left of the sternum at the level of the angle of Louis, in the mid-clavicular line.
3. After the introduction (through skin, subcutaneous tissue, muscle and fascia down to the parietal pleura) of local anesthesia, a 2 mm incision is made through the skin, subcutaneous tissue, muscle and fascia, but not through the parietal pleura.
4. Insert the needle tip, with the catheter as its sheath, through the incision site into the pleural cavity at the level of the second anterior intercostal space just above the upper border of the third rib. Advance the needle and catheter with one hand while the other hand maintains pressure upon the stylet within the needle cannula.
5. If a large pneumothorax is present (as determined by fluoroscopy or by roentgenograms) and there is certainty that the catheter and cannula are well within the pleural cavity, immediately push the catheter forward by sliding it over the needle cannula and into the pleural space. **NOTE:** If there is uncertainty as to whether the needle and catheter tip are well within the pleural cavity, as with a small pneumothorax, the inner needle stylet may be removed to determine whether there is a free flow of air. Once positioned within the pleural space, advance the catheter so that all catheter sideports are within the pleural space. As the catheter is advanced, the needle is slowly withdrawn.
6. Monitor the catheter tip position fluoroscopically. The distal tip of the catheter should be positioned with the tip pointing toward the extreme apex of the pleural space.
7. Stabilize the catheter near the locking sleeve with one hand, grasp the hub area of the locking sleeve and pull back. (**Fig. 1**)
8. Pull the sleeve until it snaps over the shoulder of the plastic fitting on the catheter shaft. (**Fig. 2**) **NOTE:** If the locking sleeve is not fully snapped into position, the loop will not be properly secured, and leaking will occur at the locking sleeve.
9. Attach the catheter to the chest with adhesive tape and skin suture.
10. **ALL CONNECTIONS MUST BE SECURE AND AIRTIGHT.** Perform inspections of the catheter and connections regularly. If the Chest Drain Valve is removed, the stopcock should be turned to the off position to prevent re-entry of air into the pleural space.
11. Obtain a chest roentgenogram within 24 to 48 hours before catheter removal.

Catheter Removal

1. Partially insert the appropriate stiffening cannula into catheter for support. (**Fig. 3**)
2. Grasp the locking sleeve just beyond the hub and advance forward over the shoulder of the plastic fitting and proximal catheter shaft (**Fig. 4**) until it stops. This will release the loop. **NOTE:** Advance the stiffening cannula into the catheter so its distal tip is only a few centimeters beyond the locking sleeve portion of the catheter.
3. To remove catheter, have patient hold his/her breath in full inspiration. Withdraw the catheter and place petrolatum-impregnated gauze dressing over the chest wall opening; dress the gauze with adhesive tape.
NOTE: For catheter replacement, total removal of the Simp-Loc mechanism is recommended. With Simp-Loc in locked position, cut completely through the locking sleeve near the hub, advance the remaining locking sleeve forward over the plastic fitting to expose suture ends, and pull one suture end completely out of the catheter. If resistance is encountered, a slight tugging motion may be required to free the suture from the catheter. Standard wire guide exchange procedure for catheter replacement can now be performed.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.