**COOK**

**MEDICAL**

**CE**  
0088

**EN**  
**3**

**Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube**

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Bruksanvisning

C \_ T \_ N J F T U \_ R E V 8

## TIGER 2 SELF-ADVANCING NASAL JEJUNAL FEEDING TUBE

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube is a 14 French, 155 cm long polyurethane feeding tube featuring innovative flaps that facilitate placement into the jejunum by allowing the natural movement of the enteral system to gently drag the catheter into the small bowel. The catheter features centimeter markings and has five sideports to optimize administration of medications and nutrition.

### INTENDED USE

The Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube is intended to provide short-term enteral access for delivery of nutrition and/or medications to the small bowel.

### CONTRAINDICATIONS

Use of the Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube is contraindicated in patients with the following conditions:

- Esophageal varices
- Gastric varices
- INR (international normalized ratio) > 1.3 (at time of insertion and/or expected at time of removal)
- Anticoagulated patients (anticoagulated at time of insertion and/or expected to be anticoagulated at time of removal)
- Pathologic coagulopathies
- History of bleeding disorder(s)
- Small or large bowel obstruction(s)
- Ischemic bowel
- Peritonitis
- Esophageal stricture or obstruction
- Gastric obstruction
- Recent nasal, oral, esophageal, or gastric surgery or trauma
- Deviated septum
- Inability to pass the feeding tube through the nares
- Uncooperative patient

### WARNINGS

- The Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube is not intended for intravascular use.
- Improper placement of the feeding tube into the trachea or lungs can lead to serious injury. Always check the final X-ray prior to tube feeding.

### PRECAUTIONS

- A thorough understanding of the technical principles, clinical applications and risks associated with nasal jejunal feeding tube placement is necessary before using the device. Tiger 2 Self-Advancing Nasal Jejunal Feeding Tubes should only be used by or under the supervision of personnel trained in standard nasal gastric tube placement procedure.
- When removing the Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube through the nose or mouth, proceed very slowly.
- The feeding adapter is manufactured from PVC (polyvinylchloride), which contains DEHP (di(2-ethylhexyl)phthalate). The catheter itself does **not** contain DEHP.
- The Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube is not intended for use beyond 30 days.
- The torque cable (supplied separately) should only be used to facilitate placement of the Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube through the nose or mouth and into the stomach. The Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube should only be advanced 50-70 cm into the stomach.
- Do not insert the torque cable into a Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube that is already *in situ*. Doing so may cause damage to internal tissues/structures.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

### POTENTIAL ADVERSE EVENTS

Potential adverse events associated with placement and use of a Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube include, but are not limited to:

- Bleeding
- Clogged or leaking feeding tube
- Sinusitis
- Premature displacement of the feeding tube
- Aspiration
- Nasal irritation
- Manual removal of the feeding tube by the patient
- Sore throat
- Bowel erosion and/or intestinal perforation

- Pneumothorax from inadvertent pleural placement
- Pulmonary injury
- Improper placement of the feeding tube into the lung or brain
- Coughing or dyspnea, indicating potential bronchial placement
- Feeding tube coiled in esophagus or posterior pharynx
- Esophageal tear from trauma of tube passing
- Feeding tube dislodgement during therapy, necessitating removal and new feeding tube placement
- Aspiration of stomach contents despite appropriate placement

## INSTRUCTIONS FOR USE

1. Upon removal of the device from the package, visually inspect with particular attention to kinks, bends or breaks. If an abnormality is detected that would prohibit proper working condition or cause an adverse patient outcome, do not use the product.
2. Measure from the patient's nose to his or her sub-xyphoid before tube placement.
3. Apply lubricant to the distal tip of the feeding tube.
4. Advance the lubricated feeding tube through the patient's nose or mouth 50-70 cm into the stomach. This distance is dependent on the patient's anatomical measurement. The operator may ask the patient to swallow as the tube is passed. Stop and confirm that the tube is in the esophagus/stomach and not in the mainstem bronchus.

**NOTE: If the patient is awake, the operator may give him or her water through a straw to help ease tube passage during insertion. Discomfort or gagging may be expected. If the patient begins coughing during this step, the Tiger 2 may be incorrectly placed in the airway. If this occurs, the procedure should be immediately stopped and the feeding tube should be removed.**

**NOTE:** The use of a torque cable (TSBD-65-170-NJFT) will provide additional support inside the Tiger 2 for advancement in difficult cases.

5. The use of insufflation and auscultation may confirm that the catheter is in the stomach. However, do not use air auscultation as definitive confirmation of final tube placement.
  6. The feeding tube should be left in place for 30 minutes to 1 hour to allow the patient's peristalsis to advance the catheter.
  7. Manually advance the feeding tube 10 cm every 30 minutes to 1 hour until the 100 cm mark is reached. If peristalsis is very weak, the feeding tube can be advanced 10 cm every two hours. If the patient's stomach is anatomically unusual, advance the tube in 5 cm increments.
- NOTE:** Pharmacologic agents may be used to increase peristalsis in accordance with institutional protocol.
8. Obtain a final confirmation X-ray that visualizes the entire course of the inserted tube prior to using the tube for the initial feeding or medication administration. This will help ensure that the tube is in the desired position (either the stomach or small bowel).
  9. Use hypoallergenic tape to secure the feeding tube in place.

## Tube Removal

**CAUTION: When removing the Tiger 2 Self-Advancing Nasal Jejunal Feeding Tube through the nose or mouth, proceed *very slowly*.**

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.