EN 4	Peripherally Inserted Central Venous Catheter Sets and Trays With Peel-Away® Introducers Instructions For Use
DA 10	Centralvenekatetersæt og bakker til perifer indføring med Peel-Away®-indførere Brugsanvisning
DE 15	Peripher eingeführte zentrale Venenkatheter-Sets mit Peel-Away® Einführschleuse Gebrauchsanweisung
EL 20	ΣΕΤ και δίσκοι κεντρικών φλεβικών καθετήρων που εισάγονται περιφερικά με εισαγωγείς Peel-Away® Οδηγίες χρήσης
ES 26	Equipos y bandejas de catéteres venosos centrales de inserción periférica con introductores Peel-Away® Instrucciones de uso
FR 31	Sets et plateaux de cathéter veineux central inséré en périphérie avec introducteurs Peel-Away® Mode d'emploi
IT 37	Set standard e set completi con catetere venoso centrale ad inserimento periferico con introduttore Peel-Away® Istruzioni per l'uso
NL 43	Perifeer ingebrachte centraalveneuze kathetersets en -pakketten met Peel-Away® introducers Gebruiksaanwijzing
PT 48	Conjuntos e tabuleiros e de cateteres venosos centrais de introdução periférica com introdutores Peel-Away® Instruções de utilização
SV 54	Perifert införda centrala venkateterset och -brickor med Peel-Away®-införare Bruksanvisning


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6


Fig. 7

PERIPHERALLY INSERTED CENTRAL VENOUS CATHETERS

WITH PEEL-AWAY® INTRODUCER

DEVICE DESCRIPTION

Silicone Peripherally Inserted Central Venous Catheters feature a soft silicone suture wing for secure fixation to the skin. Reinforced silicone extensions reduce catheter damage during catheter manipulations. Catheters can be trimmed to fit the patient's anatomy. Optional plastic clamp(s) help prevent air aspiration due to inadvertent hub dislodgment.

Single-Lumen Catheters

French Size	Lumen Volume	Lumen Diameter	Length
3.0	0.5 cc	.022 inch	50 cm
4.0	0.6 cc	.026 inch	60 cm
5.0	0.8 cc	.036 inch	60 cm

Double-Lumen Catheters

French Size	Small Lumen Diameter	Small Lumen Volume	Large Lumen Diameter	Large Lumen Volume	Length
6.0	.018 inch	0.4 cc	.026 inch	0.5 cc	49 cm
7.0 ¹	.019/.045 inch	0.7 cc	.025/.050 inch	0.8 cc	50 cm

¹D-Shaped Lumens

INTENDED USE

Silicone Peripherally Inserted Central Venous Catheters (PICC) have been designed for the intravenous administration of nutrient fluids, chemotherapeutic agents and other drugs for therapy.

Peripherally Inserted Central Venous Catheters are suggested for:

1. Delivery of whole blood or blood products
2. Drug administration
3. Blood sampling

WARNINGS

- Only medical practitioners licensed by law, trained and experienced in proper positioning of catheters in the central venous system using percutaneous entry (Seldinger) technique should place this catheter. **EXTREME CAUTION** must be used in placement and monitoring.

- Silicone Peripherally Inserted Central Venous Catheters (PICC) are not designed for power injection of contrast material. **Do not power inject contrast medium through catheter.** Catheter rupture may result. Use of 10 cc syringe or larger will reduce the risk of catheter rupture.
- **Peripherally Inserted Central Venous Catheters play an important role in treatment of critically ill patients; however, catheter tips can erode or perforate vascular walls. Every effort must be made to properly place and monitor catheters.**
- Catheter tip position should be verified by x-ray and monitored on a routine basis. Periodic lateral view x-ray is suggested to assess tip location in relation to vessel wall. Tip position should appear to be parallel to vessel wall.

PRECAUTIONS

- Select puncture site and length of catheter needed by assessing patient anatomy and condition.
- Catheter size should be as small as the use will allow.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify attending physician immediately.
- Patient movement can cause catheter tip displacement. Catheters placed via an antecubital vein have shown forward tip movement of up to 10 cm with motion of the extremity.
- Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

The potential exists for serious complications including the following:

- Air embolism
- Arterial puncture
- Bleeding
- Brachial plexus injury
- Cardiac arrhythmia
- Cardiac tamponade
- Catheter damage
- Catheter erosion through the skin
- Catheter embolism
- Catheter occlusion
- Catheter related sepsis
- Endocarditis
- Exit site infection
- Exit site necrosis
- Extravasation
- Fibrin sheath formation

- Failure to advance
- Hematoma
- Insertion site drainage
- Intolerance reaction to implanted device
- Laceration of vessels or viscus
- Myocardial erosion
- Nerve injury
- Perforation of vessels or viscus
- Phlebitis
- Spontaneous catheter tip malposition or retraction
- "Stuck" catheter
- Thromboembolism
- Venous thrombosis
- Vessel erosion
- Risks normally associated with local or general anesthesia, surgery, and post-operative recovery

PRODUCT RECOMMENDATIONS

Catheter Size And Puncture Site


Preliminary reports indicate that catheter size can influence clotting; larger diameter catheters have more tendency to promote clots. As reported by Amplatz, Gianturco and others, clot formation has less relation to type of catheter material than to size of catheter.

The angle of the catheter tip to the vessel wall should be checked carefully. Blackshear reviewed the medical literature of catheter perforations, which have confirming x-rays, and found that an incident angle of the catheter to vessel wall greater than 40° was more likely to perforate.

The following variables must also be considered in selecting appropriate catheter and length:

1. Patient history
2. Size and age of patient
3. Access site available
4. Unusual anatomical variables
5. Proposed use and duration of treatment plan

Catheter Tip Positioning

Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.

Suggested Catheter Maintenance

Catheter entry site must be prepared and maintained in a manner consistent with standard procedure for central venous catheterization. After catheter placement and prior to use, tip position and lumen patency should be confirmed by free aspiration of venous blood. **If blood is not freely aspirated, catheter tip position should be immediately reevaluated by physician.** If catheter is not to be used immediately, its lumen should be maintained by continuous saline or heparinized saline drip or locked with a catheter locking solution. **NOTE:** If CLC-2000, MicroClave or other needleless adapters approved for saline only lock are used, saline only catheter lock may be used. Catheter heparinization should be determined by institutional protocol and clinical judgement. Heparin concentrations of 10 Units/ml to 100 Units/ml have been reported adequate to maintain lumen patency. Catheter lock should be reestablished after every use or at least every 24 hours if unused. Before using catheter lumen already locked with heparin, lumen should be flushed with twice the indicated lumen volume using normal saline. Lumen should be flushed with normal saline between administration of different infusates. After use, lumen should again be flushed with twice the indicated lumen volume using normal saline before reestablishing catheter lock. Strict aseptic technique must be adhered to while using and maintaining catheter.

INSTRUCTIONS FOR USE

Activating AQ® Hydrophilic Coating

1. Flush catheter with heparinized saline solution or sterile water. **NOTE: The catheter may be trimmed if a shorter length is required.**
2. Attach a syringe with heparinized saline solution or sterile water to the Luer lock fitting of the catheter obturator holder.
3. Inject enough solution to wet the obturator surface entirely. This will activate the AQ® coating, making the obturator surface very lubricious.
4. Remove the obturator from its holder and insert it into the catheter, locking it into position. If the catheter has been trimmed, only advance the obturator to the distal end of the catheter. **NOTE: If the surface of the obturator becomes dry after removal from the holder, wetting with additional heparinized saline or sterile water will renew the hydrophilic effect.**
5. The catheter/obturator assembly can now be introduced as described in the following information.

Catheter Placement

1. After prepping the access site, introduce the thinwall percutaneous entry needle into the vessel. (**Fig. 1**)
2. Slide the Safe-T-J® wire guide straightener (positioned on the distal tip of the wire guide) over the "J" portion of the wire guide. Pass the straightened wire guide through the needle; advance the wire guide 5-10 cm into the vessel. (**Fig. 2**)
3. Leaving the wire guide in place, withdraw the needle and safety wire guide straightener. (**Fig. 3**) If necessary, enlarge the puncture site with a #11 scalpel blade.
4. Introduce the sheath-introducer assembly over the wire guide. With a twisting motion, advance the assembly into the vessel. (**Fig. 4**)
5. Leaving the sheath in place, remove the introducer and wire guide. (**Fig. 5**) (To prevent inadvertent air aspiration, place thumb or finger over the cuffed proximal end of the sheath after removing the introducer and wire guide.)
6. Introduce the catheter/obturator into the sheath; advance the catheter into position. (**Fig. 6**)
7. By grasping the two knobs of the sheath and pulling outward and upward at the same time, peel the sheath away from the catheter. (**Fig. 7**) **NOTE: Remove obturator after catheter is in final position.**
8. Secure catheter to skin and dress in standard fashion.
9. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.