|                 |   |
|-----------------|---|
| <b>en</b><br>2  | <b>Fuhrman Pleural Drainage Set</b><br>Instructions for Use                   |
| <b>bg</b><br>4  | <b>Комплект за плеврален дренаж Fuhrman</b><br>Инструкции за употреба         |
| <b>cs</b><br>7  | <b>Souprava pro Fuhrmanovu pleurální drenáž</b><br>Návod k použití            |
| <b>da</b><br>10 | <b>Fuhrman pleuradræningsæt</b><br>Brugsanvisning                             |
| <b>de</b><br>13 | <b>Fuhrman Pleuradrainage-Set</b><br>Gebrauchsanweisung                       |
| <b>el</b><br>16 | <b>ΣΕΤ παροχέτευσης υπεζωκοτικής<br/>κοιλότητας Fuhrman</b><br>Οδηγίες χρήσης |
| <b>es</b><br>19 | <b>Equipo de drenaje pleural Fuhrman</b><br>Instrucciones de uso              |
| <b>fi</b><br>22 | <b>Fuhrman-keuhkopussidreneerausetti</b><br>Käyttöohjeet                      |
| <b>fr</b><br>24 | <b>Set de drainage pleural Fuhrman</b><br>Mode d'emploi                       |
| <b>hu</b><br>27 | <b>Fuhrman pleurális drenázskészlet</b><br>Használati utasítás                |
| <b>it</b><br>30 | <b>Set per drenaggio pleurico di Fuhrman</b><br>Istruzioni per l'uso          |
| <b>nl</b><br>33 | <b>Fuhrman-pleuradrainageset</b><br>Gebruiksaanwijzing                        |
| <b>no</b><br>36 | <b>Fuhrman pleuralt drenasjesett</b><br>Bruksanvisning                        |
| <b>pl</b><br>39 | <b>Zestaw do drenażu jamy opłucnowej<br/>Fuhrman</b><br>Instrukcja użycia     |
| <b>pt</b><br>42 | <b>Conjunto de drenagem pleural Fuhrman</b><br>Instruções de Utilização       |
| <b>sk</b><br>44 | <b>Súprava na pleurálnu drenáž Fuhrman</b><br>Návod na použitie               |
| <b>sl</b><br>47 | <b>Komplet za plevralno drenažo Fuhrman</b><br>Navodila za uporabo            |
| <b>sv</b><br>50 | <b>Fuhrman-set för pleuradränage</b><br>Bruksanvisning                        |
| <b>tr</b><br>52 | <b>Fuhrman Plevral Drenaj Seti</b><br>Kullanma Talimatı                       |

## FUHRMAN PLEURAL DRAINAGE SET

**Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.**

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

The Fuhrman Pleural Drainage Set consists of a pigtail catheter with a Luer fitting on the proximal hub for connectivity, and components to aid in placement and use, including access needle, dilator, three-way stopcock, multipurpose tubing adapter and wire guide.

### Performance Characteristics

A variety of drainage sets are available across the device range to accommodate various patient anatomies. The catheters are offered in French sizes between 5 Fr and 12 Fr with an effective shaft length of 15 cm. The catheters have the following features:

- Radiopacity – the shaft material allows for visibility under fluoroscopy
- Hub – allows for standard Luer connection to accessory drainage components
- Depth markings – show user the depth in which the pigtail has been placed (only available on WCE sets/trays)

### Device Compatibility

The Fuhrman Pleural Drainage Set is compatible with:

- Saline
- Standard Luer syringe
- Mechanical drainage device (e.g., Heimlich valve, water seal suction, etc.)

### Patient Population

The Fuhrman Pleural Drainage Set should only be used in patients who require the placement of a drainage catheter into the pleural space as deemed appropriate by their physician or healthcare professional.

### Intended Users

The device is intended for use by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

### Contact with Body Tissue

The access needle, wire guide, dilator, and catheter will cross the skin, subcutaneous tissue, intercostal muscle, and parietal pleura to reach the pleural space. The three-way stopcock and multipurpose tubing adapter do not enter the patient's body.

### Operating Principle

The catheter is a single pigtail device with outer diameters ranging from 5 to 12 Fr and an effective length of 15 cm. The catheter can drain air or fluid from the pleural space. To accomplish this, the catheter is inserted into the pleural space via the Seldinger technique (access needle, wire guide, and dilators are provided with the set). The catheter sideports and endhole allow for pleural drainage, which is achieved through gravity or suction. For suction, the drainage catheter is compatible with drainage systems, such as a water seal suction apparatus, an underwater seal bottle or a similar mechanical drainage device, so long as the maximum negative pressure does not exceed 20 cmH<sub>2</sub>O (1961 Pa). Do not connect the drainage catheter device directly to wall suction.

### INTENDED USE

After percutaneous insertion via the Seldinger technique, the pigtail catheter of the Fuhrman Pleural Drainage Sets serves as a conduit to drain air or fluid from the pleural space to outside the body.

### INDICATIONS FOR USE

The indications for use of the pleural drainage catheters include pneumothorax and pleural effusion.

### DEVICE AND CLINICAL BENEFITS

- Complete or partial drainage of pleural air
- Complete or partial drainage of pleural fluid

### CONTRAINDICATIONS

None.

## WARNINGS

- Over-insertion of the access needle, wire guide, dilator, or catheter may result in serious or life-threatening injury. Clinical judgment and consideration of patient size and anatomy should guide appropriate catheter selection (e.g., French size and length) and insertion.
- The wire guide should always advance without impedance to avoid perforation of any surrounding structures.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

## PRECAUTIONS

- Imaging guidance is suggested.
- Ensure stopcock is in the ON position to the valve to prevent introduction of air.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

## MR SAFETY INFORMATION

The Fuhrman Pleural Drainage Catheter is MR Safe.

The presence of the Fuhrman Pleural Drainage Catheter may produce an image artifact.

## POTENTIAL ADVERSE EVENTS

- Air leak
- Bleeding/hemorrhage
- Death
- Dyspnea
- Fistula
- Hemothorax
- Hypoxia/hypoxemia
- Infection
- Injury to surrounding tissues/organs/structures
- Pain/discomfort
- Pleural effusion
- Pneumothorax
- Re-expansion pulmonary edema
- Subcutaneous emphysema
- Tube breakage/fracture
- Tube kinking
- Tube leakage
- Tube malfunction (dislocation, dislodgement, displacement)
- Tube occlusion/obstruction/blockage
- Vasovagal reaction or hypotension

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and out of sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

## INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

## INSTRUCTIONS FOR USE

1. Prepare the chest using standard aseptic technique and infiltrate the insertion area with local anesthesia.
2. Direct the tip of the needle towards the desired catheter placement and insert the needle under the skin.
3. Gently advance the needle.

**WARNING:** Over-insertion of the needle, wire guide, dilators, and/or catheter may result in serious harm to the patient.

**NOTE:** Due to anatomic variations in chest wall thickness, insertion depth of access needle, wire guide, dilator, and catheter will vary from patient to patient.

4. Attach a syringe to the needle. Verify placement into the pleural space by aspiration of air or fluid into syringe. Remove syringe from the needle once placement is verified.
5. Introduce the wire guide through the needle and gently advance it into the pleural cavity.  
**NOTE:** The wire guide should advance without resistance.  
**NOTE:** Maintain control of wire guide throughout the procedure.
6. Remove the needle, leaving the wire guide in place. Then, if needed, make a small incision with the scalpel. Advance the supplied dilator over the wire guide to dilate the tract to facilitate catheter introduction.
7. Remove the dilator, leaving the wire guide in place, then introduce the catheter over the wire guide.
8. Advance the catheter into position. Remove the wire guide and attach the multipurpose adapter to the Luer fitting on the catheter hub.  
**NOTE:** If the use of three-way stopcock is desired, the three-way stopcock may be attached directly to the Luer fitting on the catheter hub and the multipurpose adapter may then be attached to the three-way stopcock.
9. Secure catheter in position at the entry site by using a bio-occlusive dressing or suturing if desired.
10. Consider post-procedure x-ray to verify catheter placement.
11. If necessary, the catheter may be flushed with saline.
12. Catheter care and maintenance should follow the facility's protocol.
13. Duration of use should be less than 30 days.

## DISPOSAL OF DEVICES

After the procedure, this device, which contains sharps, may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

## REFERENCES

These instructions for use are based on experience from physicians and/or published literature. Refer to your local Cook sales representative for information on available literature.

## PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions contraindications, measures to be taken, and limitations of use that the patient should be aware of.

## SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used in.

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© COOK 2025
### Cook Medical Europe Ltd.
### William A. Cook Australia Pty. Ltd.