COOK

MEDICAL

CE
0088

EN
6

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C - T - P T B S - R E V 7


Fig. 1

1. Dolphin BT™ Ciaglia Balloon-Assisted Tracheostomy Introducer
2. Cook Inflation Device

1. Dolphin BT™ Ciaglia ballon-assisteret trakeostomi-indfører
2. Cook-inflationspistol

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
Fig. 8


Fig. 9

DOLPHIN BT™ CIAGLIA BALLOON-ASSISTED TRACHEOSTOMY INTRODUCER SET/TRAY

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Dolphin BT Ciaglia Balloon-Assisted Tracheostomy Introducer is a device used to facilitate percutaneous entry into the trachea for placement of a tracheostomy tube. A separate, sterile tracheostomy tube is also included in an optional set.

The set consists of a balloon-tipped catheter loading dilator assembly, Cook Inflation Device, (**Fig. 1**) wire guide, introducer needle, FEP sheath introducer needle, needle holder cup, 14 French dilator, gauze pads, disposable syringe, measuring tape, disposable scalpel, lubricating jelly and large full-body drape with clear plastic window.

A tray is also available in the USA only, which contains the set components and other items necessary for a bedside procedure, including lidocaine, 22, 25 and 18 gage needle, double swivel connector, chlorhexidine/alcohol prep solution, suture with needle, CSR wrap and prep tray.

INTENDED USE

The Dolphin BT Ciaglia Balloon-Assisted Tracheostomy Introducer is intended for controlled elective subcricoid insertion of a tracheostomy tube.

Tracheostomy tube placement, performed by a physician, using the technique as described in the following text, is done in a controlled setting (i.e., ICU) with the assistance of an anesthetist or respiratory therapist and nursing personnel.

CONTRAINDICATIONS

- Emergency tracheostomy tube placement
- Pediatric applications
- Patients with enlarged thyroids
- Non-palpable cricoid cartilage
- Non-intubated patients
- PEEP (Positive End Expiratory Pressure) value greater than or equal to 20
- Uncorrected coagulopathy

WARNINGS

- Excessive force applied to the manifold of the Blue Dolphin's loading dilator assembly can result in breakage and consequently procedural failure and loss of adequate control of the patient's airway.
- Only physicians trained and experienced in balloon-assisted tracheostomy technique should use this device.
- Performance of balloon-assisted tracheostomy under bronchoscopic guidance is recommended to reduce chance of paratracheal insertion and to determine the intratracheal position of the needle, wire, balloon, dilator and tracheostomy tube.
- Before balloon-assisted tracheostomy is attempted, the patient's airway must be secured with an endotracheal tube.
- Maintain alignment via safety positioning marks of wire guide and catheter during dilatation to prevent trauma to posterior wall of the trachea.
- Anatomic anomalies may make the procedure difficult to perform. The presence of anomalous blood vessels may cause excessive bleeding during procedure.
- Tracheostomy tube position should be verified by a chest X-ray or other method.
- In small or short patients (40 to 50 kilograms or less and/or individuals 4 feet or less in height) care must be taken not to injure the airway at the level of the carina.
- Do not exceed 11 atm. Over-inflation may cause rupture of the balloon, with resultant damage to the tracheal wall. Use of a pressure gauge is required to monitor inflation pressures.
- The catheter lumen marked "WIRE" should not be used for injection.

PRECAUTIONS

- Do not resterilize.
- The balloon is constructed of a heat-sensitive material. Do not heat or attempt to shape the catheter tip.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to the expiration date specified on the package.
- The catheter is not intended for the delivery of stents.
- Follow strict adherence to aseptic technique during placement of this device.
- Always confirm access into trachea by air bubble aspiration.

- Generous lubrication to the surface of the tracheostomy tube itself will facilitate smooth placement.
- A red sharps needle holder has been provided in this tray. Push any sharps into the foam after use. Do not reuse the sharps as particulate matter may adhere. Entire cup must be discarded into an appropriate sharps container after completion of the procedure.
- A greater degree of blunt dissection may be required in patients with strong connective tissue, pre-tracheal fascia or subcutaneous tissue for successful tracheostomy tube insertion.
- DO NOT pre-inflate the balloon to check for air leaks. This may potentially cause inadequate inflation during the procedure.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

Patient Preparation

1. Raise the inspired oxygen concentration to 100% and adjust the mechanical ventilator to ensure adequate ventilation during the procedure.
2. Monitor continuous pulse oximetry, airway pressures and exhaled tidal volumes to ensure adequate ventilation.
3. Provide adequate sedation and analgesia to control patient comfort.
4. Place the patient in the tracheostomy position. Position a pillow under the shoulders to permit full extension of head and neck. The head of the patient's bed may be elevated 30-40 degrees at the physician's discretion. If possible, perform a screening ultrasound of the neck to identify aberrant vessels that may be within the predicted path of the tracheostomy tube.
5. Prep and drape the anterior neck area.
6. Prepare the included Cook Inflation Device in standard manner and purge to remove all air from the device and tubing. Fill the device with a minimum of 20 ml of sterile saline.
7. After generously lubricating the loading dilator, load the tracheostomy tube onto the loading dilator. (**Fig. 5**)
8. Ensure that all tubing connections on the balloon catheter are tight. Attach the prepared Cook Inflation Device to the balloon port on the balloon catheter assembly. Ensure that there is no air in the Cook Inflation Device before connecting to the balloon port.
9. Slide the lock lever to the right (engage) to prevent inadvertent balloon inflation.

The Procedure

10. Palpate landmark structures (thyroid notch, cricoid cartilage) to ascertain proper location for intended tracheostomy tube placement. Ideally, tube placement is made at the level between the first and second tracheal cartilages or between the second and third tracheal cartilages whenever feasible. (**Fig. 2**)
11. Introduce local anesthesia into the desired tracheostomy site.
12. Make a 2 cm vertical or horizontal skin incision over the intended site of tracheostomy tube placement. No further sharp dissection is necessary; however, some blunt dissection can be performed using a hemostat if preferred. **NOTE:** Skin incision must be full 2 cm for the balloon to properly inflate and subsequently advance into the trachea. **NOTE:** A greater degree of blunt dissection may be required in patients with strong connective tissue, pre-tracheal fascia or subcutaneous tissue for successful tracheostomy tube insertion.
13. Instruct the anesthetist/respiratory therapist to loosen the fixation tapes of the endotracheal tube, deflate the cuff, and withdraw the tube approximately 1 cm under direct visualization so that the endotracheal tube cuff is just below the vocal cords. It may be necessary to make changes in tidal volume, frequency, etc., to compensate for air leaks. **Continuous oximetry monitoring should be employed.** Because the endotracheal tube is not secure in this location, it must be held in place.
14. Make a puncture in the trachea with the 18 gage introducer needle or FEP sheath needle mounted on a partially fluid-filled syringe. The needle stick should be made at the level between the first and second tracheal cartilages or between the second and third tracheal cartilages whenever feasible. Note that tracheal puncture is performed by directing the needle, in the midline, posterior and caudad. (**Fig. 3**)
15. When advancing the needle forward, verify entrance into the tracheal lumen by aspiration on the fluid-filled syringe resulting in air bubble return. **NOTE:** It is important not to impale the endotracheal tube with the needle. This can be checked by having the anesthetist/respiratory therapist gently move the endotracheal tube in and out 1 cm. If impaled, the needle is seen and felt to also move and it will be necessary to withdraw the needle, adjust the position of the endotracheal tube, and reinsert the needle.
16. Remove the syringe (and needle, if using the FEP sheath) and introduce the wire guide several centimeters through the FEP sheath or introducer needle and into the trachea. **NOTE:** The wire guide will not pass through the needle of the FEP sheath introducer needle assembly. If resistance is met, do not force the wire guide. Instead, re-confirm correct catheter

- placement bronchoscopically, then advance the wire guide into the trachea to the distal skin level marking on the wire guide.
17. Remove the FEP sheath or introducer needle while maintaining the wire guide position within the tracheal lumen. (**Fig. 4**)
 18. Maintaining the wire guide position, advance the short 14.0 French introducing dilator over the wire guide. Using a slight twisting motion, dilate the initial access site into the trachea. Remove the dilator while maintaining the wire guide in position with the skin level mark on the wire guide at its proper level.
 19. With the balloon fully deflated, advance the balloon catheter and tracheostomy tube assembly as a unit over the wire guide while maintaining the position of the distal wire guide marking at skin level. (**Fig. 6**) Align the proximal end of the balloon catheter at the mark on the proximal portion of the wire guide. This will ensure the distal end of the balloon catheter is correctly aligned on the wire guide, preventing possible trauma to the posterior tracheal wall during subsequent manipulations.
 20. With the balloon still deflated, advance the balloon catheter assembly with the wire guide until the balloon is halfway into the trachea. Ensure that the proximal wire guide marking is in alignment with the proximal tip of the balloon catheter hub. **NOTE:** The black mark on the balloon catheter indicates the mid-point of the balloon. Correct placement of the black mark should be at the anterior tracheal wall.
 21. While maintaining the visual reference points and positioning relationships of the wire guide and balloon catheter assembly, inflate the balloon. **NOTE:** Digital stabilization of the balloon catheter should be maintained during initial inflation of the balloon. Slide the lock lever on the Cook Inflation Device to the left (disengage) to allow free movement of the plunger. Inflate the balloon firmly and then slide the lock lever back to the right (engage). This will maintain the pressure and allow higher inflation pressures to be reached. Using the palm grip, rotate clockwise until the desired pressure is reached. **NOTE:** To avoid balloon rupture, do not exceed 11 atm. (**Fig. 7**)
 22. Maintain maximum inflation pressure for 10-20 seconds. **NOTE:** It is very important to rotate the palm grip as necessary in order to maintain the maximum 11 atm during the full 10-20 second inflation. **NOTE:** Assure adequate ventilation during maximum balloon inflation by adjusting ventilatory settings accordingly. The 2 cm skin incision allows the balloon on the catheter to slightly over-dilate the tracheal entrance site to a size appropriate for passage of the tracheostomy tube of choice. **NOTE:** Balloon can allow for repeated inflations if necessary.
 23. Slide the lock lever to the left (disengage) and immediately deflate the balloon until all infused saline is withdrawn. Slide the lock lever to the right (engage) to prevent inadvertent balloon inflation.
 24. Advance the deflated balloon catheter/loading dilator/tracheostomy tube/wire guide assembly as a unit into the trachea. (**Fig. 8**) The assembly should be directed perpendicular to the axis of the trachea during insertion for uniform dilation between tracheal cartilage. **NOTE:** The wire guide must always lead the balloon catheter assembly to prevent possible trauma to the posterior tracheal wall during dilation and advancement. Care should be taken to keep the balloon assembly properly aligned with the mark on the proximal portion of the wire guide. This will assure that the tip of the balloon assembly does not advance beyond the distal tip of the wire guide within the trachea. Dual cannula tracheostomy tubes may also be placed using this technique. The inner cannula must be removed for introduction. Follow tracheostomy tube manufacturer's instructions for testing of balloon cuff and inflation system prior to insertion.
 25. Once the tracheostomy tube is within the tracheal lumen, the assembly may be directed caudad. Advance the tracheostomy tube to its flange. (**Fig. 9**)
 26. When the tracheostomy tube is in place, slowly remove the balloon catheter/loading dilator/wire guide assembly. **If resistance is met during withdrawal of balloon, ensure that all saline is returned to Cook Inflation Device.** If continued resistance is met, consider removing entire assembly, including tracheostomy tube.
 27. If using a dual cannula tracheostomy tube, insert the inner cannula at this point. Connect the tracheostomy tube to the ventilator and inflate the balloon cuff. Prior to removal of the endotracheal tube, test ventilation through tracheostomy tube. **NOTE:** Insert the bronchoscope into the newly inserted tracheostomy tube to confirm correct airway placement.
 28. Perform suction of the tracheostomy tube to determine if any significant bleeding or possible obstruction exists which has not been noted to this point.
 29. If necessary, one suture may be used to close the bottom of the initial incision.

Retrieval Technique in the Event of Balloon Rupture

In the event of intraoperative balloon rupture, follow the steps listed below.

1. Attempt to withdraw the balloon-tipped catheter and loading dilator assembly while leaving the wire guide in place. If some level of dilatation was accomplished before balloon rupture occurred, removal of the ruptured balloon-tipped catheter and loading dilator assembly may be successful. If the balloon-tipped catheter begins to stretch or if resistance is noted, attempted withdrawal should be stopped immediately.

2. Remove the portion of the balloon remaining outside of skin level (proximal portion) from the catheter by trimming the balloon material away, leaving the catheter itself intact. **NOTE:** Care must be taken not to damage the catheter during the trimming process.
3. Advance the balloon-tipped catheter and loading dilator assembly over the *in situ* wire guide into the stoma that the ruptured balloon failed to dilate. Using the loading dilator, dilate the stoma to a size that will permit removal of the ruptured balloon.
4. Once the stoma is dilated to the size of the loading dilator, the ruptured balloon catheter may be removed over the *in situ* wire guide.
5. Leaving the wire guide in place allows a second balloon-tipped catheter and loading dilator assembly to be placed. Alternatively, a Ciaglia Blue Rhino Percutaneous Tracheostomy Introducer may be utilized over the *in situ* wire guide. **NOTE:** The wire guide included in the Dolphin BT set has a smaller diameter than the wire guide provided in the Ciaglia Blue Rhino Percutaneous Tracheostomy Introducer set.
6. Do not attempt to place the tracheostomy tube through the partially formed stoma until full dilatation is successfully achieved with either a second Dolphin BT device or a Ciaglia Blue Rhino device, which is suitable for the selected tracheostomy tube.

Post Placement

Elevate the head of the patient's bed 30-40 degrees if possible.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

Zgoda M, Berger R: "Balloon Facilitated Percutaneous Tracheostomy: A Novel Technique," *Chest*; October 2003; Supplement 130S-131S.

Zgoda M, Berger R: "Balloon-Facilitated Percutaneous Dilational Tracheostomy Tube Placement : Preliminary Report of a Novel Technique," *Chest*; November 2005; 128; 3688-3690.

Zgoda M, Berger R: "Tracheostomy in the Critically Ill Patient: Who, When, and How? A Review." *Clinical Pulmonary Medicine*; March 2006; 13(2); 111-120.

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Bloomington, IN 47404 U.S.A.