- EN**
11 **Blue Rhino® G2-Multi Percutaneous Tracheostomy Introducer**
Instructions for Use
- CS**
15 **Perkutánní tracheostomický zavaděč Blue Rhino® G2-Multi**
Návod k použití
- DA**
19 **Blue Rhino® G2-Multi perkutan trakeostomiindfører**
Brugsanvisning
- DE**
23 **Blue Rhino® G2-Multi Einführhilfe für die perkutane Tracheostomie**
Gebrauchsanweisung
- EL**
28 **Εισαγωγέας διαδερμικής τραχειοστομίας Blue Rhino® G2-Multi**
Οδηγίες χρήσης
- ES**
32 **Introductor para traqueostomía percutánea Blue Rhino® G2-Multi**
Instrucciones de uso
- FR**
37 **Introduceur pour trachéotomie percutanée Blue Rhino® G2-Multi**
Mode d'emploi
- HU**
41 **Blue Rhino® G2-Multi perkután tracheosztómiás bevezető**
Használati utasítás
- IT**
46 **Introduttore per tracheostomia percutanea Blue Rhino® G2-Multi**
Istruzioni per l'uso
- NL**
50 **Blue Rhino® G2-Multi introducer voor percutane tracheostomie**
Gebruiksaanwijzing
- NO**
55 **Blue Rhino® G2-Multi perkutan trakeostomiinnføringsenhet**
Bruksanvisning
- PL**
59 **Introduktor do tracheostomii przezskórnej Blue Rhino® G2-Multi**
Instrukcja użycia
- PT**
64 **Introdutor para traqueostomia percutânea Blue Rhino® G2-Multi**
Instruções de utilização
- SV**
68 **Blue Rhino® G2-Multi införare för perkutan trakeostomi**
Bruksanvisning


Fig. 1

Fig. 4


Fig. 5


Fig. 6


Fig. 13

BLUE RHINO® G2-MULTI PERCUTANEOUS TRACHEOSTOMY INTRODUCER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Blue Rhino G2-Multi Percutaneous Tracheostomy Introducer Set/Tray consists of these primary components: an introducer needle, J-tipped wire guide, introducer dilator, guiding catheter, loading dilators, and single-staged Blue Rhino G2-Multi dilator. Dilation takes place in one step with the Blue Rhino G2-Multi dilator using the Seldinger technique. Proper alignment of the components is illustrated in (Fig. 1).

INTENDED USE

The Blue Rhino G2-Multi Percutaneous Tracheostomy Introducer Set/Tray is intended for percutaneous dilational tracheostomy for management of the airway in adults only. Tube placement, using the technique described herein, should be performed in a controlled setting (e.g., ICU or operating room) with the assistance of trained personnel.

CONTRAINDICATIONS

- Emergency tracheostomy tube placement
- Patients with enlarged thyroids
- Nonpalpable cricoid cartilage
- Pediatric applications
- Non-intubated patients
- PEEP (Positive End Expiratory Pressure) value greater than or equal to 20
- Uncorrected coagulopathy
- Pre-existing infection at the tracheostomy site
- Pre-existing malignancy at the tracheostomy site
- Unstable cervical spine fracture
- Previous surgery at the tracheostomy site (e.g., thyroidectomy)

WARNINGS

- Only physicians trained and experienced in percutaneous tracheostomy techniques should use this device.
- Exercise care to ensure that the components used in each step are properly positioned within the trachea. Improper placement of the components may lead to potentially life-threatening injury.
- Improper dilation technique or tracheostomy tube placement can lead to delayed complications (e.g., Corkscrew stenosis).
- Aseptic technique must be strictly adhered to during placement of this device.
- Prior to attempting percutaneous tracheostomy, the patient's airway must be secured with an endotracheal tube.
- Anatomic anomalies may make the procedure difficult to perform. The presence of anomalous blood vessels may cause excessive bleeding during the procedure.
- In small or short patients (weighing 40-50 kg or less, and/or 4 ft. or less in height), care must be taken not to injure the airway at the level of the carina.

PRECAUTIONS

- Bronchoscopic guidance is **strongly** recommended during placement of this device to reduce the likelihood of paratracheal insertion and to determine the intratracheal position of the needle, wire guide, dilators and tracheostomy tube.

- An ultrasound evaluation of the patient's neck prior to the procedure may aid in identification of anatomical variances.
- This product is intended for use by physicians trained and experienced in percutaneous tracheostomy techniques. Standard techniques for percutaneous placement of tracheostomy tubes should be employed.
- Do not resterilize components.
- Always confirm needle access into the trachea by air bubble aspiration or bronchoscopic visualization.
- Maintain safety positioning marks of the wire guide, guiding catheter and dilator during dilating procedure to prevent trauma to posterior wall of the trachea.
- The tracheostomy tube should fit snugly to the loading dilator for insertion.
- The generous lubrication of the loading dilator surface will enhance fit and placement of the tracheostomy tube.
- Loading dilators are designed to be inserted within a tracheostomy tube only. Loading dilators should not be used for creation of a stoma.
- Tracheostomy tube position should be verified by a chest X-ray.
- Take care not to advance the Blue Rhino G2-Multi Dilator too far into the patient's airway.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects. **NOTE:** The bronchoscopy adapter is the only component in this kit that contains phthalates.

POTENTIAL ADVERSE EVENTS

- Perforation of the trachea
- Failed tracheostomy tube placement
- Hypoxia

PRODUCT RECOMMENDATIONS

Recommended Loading Dilator sizes and associated tracheostomy tubes

Loading Dilator	Compatible ISO Tracheostomy Tube Sizes
6.5 mm	6.5 mm
7.0 mm	7.0 mm
7.5 mm	7.5 mm
8.0 mm	8.0 mm
8.5 mm	8.5 mm
9.0 mm	9.0 mm
10.0 mm	10.0 mm

Compatibility testing was performed with Shiley®¹ Flex and Shiley®¹ Evac tracheostomy tubes. When using another appropriately sized tracheostomy tube, ensure that the tracheostomy tube's tip fits snugly on the dilator.

¹ "Shiley" is a registered trademark of Covidien LP

INSTRUCTIONS FOR USE

Patient Preparation

1. Following the tracheostomy tube manufacturer's instructions, test the balloon cuff and inflation system.
2. Place the patient in the tracheostomy position. (**Fig. 4**) Position a pillow under the shoulders to permit full extension of head and neck. The head of the patient's bed may be elevated 30-40 degrees at the physician's discretion.

3. Use ventilator changes and sedation to control patient respiration. A PEEP (Positive End Expiratory Pressure) level of 5-10 is recommended.
4. Instruct the clinician managing the airway to loosen the fixation tapes of the in-place endotracheal tube and deflate the cuff, making necessary changes in tidal volume, frequency, etc., to evaluate compensation needed for air leak. Reinflate the endotracheal tube cuff. **Continuous oximetry monitoring should be employed.**
5. Prep and drape the anterior neck area.
6. Generously lubricate the surface of the appropriately sized loading dilator and load the tracheostomy tube onto the dilator. Ensure that the tracheostomy tube's tip fits snugly on the dilator. (Fig. 2) Ensure that the balloon is completely deflated. Thoroughly lubricate tracheostomy tube assembly.

Tracheostomy Procedure

NOTE: Dual cannula tracheostomy tubes may also be placed using the following technique. The inner cannula must be removed for introduction. Always check the fit of the dilator to the tracheostomy tube prior to insertion.

1. Palpate landmark structures (thyroid notch, cricoid cartilage) to ascertain proper location for tracheostomy tube placement. **Access and, ultimately, tube placement is ideally made at the level between the first and second tracheal cartilages or between the second and third tracheal cartilages whenever feasible.** (Fig. 3)
2. After introducing local anesthesia, make a 1.5-2.0 cm skin incision (vertical or horizontal) at the chosen insertion site. (Fig. 5)
3. If desired, use a curved mosquito clamp to gently dissect vertically and transversely down to the anterior tracheal wall. (Fig. 6) With a finger tip, dissect the front of the trachea, in the midline, free of any tissues and identify the cricoid cartilage. Displace the isthmus of the thyroid downward, if present. **NOTE:** An adequate skin incision and blunt dissection of the subcutaneous tissue can minimize the need for excessive force and torque throughout the procedure. Excessive force and rotation may lead to long-term complications (e.g., stenosis).
4. Deflate the endotracheal tube cuff and withdraw to an appropriate distance above the insertion site, yet still within the trachea. Re-inflate the cuff once the proper position of the endotracheal tube has been reached.
5. Attach a syringe half-filled with fluid to the introducer needle and seek the tracheal air column by directing the needle, in the midline, posterior. **Verify entrance into the tracheal lumen via aspiration on the syringe resulting in air bubble return.** (Fig. 7) **Alternatively, if using bronchoscopy, visualize the needle entering the trachea.**

NOTE: It is important that the needle not impale the endotracheal tube. To ensure that the endotracheal tube is not impaled, gently move it in and out 1 cm. If the tube is impaled, the needle will be seen and felt to also move. If this occurs, it will be necessary to withdraw the needle, pull back the endotracheal tube, and then reinsert the needle. **NOTE:** Proper positioning and alignment may help minimize complications (e.g., stenosis).

6. With the needle tip positioned in the trachea, local anesthesia may be injected (if necessary).
7. When free flow of air is obtained, with no impalement of the endotracheal tube, remove the inner needle of the introducer needle assembly and advance the outer FEP sheath several millimeters. **NOTE:** If using an introducer needle without a sheath, proceed to step 9.
8. Attach a syringe to the FEP sheath and re-confirm position within the tracheal lumen by visualizing free flow of air into the syringe when aspirated. (Fig. 8) **Alternatively, re-confirm position by visualizing the FEP sheath in the trachea with the bronchoscope. Remove the syringe.**
9. Introduce the J-tipped wire guide several cm into the trachea. (Fig. 9) **NOTE:** The wire should advance freely, without resistance. If resistance is encountered, do not force wire guide. Confirm correct FEP sheath or introducer needle placement via bronchoscopy, then advance wire guide into the trachea until the distal mark on the wire guide reaches skin level.
10. Remove the FEP sheath or introducer needle while maintaining wire guide position within the tracheal lumen. (Fig. 9)
11. Maintaining the wire guide's position at the skin level mark, dilate the initial tracheal access site by advancing the short, 14 French introducer dilator over the wire guide with a slight twisting motion. (Fig. 10)
12. Remove the dilator while maintaining wire guide position.
13. Activate the hydrophilic coating by immersing the distal end of the Blue Rhino G2-Multi dilator in sterile water or saline.

14. Advance the Blue Rhino G2-Multi dilator and the guiding catheter as a unit over the wire guide, while maintaining wire guide position. **(Fig. 11) NOTE:** Align the proximal end of the guiding catheter at the mark on the proximal portion of the wire guide. **(Fig. 1)** This will ensure that the distal end of the guiding catheter is properly positioned back on the wire guide, preventing possible trauma to the posterior tracheal wall during subsequent manipulations. **NOTE:** Bronchoscopic guidance may also prevent possible trauma to the posterior tracheal wall.
15. Begin to dilate the tracheal access site by advancing the guiding catheter and Blue Rhino G2-Multi dilator into the trachea. **To properly align the dilator on the wire guide/guiding catheter assembly, position the proximal end of the dilator at the single positioning mark on the guiding catheter. This will ensure that the distal tip of the dilator is properly positioned at the safety ridge on the guiding catheter to prevent possible trauma to the posterior tracheal wall during introduction.** While maintaining the visual reference points and positioning relationships of the wire guide, guiding catheter and dilator, advance them **as a unit** to the skin level mark on the Blue Rhino G2-Multi dilator. **(Fig. 11) NOTE:** Proper positioning and alignment may help minimize complications (e.g., stenosis).
16. Advance and pull back the dilating assembly several times to effectively dilate the tracheal access site. **NOTE:** The wire guide must always lead the dilator and the guiding catheter assembly to prevent possible trauma to the posterior tracheal wall during dilation. Care should be taken to keep the guiding catheter assembly properly aligned with the mark on the proximal portion of the wire guide. This will ensure that the tip of the guiding catheter assembly does not advance beyond the distal tip of the wire guide within the trachea.
17. Remove the Blue Rhino G2-Multi dilator, leaving the wire guide/guiding catheter assembly in position. **Respiratory air leak through the tracheostomy stoma should be noted to confirm intratracheal location of the wire guide and guiding catheter.**
18. Advance the tracheostomy tube (loaded on the dilator) over the wire guide/guiding catheter assembly to the safety ridge of the guiding catheter, then advance wire guide, guiding catheter, loading dilator and tracheostomy tube **as a unit** into trachea. **(Fig. 12)**
NOTE: The assembly should be directed perpendicular to the axis of the trachea during insertion for uniform dilation between tracheal cartilages. Once the tracheostomy tube is within the tracheal lumen, the assembly may be directed caudad. **NOTE:** Proper positioning and alignment may help minimize complications (e.g., stenosis).
19. Advance the tracheostomy tube to its flange. Remove the dilator, guiding catheter and wire guide, leaving the tracheostomy tube in place. **(Fig. 13) NOTE:** At this point, the bronchoscope may be inserted into the tracheostomy tube to confirm correct placement. **NOTE:** If using a dual cannula tracheostomy tube, insert the inner cannula at this point.
20. Inflate the tracheostomy tube balloon cuff. Connect the tracheostomy tube to the ventilator. Confirm position of the tracheostomy tube via standard methods (e.g., capnography, breath sounds, etc.).
21. Deflate and remove the endotracheal tube.
22. Perform suction to determine if any significant bleeding or possible obstruction exists that has not been noted to this point.
23. If necessary, one suture may be taken at the bottom of the initial incision.

Post-Placement

Follow hospital protocol for post-tracheostomy care and maintenance.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2020
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd