- EN**
5-6
Thal-Quick Chest Tube Sets and Trays
Instructions for Use
- CS**
6-7
Sady a tácky hrudních sond Thal-Quick
Návod k použití
- DA**
7-9
Thal-Quick thoraxtubesæt og bakker
Brugsanvisning
- DE**
9-10
Thal-Quick Thoraxdrainage-Sets
Gebrauchsanweisung
- EL**
10-12
ΣΕΤ και δίσκοι θωρακικού σωλήνα Thal-Quick
Οδηγίες χρήσης
- ES**
12-14
Equipos y bandejas de catéteres torácicos Thal-Quick
Instrucciones de uso
- FR**
14-15
Sets et plateaux de drain thoracique Thal-Quick
Mode d'emploi
- IT**
15-17
Set standard e set completi con tubo toracico Thal-Quick
Istruzioni per l'uso
- KO**
17-18
Thal-Quick 흉관 세트 및 트레이
사용 지침
- NL**
18-19
Thal-Quick thoraxdrainsets en -pakketten
Gebruiksaanwijzing
- NO**
20-21
Thal-Quick thoraxslangesett og brett
Bruksanvisning
- PT**
21-22
Conjuntos e tabuleiros de tubos torácicos Thal-Quick
Instruções de utilização
- SV**
23-24
Thal-Quick toraxdränset och -brickor
Bruksanvisning


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5

THAL-QUICK CHEST TUBE SETS AND TRAYS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Thal-Quick Chest Tube Set (**TQTS-**) consists of introducer needle, wire guide, dilators, coaxial chest tube inserter and clear catheter with radiopaque strip (**TQT-**).

The letters "Y" or "YJ" in the order number prefix, (e.g., **TQTSYJ-**) indicate a tray. Trays include the same items as sets, plus additional accessories for use in placement of a chest tube. The Thal-Quick is also available in a double-lumen configuration. The Double-Lumen Thal-Quick Chest Tube is a 24 French chest tube with a 17 gage second lumen that can be used for infusions or aspirations.

INTENDED USE

The Thal-Quick Chest Tube is intended for percutaneous introduction of a chest tube for pleural fluid drainage.

CONTRAINDICATIONS

None known

WARNINGS

Over-insertion of the coaxial chest tube inserter and clear catheter assembly may result in serious or life-threatening injury. Clinical judgment, and consideration of patient size and anatomy should guide appropriate chest tube selection (e.g., size French, length, number of sideports) and insertion.

NOTE: To maintain an air tight drainage system and avoid further complications, it is important that all sideports of the chest tube are positioned within the pleural space. The centimeter markings on the chest tube denote the distance from the most proximal side port, see **Fig. 5** (i.e., this is not the distance from tip of catheter). These markings can be used as a reference while establishing appropriate device positioning.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in percutaneous pleural drainage techniques. Standard techniques for placement of chest tubes should be employed.
- Manipulation of products requires fluoroscopic, CT scan or ultrasound guidance.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Prepare the chest tube insertion site with antiseptic solution and drape in a standard fashion.
NOTE: A mid-axial approach is recommended.
2. After injection of local anesthetic, make a small skin incision slightly larger than the diameter of the chest tube.
3. Attach the introducer needle to a syringe filled with local anesthetic (e.g., lidocaine), and administer additional local anesthetic to the skin and underlying tissues.
4. Advance the introducer needle over the superior border of the rib and into the pleural space. Fluid or air should be aspirated to verify an intrapleural position.
WARNING: Over insertion of the needle and/or dilators may result in serious harm to the patient.
NOTE: The needle should be introduced and directed with appropriate orientation inferiorly or superiorly.
5. When the appropriate drainage site has been identified, advance the soft "J" end of the wire guide through the needle and into the pleural space. (**Fig. 1**)
NOTE: The wire guide should pass through the needle and advance into the pleural space without resistance.
NOTE: The skin level marker on the distal end of the wire guide can be used as an insertion depth marker on patients with normal anatomy.
NOTE: Due to anatomical variations in chest wall thickness, insertion depth of introducer needle, wire guide and dilators will vary from patient to patient.
6. Remove the needle, leaving the wire guide in place.
7. Liberally inject additional local anesthetic into the intercostal muscles surrounding the wire guide.
8. While maintaining the wire guide position, dilate the tract and opening into the pleural space by advancing, in sequence small to large, the supplied dilators over the wire guide. Introduction into the pleural space is facilitated by rotating and advancing the dilators in line with the wire guide to prevent its kinking. (**Fig. 2**)

WARNING: Over insertion of the needle and/or dilators may result in serious harm to the patient.

NOTE: It is important to have enough length of wire guide exiting the access site to facilitate controlled introduction of the dilators. The dilators only have to enter the pleural space to their full diameter and should never be advanced beyond the distal end of the wire guide.

9. With the wire guide still positioned within the pleural space, advance the chest tube inserter/chest tube assembly over the wire guide and into the pleural space. (**Fig. 3**)

NOTE: If resistance is encountered during chest tube assembly insertion, determine the cause of resistance and take necessary action to relieve resistance before proceeding.

NOTE: It is important to advance the chest tube assembly into the pleural space in the same line as the wire guide. This will make introduction easier and avoid kinking of the wire guide.

NOTE: The distal end of the chest tube inserter should not be advanced beyond the distal end of the wire guide.

10. Remove the wire guide and chest tube inserter, leaving the chest tube in place. (**Fig. 4**)

NOTE: It is important that all sideports of the chest tube are positioned within the pleural space.

11. The chest tube can now be sutured to the skin and is ready for use.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd