EN
3

Thoracentesis Sets

Instructions for Use

DA
3

Pleurapunktursæt

Brugsanvisning

DE
4

Pleurapunktions-Sets

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EL
5

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C - T - T H O R A - R E V 1

THORACENTESIS SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Thoracentesis Sets include an access needle, wire guide, centesis catheter (with or without a stiffening cannula) and connecting tubes.

INTENDED USE

The Thoracentesis Set is intended to remove fluid from the pleural space.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The products are intended for use by physicians trained and experienced in thoracentesis. Standard thoracentesis techniques should be employed.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.
- Lung puncture may result in an air embolus, which could lead to ischemia or infarction of major organs including the brain or cardiac system.

INSTRUCTIONS FOR USE

1. Prep lateral chest using standard aseptic technique.
2. Infiltrate the seventh and eighth intercostal spaces in the mid-axillary line with local anesthesia.
3. Remove the inner stylet from the 18 gage introducer needle and attach a syringe to the needle.
4. Introduce the needle 45 degrees cephalad in the mid-axillary line of the seventh intercostal space.
5. After entering the pleural space or after aspirating fluid, remove the syringe from the needle.
6. Insert the wire guide through the needle and advance the wire guide approximately 10–15 cm into the pleural space.
NOTE: The wire guide should advance without resistance.
7. Leaving the wire guide in place, remove the needle.
8. Introduce the catheter, with stiffening cannula in place, over the wire guide and advance into the pleural cavity.
9. Holding the catheter steady, withdraw the wire guide while at the same time twisting and withdrawing the stiffening cannula.
10. Attach the connecting tube to the catheter and aspirate.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2013
Bloomington, IN 47404 U.S.A.