EN 3	Cook TPN® Catheters and Sets Instructions for Use
DA 7	Cook TPN®-katetere og -sæt Brugsanvisning
DE 11	Cook TPN® Katheter und Sets Gebrauchsanweisung
EL 15	Καθετήρες Cook TPN® και Σετ Καθετήρων Οδηγίες Χρήσης
ES 19	Equipos y Catéteres Cook TPN® Instrucciones de Uso
FR 23	Cathéters et sets TPN® Cook Mode d'emploi
IT 27	Catetere e set Cook TPN® Istruzioni per l'uso
NL 31	Cook TPN® katheters en sets Gebruiksaanwijzing
PT 35	Cateteres e Conjuntos Cook TPN® Instruções de utilização
SV 39	Cook TPN®-katetrar och set Bruksanvisning

Cook TPN® Catheters and Sets

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION


Fig. 1

Single-Lumen Catheter			
French Size	Lumen Diameter (inch)	Lumen Volume (ml)*	Length (cm)
3.0	.022	.4	65
4.0	.026	.5	65
5.0	.036	.7	65
6.5	.040	1.1	90
9.5	.062	2.2	90

*All lumen volumes are for untrimmed length. If catheter is trimmed, volumes should be adjusted accordingly.

The Cook TPN Single-Lumen Catheter Set consists of a single-lumen catheter with a clamp and a synthetic fiber cuff for tissue in-growth (**Fig. 1**), a wire guide, a needle and a Peel-Away® sheath introducer for percutaneous introduction of the catheter. The set also includes a needle-less leak-proof adapter and a syringe.


Fig. 2

Double-Lumen Catheter					
French Size	Small Lumen Diameter (inch)	Small Lumen Volume (ml)*	Large Lumen Diameter (inch)	Large Lumen Volume (ml)*	Length (cm)
5.0D	.019	.5	.019	.5	65
6.0D	.018	.6	.026	.7	65
7.0D ¹	.019 x .045	.7	.025 x .050	.9	65
9.0D ¹	.025 x .065	1.6	.036 x .074	2.1	90
10.0D	.043	1.7	.043	1.7	90
12.0D ¹	.034 x .070	2.0	.047 x .082	2.7	90

*All lumen volumes are for untrimmed length. If catheter is trimmed, volumes should be adjusted accordingly.

¹ D-Shaped Lumens

The TPN Double-Lumen Catheter Set has a double-lumen catheter (**Fig. 2**), which incorporates two separate, non-communicating vascular access lumens within a single catheter body. Exit ports of individual lumens are separated. Set contents are as in the single-lumen set, including a Luer lock plug.

ENGLISH


Fig. 3

Triple-Lumen Catheter					
French Size	Small Lumen Diameter (inch)	Small Lumen Volume (ml)*	Large Lumen Diameter (inch)	Large Lumen Volume (ml)*	Length (cm)
8.0T	.029	1.2	.043	1.8	90

*All lumen volumes are for untrimmed length. If catheter is trimmed, volumes should be adjusted accordingly.

The Cook TPN Triple-Lumen Catheter Set has a triple-lumen catheter (Fig. 3), which incorporates three separate, non-communicating vascular access lumens within a single catheter body. Set contents are as in the sets above.

The catheter sets can be ordered with different suffixes added to the order number (e.g., -REDO). The Redo TPN catheter has circumferential suture rings to allow the catheter to be secured to the skin with suture during tissue in-growth into cuff.

INTENDED USE

Cook TPN Catheters are intended for the administration of parenteral nutrition, chemotherapy, other therapeutic drugs and blood.

The catheters are intended for use by physicians trained and experienced in the placement of central venous catheters using percutaneous entry (Seldinger) technique. Standard Seldinger technique for placement of percutaneous vascular access sheaths, catheters, and wire guides should be employed.

CONTRAINDICATIONS

None known

WARNINGS

- Possible allergic reactions should be considered.
- Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system. Tip position should be verified by X-ray and monitored on a routine basis. Periodic lateral view X-ray is suggested to assess tip location in relation to vessel wall. Tip position should appear to be as parallel to the vessel wall as possible, as angles larger than 40° are more likely to perforate.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify the attending physician immediately.

PRECAUTIONS

- Volumes indicated are for untrimmed catheters. If catheter is trimmed, volumes should be adjusted accordingly.
- Manipulation of products requires fluoroscopic control.
- The catheter size should be as small as the use will allow, as larger diameter catheters have more tendency to promote clots.
- Select the puncture site and length of catheter needed by assessing patient anatomy and condition.
- The use of ECG and/or fluoroscopy is suggested for accurate catheter placement.
- Extreme caution must be used in placement and monitoring.
- To guarantee extrapericardial location, the catheter tip should not be advanced beyond the innominate vein or initial segment of the superior vena cava.

- Placement into the subclavian vein may result in compression of the catheter by the clavicle and first rib. Excessive compression may result in catheter damage, including rupture or catheter embolus. Compression of the catheter and consequent risk may be minimized by placement lateral to the junction of the clavicle and first rib or use of alternative venous sites.
- Silicone catheters are not designed for power injection. Catheter rupture may occur. Use of a 10 ml syringe or larger will reduce the risk of catheter rupture.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

PRODUCT RECOMMENDATIONS

Catheter Selection

The following variables must be considered in selecting the appropriate catheter and length:

- Patient history
- Size and age of patient
- Access site available
- Unusual anatomic variables
- Proposed use and duration of treatment

INSTRUCTIONS FOR USE

1. Prepare and maintain the catheter entry site in a manner consistent with standard procedures for catheter placement.
2. Make a small incision in the skin just beneath the right clavicle in the mid-clavicular line and puncture the subclavian vein with the needle attached to the syringe.
3. Remove the syringe, introduce the wire guide through the needle and advance it into the vessel.
4. While maintaining wire guide position, withdraw the needle.

Subcutaneous Catheter Placement with Curved Hemostat

1. Select the catheter exit site, usually at the medial aspect of the breast at the nipple line, and make a transverse 1 cm incision.
2. A subcutaneous tunnel is created from the wire guide entry to the exit site using a curved hemostat.
3. Grasp the catheter tip with hemostats and pull it through the tunnel while positioning the cuff approximately 1 cm inside the exit wound.
4. The catheter may then be trimmed to proper length and flushed with heparinized saline. **NOTE:** In order to guarantee extrapericardial location, the catheter tip should not be advanced beyond the innominate vein or the initial segment of the superior vena cava.
5. Introduce the Peel-Away sheath introducer with a gentle twisting motion over the wire guide into the superior vena cava.
6. Leaving the sheath in place, remove the introducer and wire guide. To prevent inadvertent air aspiration, place thumb or finger over the cuffed proximal end of the sheath after removing the introducer and wire guide.
7. Introduce the catheter into the sheath and advance into position.
8. By grasping the two knobs of the sheath and pulling outward and upward at the same time, peel the sheath away from the catheter.
9. Verify the catheter tip position radiographically and close all wounds. The catheter can then be sutured to the skin and dressed in a standard fashion.

ENGLISH

SUGGESTED CATHETER MAINTENANCE

After catheter placement and prior to use, tip position and lumen patency should be confirmed by free aspiration of venous blood. **If blood is not freely aspirated, catheter tip position should be immediately reevaluated by the physician.** If catheter is not to be used immediately, its lumen should be maintained by continuous saline or heparinized saline drip or locked with heparinized saline solution. Normal saline lock is permissible if utilizing the CLC2000 injection cap. Catheter heparinization should be determined by institutional protocol and clinical judgement. Heparin concentrations of 10 Units/ml to 100 Units/ml have been reported adequate to maintain lumen patency. Catheter lock should be reestablished after every use or at least every 24 hours if unused. Before using catheter lumen already locked with heparin, lumen should be flushed with twice the indicated lumen volume using normal saline. Lumen should be flushed with normal saline between administration of different infusates. After use, lumen should again be flushed with twice the indicated lumen volume using normal saline before reestablishing heparin or saline lock. Strict aseptic technique must be adhered to while using and maintaining catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.