- EN**
3 **Pericardiocentesis Sets**
Instructions for Use
- CS**
4 **Sady pro perikardiocentézu**
Návod k použití
- DA**
5 **Pericardiepunktursæt**
Brugsanvisning
- DE**
6 **Perikardpunktions-Sets**
Gebrauchsanweisung
- EL**
7 **ΣΕΤ ΠΕΡΙΚΑΡΔΙΟΚΕΝΤΗΣΗΣ**
Οδηγίες χρήσης
- ES**
8 **Equipos de pericardiocentesis**
Instrucciones de uso
- FR**
9 **Sets de péricardiocentèse**
Mode d'emploi
- IT**
10 **Set per pericardiocentesi**
Istruzioni per l'uso
- NL**
12 **Pericardiocentesesets**
Gebruiksaanwijzing
- PT**
13 **Conjuntos de pericardiocentese**
Instruções de utilização
- SV**
14 **Perikardiocentes-set**
Bruksanvisning

PERICARDIOCENTESIS SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Pericardiocentesis Sets may include an introducer needle (with or without an inner stylet), J-tipped wire guide, dilator, scalpel, centesis catheter (with or without a stiffening cannula) connection tubing, and multipurpose tubing adapter. Some Pericardiocentesis Sets also include an alligator clip cable for electrocardiographic monitoring during needle insertion.

INTENDED USE

Pericardiocentesis Sets are intended to remove fluid from the pericardial sac.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- These products are intended for use by physicians trained and experienced in pericardiocentesis. Standard pericardiocentesis techniques should be employed.
- A TFE-coated wire guide must be used with the catheter in this set.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.
- Lung puncture may result in an air embolus, which could lead to ischemia or infarction of major organs, including the brain or cardiac system.
- Only the wire guide included in this set should be used with this catheter.

INSTRUCTIONS FOR USE

1. Prep the left chest using standard aseptic technique.
2. Infiltrate the left xiphocostal area with local anesthesia.

NOTE: Sets that contain an alligator clip cable can be used for constant electrocardiographic (ECG) monitoring. To do so, attach an insulated V lead to one end of the alligator clip cable. **CAUTION: Ensure that the electrocardiograph is properly grounded prior to attaching the alligator clip cable to the lead.** Attach the other end of the alligator clip cable to the metal hub of the introducer needle.

3. With the needle directed toward the right shoulder and elevated 30 degrees from the horizontal plane at the left xiphocostal margin, advance the needle until it enters the pericardial space.

NOTE: If ECG monitoring is being performed, the ECG will demonstrate a current of injury as the needle makes contact with the epicardium. It may also be possible to feel cardiac movement transmitted through the needle.

NOTE: If the needle has an inner stylet, it can be removed prior to introducing the needle, if continuous aspiration is desired while advancing the needle. If the alligator clip cable is in use, it should remain attached to a metal component of the needle.

4. Once the needle is in position, remove the inner stylet (if present) from the needle.
5. Slide the Safe-T-J wire guide straightener (positioned on the distal end of the wire guide) over the "J" portion of the wire guide and seat the straightener in the hub of the needle.
6. Insert the wire guide through the needle and advance the wire guide into the pericardial sac. **NOTE:** The wire guide should advance without resistance.
7. Leaving the wire guide in place, remove the needle.
8. Create a small skin incision. If a dilator is included in set, the dilator should be advanced over the wire guide and removed prior to insertion of the catheter.
9. Introduce the catheter (with stiffening cannula, if applicable) over the wire guide and advance until the tip of the catheter is posterior to the muscle fascia.
10. If the set contains a catheter with stiffening cannula, remove the cannula with a twisting motion, while holding the catheter and wire guide in place.
11. Advance the catheter into the pericardium.

NOTE: The wire guide should always extend beyond the catheter tip.

12. Remove the wire guide.
13. If applicable, attach the appropriate connecting tube and prepare to aspirate.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt the product is sterile. Store in a dark, dry cool, place. Avoid extended exposure to light. Upon removal from package inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd