EN
4

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
Fig. 1


Fig. 2

UNCOATED AND HEPARIN-COATED CENTRAL VENOUS CATHETERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Uncoated and heparin-coated central venous catheters are available in either a single-lumen or multi-lumen configuration. Multi-lumen catheters incorporate separate, non-communicating vascular access lumens within a single catheter body.

Some catheters are available with heparin coating for short-term protection against thrombus formation. The central venous catheter with heparin coating is designated by the suffix -BH in the reorder number. (Example: C-UTLMY-701J-RSC-BH)

Minipuncture™ Access Sets and Trays (denoted by -MNP)

Used for the insertion of 5.0 French or larger central venous catheters (single and multi-lumen) through the use of a 20 gage extra-thinwall access needle and a specially designed .025 inch diameter Amplatz Extra Stiff Wire Guide using percutaneous entry (Seldinger) technique. All Minipuncture central venous catheters have a removable inner cannula, which must be removed after introduction.

Micropuncture® Access Sets and Trays (denoted by -MP)

Used for the insertion of central venous catheters (single and multi-lumen) through the use of a 22 gage access needle and a .018 inch diameter wire guide using percutaneous entry (Seldinger) technique. Most catheters for adult applications have a removable inner cannula, which must be removed after introduction.

INTENDED USE

The central venous catheter is designed for treatment of critically ill patients and is suggested for:

1. Continuous or intermittent drug infusions
2. Central venous blood pressure monitoring (CVP)
3. Acute hyperalimentation
4. Blood sampling
5. Delivery of whole blood or blood products
6. Simultaneous, separate infusion of drugs

CONTRAINDICATIONS

Do not use heparin-coated central venous catheter in patients with a history of heparin-induced thrombocytopenia.

WARNINGS

- Complications arising from the use of this device can result in serious injury or death; catheter tip can erode or perforate vascular walls.
- If heparin-induced thrombocytopenia is suspected, the heparin-coated catheter should be immediately removed.
- Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of central venous system. Tip position should be verified by x-ray and monitored on a routine basis. Periodic lateral view x-ray is suggested to assess tip location in relation to vessel wall. Tip position should appear to be parallel to vessel wall.
- To avoid vascular injury, do not use excessive force when advancing dilators. Use the smallest size dilator catheter placement will allow. Wire guide must always lead dilator by several centimeters. Do not advance dilator more than a few centimeters into the vessel.
- **Do not power inject contrast medium through catheter.** Catheter rupture may result. Use of 10 ml or larger syringe will reduce the risk of catheter rupture.
- To distend great vessels and to prevent inadvertent air aspiration during catheter insertion, patient should be placed in Trendelenburg position.
- Development of a hypersensitivity reaction should be followed by removal of the catheter and appropriate treatment at the discretion of the attending physician.

PRECAUTIONS

- The catheter is intended for use by physicians trained and experienced in the placement of central venous catheters using percutaneous entry (Seldinger) technique. Standard Seldinger technique for placement of percutaneous vascular access sheaths, catheters and wire guides should be employed during the placement of a central venous catheter.
- Do not re-sterilize catheter.
- Do not cut, trim or modify catheter or components prior to placement or intraoperatively.
- Patient movement can cause catheter tip displacement. Use should be limited to controlled hospital

situations. Catheters placed via an antecubital vein have shown forward tip movement of up to 10 cm with motion of the extremity. Catheters placed from either a jugular or subclavian vein have demonstrated forward tip movement of 1-3 cm with neck and shoulder motion.

- Catheter should not be used for long-term indwelling applications.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify attending physician immediately.
- Catheter should not be used for chronic hyperalimentation.
- Select puncture site and length of catheter needed by assessing patient anatomy and condition.
- Left subclavian and left jugular veins should be used only when other sites are not available
- Use of ECG, ultrasound and/or fluoroscopy is suggested for accurate catheter placement.
- For heparin coated devices, standard flushing procedures are recommended.

PRODUCT RECOMMENDATIONS

Catheter Size and Puncture Site

Preliminary reports indicate that catheter size can influence clotting; larger diameter catheters have more tendencies to promote clots. As reported by Amplatz and others, clot formation has less relation to type of catheter material than to size of catheter.¹ The angle of the catheter tip to the vessel wall should be checked carefully. Blackshear reviewed the medical literature of catheter perforations, which have confirming x-rays, and found that an incident angle of the catheter to the vessel wall greater than 40 degrees was more likely to perforate.²

Another critical factor that can cause a catastrophic event is the choice of puncture site. Findings by Tocino and Watanabe indicate that the left subclavian and left jugular veins should be avoided when practical.³ Eighty percent of the perforations or erosions were found when these vessels were used. In addition, they have observed that the tip curve of a wedged catheter can be detected with lateral view x-ray.

The above discussion is meant to be a guide for catheter size and puncture site. As more data become available, other causal factors may become evident, but present information suggests that:

- The catheter size should be as small as the use will allow.
- Left subclavian and left jugular veins should be used only when other sites are not available.

The following variables should be considered when selecting appropriate catheter and length:

- Patient history
- Size and age of patient
- Access site available
- Unusual anatomical variables
- Proposed use and duration of treatment plan

Standard Catheter Lengths Available for Pediatric Use

These suggested lengths must be viewed only as a guideline.

French Size	Length	Access Site
2.5	2.5/5 cm	Internal and external jugular veins; for patients weighing less than 10 kg
3.0	5/8 cm	Internal and external jugular veins; for patients weighing less than 10 kg
3.0	8 cm	Internal and external jugular veins
4.0/5.0	5 cm	Internal and external jugular veins; for patients weighing less than 10 kg
4.0/5.0	8 cm	Internal and external jugular veins; for patients weighing 10-40 kg
4.0/5.0	12 cm	Right subclavian vein; for patients weighing 10-40 kg
4.0/5.0	15 cm	Left subclavian vein; for patients weighing 10-40 kg

Standard Catheter Lengths Available for Adult Use

These suggested lengths must be viewed only as a guideline.

French Size	Length	Access Site
5.0/6.0/6.3	15 cm	Internal and external jugular veins
5.0/6.0/6.3	20 cm	Right subclavian vein
5.0/6.0/6.3	25 cm	Left subclavian vein
7.0/7.5/9.0/9.5	15 cm	Internal and external jugular veins
7.0/7.5/9.0/9.5	20 cm	Right subclavian vein
7.0/7.5/9.0/9.5	25 cm	Left subclavian vein

Single-Lumen Information

French Size	Length	Equivalent Gage	Minimum Lumen Volume for Heparin Lock
2.5	2.5 cm	21	.05 ml
2.5	5 cm	21	.05 ml
3.0	5 cm	20	.1 ml
3.0	8 cm	20	.1 ml
4.0	12 cm	18	.1 ml
5.0	15 cm	16	.2 ml
5.0	20 cm	16	.3 ml
5.0	25 cm	16	.3 ml
6.0	15 cm	14	.3 ml
6.0	20 cm	14	.3 ml
6.0	25 cm	14	.4 ml
6.3	15 cm	14	.4 ml
6.3	20 cm	14	.5 ml
6.3	25 cm	14	.6 ml

Double-Lumen Information

French Size	Lumens	Equivalent Gage	Minimum Lumen Volume
4.0*	#1	20	.2 ml
	#2	23	.1 ml
4.0**	#1	20	.3 ml
	#2	22	.3 ml
5.0	#1	20	.2 ml
	#2	20	.2 ml
7.0	#1	16	.6 ml
	#2	20	.2 ml
7.5	#1	14	.8 ml
	#2	21	.2 ml
9.5	#1	13	2.1 ml
	#2	18	.4 ml

*Pediatric polyurethane Double Lumen Central Venous Catheter (C-UDLM, C-DLM)

**Pediatric Double Lumen Soft Central Venous Catheter (C-SDLM)

Triple-Lumen Information

French Size	Lumens	Equivalent Gage	Minimum Lumen Volume
5.0	#1	18	.3 ml
	#2	23	.2 ml
	#3	23	.2 ml
7.0	#1	16	.5 ml
	#2	18	.3 ml
	#3	18	.3 ml
9.0	#1	14	1.1 ml
	#2	18	.4 ml
	#3	18	.4 ml

Suggested Lumen Utilization: Double-Lumen

- **#1 Distal exit port (endhole)** — whole blood or blood product delivery and sampling; any situation requiring more flow rate; CVP monitoring; medication delivery. **It is strongly recommended that this lumen be used for all blood sampling.**
- **#2 Proximal exit port** — medication delivery; acute hyperalimentation.

Suggested Lumen Utilization: Triple-Lumen

- **#1 Distal exit port (endhole)** — whole blood or blood product delivery and sampling; any situation requiring more flow rate; CVP monitoring; medication delivery. **It is strongly recommended that this lumen be used for all blood sampling.**
- **#2 Middle exit port** — medication delivery; acute hyperalimentation.
- **#3 Proximal exit port** — medication delivery.

NOTE: This catheter is not recommended for right atrial air aspiration because of inadequate sideport surface area. For this procedure, order the Bunegin-Albin Air Aspiration Set. (Order Number: **C-VAE-581J-SCJ**)


Suggested Catheter Maintenance

- Catheter entry site must be prepared and maintained in a manner consistent with standard procedure for central venous catheterization.
- To prevent clotting or possibility of air embolus, the double-lumen's #2 lumen, and the triple-lumen's #2 and #3 lumens should be filled with heparinized saline solution (100 units of heparin per cc of saline is usually adequate) prior to catheter introduction.
- After catheter is placed and prior to use, tip position and lumen patency should be confirmed by free aspiration of venous blood. **If blood is not freely aspirated, physician should immediately reevaluate catheter tip position.**

- Any unused lumens should be maintained with continuous saline or heparinized saline drip or locked with heparinized saline solution. Heparin-locked lumens should be reestablished at least every 8 hours.
- Before using any lumen already locked with heparin, lumen should be flushed with twice the indicated lumen volume using normal saline. Lumens should be flushed with normal saline between administrations of different infusates. After use, lumen should again be flushed with twice the indicated lumen volume using normal saline before reestablishing heparin lock.
- Strict aseptic technique must be adhered to while using and maintaining catheter.

Access Sites of Choice

- High Internal Jugular
- External Jugular
- Low Internal Jugular
- Supraclavicular
- Intraclavicular
- Femoral (not shown)

INSTRUCTIONS FOR USE

- If applicable, remove the Luer-lock end cap from each extension.
- Prepare the catheter for insertion by flushing each of the lumens and clamping or attaching the injection caps to the appropriate extensions. Leave the distal extension uncapped for wire guide passage.
- Introduce thinwall percutaneous entry needle into vessel. Venous blood should be easily aspirated to confirm position of needle tip within vessel.
- Slide Safe-T-J® wire guide straightener (positioned on distal tip of wire guide) over "J" portion of wire guide. Pass straightened wire guide through needle; advance wire guide 5-10 cm into vessel. If straight wire is used, always advance soft, flexible end through needle hub and into vessel. **If resistance is encountered during wire guide insertion, do not force wire guide. Withdrawal of wire guide through needle should be avoided; breakage may result.**
- While maintaining wire guide position, withdraw needle and Safe-T-J wire guide straightener.
- Enlarge puncture site with number 11 scalpel blade, if required. If dilation is required, dilator can be advanced over wire guide and removed prior to insertion of central venous catheter.
CAUTION: To avoid vascular injury, do not use excessive force when advancing dilators. Use the smallest size dilator catheter placement will allow. Wire guide must always lead dilator by several centimeters. Do not advance dilator more than a few centimeters into the vessel.
- Measure catheter to be used against patient to determine approximate length of catheter needed from puncture site to central venous tip position.
- Introduce the central venous catheter over wire guide. While maintaining wire guide position, advance catheter into vessel with a gentle twisting motion. (Fig. 1)
NOTE: Do not advance catheter tip beyond distal tip of wire guide. Always have wire guide leading during catheter placement. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.
- After catheter is in position, remove wire guide. (Fig. 2) **Venous blood should be easily aspirated.** Winged hub can now be sutured into place. If catheter is not introduced to its full length, additional suture should be carefully placed around catheter and affixed to the skin at entry site if movable suture wing is not included. This will help prevent backward or forward catheter movement. Lumens should now be flushed with 5-10 cc normal saline prior to use or establishment of heparin lock.
NOTE: A wire guide that is at least twice as long as the catheter is recommended for catheter exchange procedure.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

1. Amplatz K. "A Simple Non-thrombogenic Coating." Investigative Radiology, 1971; 6:280-289.
2. Blackshear RH, Gravenstein N. "Critical Angle of Incidence for Delayed Vessel Perforation by Central Venous Catheter: A Study of *In Vivo* Data." Annals of Emergency Medicine, 1992; 21:242.
3. Tocino IM, Watanabe A. "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition." American Journal of Roentgenology, 1986; 146:487-490.

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1. Amplatz K. "A Simple Non-thrombogenic Coating." *Investigative Radiology*, 1971; 6:280-289.
2. Blackshear RH, Gravenstein N. "Critical Angle of Incidence for Delayed Vessel Perforation by Central Venous Catheter: A Study of In Vivo Data." *Annals of Emergency Medicine*, 1992; 21:242.
3. Tocino IM, Watanabe A. "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition." *American Journal of Roentgenology*, 1986; 146:487-490.

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1. Amplatz K. "A Simple Non-thrombogenic Coating." Investigative Radiology, 1971; 6:280-289.
2. Blackshear RH, Gravenstein N. "Critical Angle of Incidence for Delayed Vessel Perforation by Central Venous Catheter: A Study of In Vivo Data." Annals of Emergency Medicine, 1992; 21:242.
3. Tocino IM, Watanabe A. "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition." American Journal of Roentgenology, 1986; 146:487-490.

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1. Amplatz K. "A Simple Non-thrombogenic Coating." *Investigative Radiology*, 1971; 6:280-289.
2. Blackshear RH, Gravenstein N. "Critical Angle of Incidence for Delayed Vessel Perforation by Central Venous Catheter: A Study of *In Vivo* Data." *Annals of Emergency Medicine*, 1992; 21:242.
3. Tocino IM, Watanabe A. "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition." *American Journal of Roentgenology*, 1986; 146:487-490.

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1. Amplatz K. "A Simple Non-thrombogenic Coating." *Investigative Radiology*, 1971; 6:280-289.
2. Blackshear RH, Gravenstein N. "Critical Angle of Incidence for Delayed Vessel Perforation by Central Venous Catheter: A Study of *In Vivo* Data." *Annals of Emergency Medicine*, 1992; 21:242.
3. Tocchio IM, Watanabe A. "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition." *American Journal of Roentgenology*, 1986; 146:487-490.

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1. Amplatz K. "A Simple Non-thrombogenic Coating." *Investigative Radiology*, 1971; 6:280-289.
2. Blackshear RH, Gravenstein N. "Critical Angle of Incidence for Delayed Vessel Perforation by Central Venous Catheter: A Study of In Vivo Data." *Annals of Emergency Medicine*, 1992; 21:242.
3. Tocchio IM, Watanabe A. "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition." *American Journal of Roentgenology*, 1986; 146:487-490.

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1. Amplatz K. "A Simple Non-thrombogenic Coating." *Investigative Radiology*, 1971; 6:280-289.
2. Blakeshear RH, Gravenstein N. "Critical Angle of Incidence for Delayed Vessel Perforation by Central Venous Catheter: A Study of *In Vivo* Data." *Annals of Emergency Medicine*, 1992; 21:242.
3. Tocino IM, Watanabe A. "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition." *American Journal of Roentgenology*, 1986; 146:487-490.

**MANUFACTURER**

COOK INCORPORATED
750 Daniels Way
Bloomington, IN 47404 U.S.A.

www.cookmedical.com
© COOK 2019

**EC REPRESENTATIVE**

Cook Medical Europe Ltd
O'Halloran Road
National Technology Park
Limerick, Ireland

2019-10
C_T_ULMBH_REV7