COOK

MEDICAL

CE
0088

EN
4

**Wayne Pneumothorax Set
For Trocar Placement**

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NO
17

**Wayne pneumothoraxsett
For trokarplassering**

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WAYNE PNEUMOTHORAX SET

For Trocar Placement

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Wayne Pneumothorax Set consists of a stainless steel needle obturator with Peel-Away® spacer, radiopaque polyurethane catheter, polyvinylchloride connecting tube, plastic three-way stopcock, and latex-free Cook Chest Drain Valve to be used for evacuation of air or fluid.

INTENDED USE

The Wayne Pneumothorax Set is intended for relief of simple, spontaneous, iatrogenic and tension pneumothorax.

CONTRAINDICATIONS

Not recommended for large fluid accumulation or hemothorax.

WARNINGS

None known


PRECAUTIONS

- This product is intended for use by physicians trained and experienced in the treatment of a pneumothorax. Standard techniques for placement of pneumothorax catheters should be employed.
- Ensure stopcock is in the ON position to the valve to prevent aspiration of air.
- Lung puncture may result in an air embolus, which could lead to ischemia or infarction of major organs, including the brain or cardiac system.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.


INSTRUCTIONS FOR USE

1. Preassemble catheter, connecting tube with stopcock, and Cook Chest Drain Valve. Attach Cook Chest Drain Valve in direction indicated by arrow on valve. **NOTE:** Chest Drain Valve may be obviated if catheter is to be connected to a water seal suction apparatus or similar mechanical suction device. Do not connect catheter directly to wall suction.
2. Prep the access site with appropriate antiseptic solution and drape in standard fashion. **NOTE:** The suggested insertion site is the fourth intercostal space at the anterior or mid-axillary line.
3. Introduce local anesthesia through skin and subcutaneous tissue down to pleura. Make an incision through skin only. **Anesthesia may be omitted in an emergency decompression.**
4. Without removing the Peel-Away spacer, fully straighten the curved catheter tip by advancing the needle obturator. (**Fig. 1**)
5. After straightening the catheter tip, remove the Peel-Away spacer and advance the needle until the point is showing at the catheter tip. **NOTE:** Care must be taken to fully straighten pigtail before advancing needle to avoid puncturing catheter. (**Fig. 2**)
6. Holding the base of the needle against the palm of the hand, push the catheter through the incision and over the rib into the pleural cavity. Once into the pleural cavity, withdraw the inner needle stylet slightly and advance the catheter slightly using the obturator.
7. The catheter may then be advanced over the obturator to the desired depth. The catheter pigtail will form as advanced. Depth may be guided by the 2.5 cm markers on the catheter. The first one begins 2.5 cm from the last sidehole. **NOTE:** If there is uncertainty as to whether the obturator and catheter tip are well within the pleural cavity as with a small pneumothorax, the inner needle stylet may be fully removed to determine whether there is a free flow of air.
8. Remove the obturator and confirm placement by valve movement and fluoroscopic or roentgenographic verification. **All connections must be secure and airtight.** Perform inspections of the catheter and connections regularly.
9. Secure catheter in position at the entry site by using a bio-occlusive dressing or suturing if desired.
10. The catheter and connected drain lines should be secured to the patient. Excessive tension on catheter connections (e.g., in instances of patient movement where the catheter is connected to a vacuum or drainage collection apparatus) may cause catheter/hub separation or accidental catheter dislodgement. To help prevent this from occurring, it is recommended to do one or both of the following:

- a. Secure the catheter hub to the patient's skin by placing tape, suture, or a catheter securement device at the location shown in the illustration below.

- b. Form a strain relief loop in the connecting tube, as shown in the illustration below, and secure the loop to the patient's skin with tape, suture, or catheter securement device.

11. To remove catheter, remove any catheter skin attachments. Have patient hold his/her breath in full inspiration. Quickly withdraw catheter and apply an occlusive dressing such as petrolatum-impregnated gauze. **A slight resistance will be noted as the curve of the catheter straightens.**

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.