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SV **Cook-Swartz dopplerflödessond med förlängningskabel**

69 Bruksanvisning

Cook-Swartz Doppler Flow Probe and Extension Cable

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Cook-Swartz Doppler Flow Probe contains at its distal end an ultrasound transducer (20 MHz) assembly that is attached to a silicone cuff which is designed to be secured around a blood vessel. Proximal to the cuff, a flexible wire connects the transducer to stress retention tabs and finally to a terminal connector. The Cook-Swartz Doppler Extension Cable connects the proximal connector of the Probe to either of the Cook Doppler Blood Flow Monitor's output channels.

INTENDED USE

For monitoring blood flow in vessels intraoperatively, and following reconstructive micro-vascular procedures, re-implantation, and free-flap transfers.

CONTRAINDICATIONS

None Known.

WARNINGS

The Cook-Swartz Doppler Flow Probe emits and detects 20 MHz signals. External 20 MHz signals may also be detected.

MR Unsafe – Do Not expose patient to an MRI procedure while Cook-Swartz Doppler Flow Probe is implanted. Substantial MRI-related heating of the Doppler Flow Probe may occur. Doppler Flow Probe must be removed prior to any MRI procedure.

The Cook-Swartz Doppler Flow Probe is not to be used with High Frequency surgical equipment.

PRECAUTIONS

Not intended for fetal use.

The Cook-Swartz Doppler Flow Probe is NOT FOR USE ON THE CENTRAL CIRCULATORY SYSTEM.

The Cook-Swartz Doppler Flow Probe should only be used with the Cook Doppler Blood Flow Monitor.

During use of all ultrasound devices, the operator should minimize the exposure of ultrasound energy to the patient using the principle of ALARA (As Low As Reasonably Achievable).

Avoid use of excessive force to remove the transducer assembly from the patient, which may cause injury to the blood vessel. If the transducer assembly can not be removed using gentle traction, the transducer assembly should be removed surgically.

Changes in audio signals produced by the Cook Doppler Blood Flow Monitor system should be immediately documented and reported to the responsible healthcare provider.

Caution: Do not remove the probe conductor wire and crystal assembly (leaving only the cuff in situ on the vessel), until vessel monitoring is completed (commonly 3-5 days). Probe conductor wire and crystal assembly placement must not exceed 29 days. Cuff alone may remain within the patient indefinitely.

Avoid the application of electrosurgical energy on or near a connected Cook-Swartz Doppler Flow Probe or extension cable as damage to the monitor may occur.

POTENTIAL ADVERSE EVENTS

Use of the Cook-Swartz Doppler Flow Probe involves potential risks normally associated with any implanted device, e.g., infection, perforation or laceration of vessels, erosion, implant rejection, or device dislodgement/migration.

Device specific risks include separation of the doppler crystal from the cuff, inability to percutaneously remove the crystal after monitoring is complete, loss of reception or transmission of ultrasound monitoring signal.

INSTRUCTIONS FOR USE

IPX7 - Ingress Protection Rating. The distal 28cm of the Cook Doppler Flow Probe which contains the crystal transducer is rated IPX7 - Water Immersible.

1. Confirm proper operation of the Cook Doppler Blood Flow Monitor system. (See Suggested Instructions for the Cook Doppler Blood Flow Monitor.)
2. Trim curved silicone cuff to accommodate vessel size. The cuff length should allow close approximation of the crystal to the vascular adventitia. Position the silicone cuff around the targeted blood vessel, ensuring that the transducer assembly is directed towards the targeted vessel.
3. Secure the free ends of the silicone cuff together around the target vessel, using appropriate sutures or clip to maintain cuff position.
4. Stabilize the location of the probe by suturing the braided wire adjacent to the skin incision site, leaving some slack to alleviate tension on the cuff transducer assembly. Enough wire length should be provided so that there is no tension on the vascular anastomosis. Loop any remaining braided wire, and suture or tape to the skin. The retention tab may be placed over the braided wire and sutured or taped to the skin providing additional strain relief.
5. Attach the metal connector end of the extension cable (DP-CAB01) into either of the two channels by pushing it firmly into place. Attach the proximal (red plastic) connector of the probe to the distal (red plastic) connector of the extension cable, aligning the black dots. Verify appropriate position of the Cook-Swartz Doppler Flow Probe by turning on the flow monitor and adjusting the volume until an adequate audible signal is obtained.

Note: The proximal connector of the Cook-Swartz Doppler Flow Probe should not be attached to the flow monitor/extension cable until the retention tabs are secured to the skin. This helps ensure that accidental tugging of the wires does not disrupt the attachment of the probe to the vessel.

6. If a strong audible signal is not identified, irrigate the crystal with saline at its interface with the blood vessel adventitia. During irrigation of the crystal, an audible signal from the monitor verifies proper function of the device.
7. For operation of the flow monitor, refer to the following suggested INSTRUCTIONS FOR USE for the Cook Doppler Blood Flow Monitor.
8. Following verification of proper function of the Cook-Swartz Doppler Flow Probe, close the incision site using standard techniques.
9. To remove the Cook-Swartz Doppler Flow Probe, first free the retention tab and braided wires to the skin by cutting the sutures (and/or removing the tape). Remove the probe by applying gentle traction to the braided wires at the skin entry site until the transducer assembly is withdrawn. (The silicone cuff remains in situ.)
10. Upon removal of the Cook-Swartz Doppler Flow Probe, examine the distal tip of the probe to ensure that the transducer assembly is present. In the unlikely event that the transducer assembly has become detached and remains in the cuff in the patient, the transducer assembly should be removed surgically.


HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

DISPOSAL

Upon completion of procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

TRANSPORT and STORAGE

Recommended Storage,
Shipping Temperature: -40°C - +60°C

Keep Dry

NON-CONDENSING

Recommended Humidity: non-condensing

Keep away
from sunlight

ACOUSTIC OUTPUT:

See following tables.

ACOUSTIC OUTPUT TABLES FOR TRACK 1

Cook-Doppler Flow Probe and Monitor System

OPERATING MODE(S)

Clinical Application	B	CD	M	PWD	CWD	A	Combined (specify)
Ophthalmic							
Other (Intraoperative)				X			
Cardiac							
Clinical Application							

OUTPUT DATA

Transducer Model: **20 MHz Swartz Doppler Flow Probe with Monitor System**
 Operating Mode: **PWD-Mode**
 Application(s): **Intraoperative**

DP-M350**Acoustic Output Reporting Table**

Transducer Model: 20 MHz 1mm			Operating Mode: Pulse Doppler (PD)				
Index Label	MI	TIS			TIB	TIC	
		Scan	Non-scan		Non-scan		
			Aaprt $\leq 1\text{cm}^2$	Aaprt $> 1\text{cm}^2$			
Maximum index value	0.0575	#	0.0175	-	0.0338	(a)	
Associated acoustic parameter	$P_{\text{e,a}}$ (MPa)	0.258					
	P (mW)		#	0.183		0.183	#
	min of $[P_a(z_s), I_{\text{ta},a}(z_s)]$ (mW)				-		
	z_s (cm)				-		
	z_{bp} (cm)				-		
	z_b (cm)					-	
	z at max. $I_{\text{pi},a}$ (cm)	0.150					
	$d_{\text{eq}}(z_b)$ (cm)					-	
	f_{awf} (MHz)	20.1	#	20.1	-	20.1	#
	Dim of A_{aprt}	X (cm)		#	0.100	-	0.100
Y (cm)			#	0.100	-	0.100	#
Other Information	t_d (μsec)	0.345					
	p_{rr} (Hz)	7.81E+4					
	p_r at max. I_{pi} (MPa)	0.286					
	d_{eq} at max. I_{pi} (cm)					0.0658	
	$I_{\text{pa},3}$ at max. MI (W/cm^2)	1.63					

Note 1: Information need not be provided for any formulation of *TIS* not yielding the maximum value of *TIS* for that mode.

Note 2: Information need not be provided regarding *TIC* for any TRANSDUCER ASSEMBLY not intended for transcranial or neonatal cephalic uses.

Note 3: *MI* and *TI* are less than 1.0 for all device settings.

(a) Intended use does not include cephalic so *TIC* is not computed.

No data reported.

SYMBOL INDEX

PWD:	Pulsed Wave Doppler
$I_{SPTA,3}$:	Derated spatial-peak, temporal-average intensity (milliwatts per square cm)
$I_{SPPA,3}$:	Derated spatial-peak, pulse-average intensity (watts per square cm)
W_O :	Ultrasound power (mW)
F_C :	Center frequency (MHz)
Z_{sp} :	Axial distance at which the reported parameter is measured (cm)
X_{-6} & Y_{-6} :	Respectively the in-plane (azimuthal) and out of plane (elevational) -6 dB dimensions in the x-y plane where Z_{sp} is found (cm)
PD:	Pulse duration (μ S)
PRF:	Pulse repetition frequency (Hz)
EBD:	Entrance beam dimensions for the azimuthal and elevational planes (cm ²)
a	Acoustic Attenuation Coefficient
Aaprt	-12 dB Output Beam Area
Deq	Equivalent Aperture Diameter
Deq	Equivalent Beam Diameter
Fawf	Acoustic Working Frequency
I_{pa}	Pulse-Average Intensity
I_{pi}	Pulse-Intensity Intergral
$I_{pi, a}$	Attenuated Pulse-Intensity Intergral
$I_{ta, a(z)}$	Attenuated Temporal-Average Intensity
<i>MI</i>	Mechanical Index
<i>P</i>	Output Power
<i>Pa</i>	Attenuated Output Power
<i>Pr</i>	Peak-Rarefactional Acoustic Pressure
<i>Pra</i>	Attenuated Peak-Rarefactional Acoustic Pressure
<i>prf</i>	Pulse Repetition Rate
<i>TIS</i>	Soft-Tissue Thermal Index
t_d	Pulse Duration
z	Distance from the source to a specified point
z_b	Depth for TIB
z_{bp}	Break-Point Depth
z_s	Depth for TIS

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© COOK 2019
Vandergrift, PA 15690 U.S.A.
COOK MEDICAL EUROPE LTD