- EN 2 **Evolution® RL Controlled-Rotation Dilator Sheath Set**
Instructions for Use
- CS 6 **Souprava dilatačních sheathů s řízenou rotací Evolution® RL**
Návod k použití
- DA 10 **Evolution® RL dilatatorsheathsæt med kontrolleret rotation**
Brugsanvisning
- DE 14 **Evolution® RL Dilatatorschleusen-Set mit kontrollierter Rotation**
Gebrauchsanweisung
- EL 18 **Σετ θηκαριών διαστολέα ελεγχόμενης περιστροφής Evolution® RL**
Οδηγίες χρήσης
- ES 22 **Equipo de vainas dilatadoras de giro controlado Evolution® RL**
Instrucciones de uso
- FR 26 **Set de gaines de dilatation à rotation contrôlée Evolution® RL**
Mode d'emploi
- HU 30 **Evolution® RL kontrollált forgatású dilatátorhüvely-készlet**
Használati utasítás
- IT 34 **Set di guaine dilatatrici a rotazione controllata Evolution® RL**
Istruzioni per l'uso
- LT 38 **„Evolution® RL“ kontroliuojamojo sukimo plėtiklio vamzdelių rinkinys**
Naudojimo nurodymai
- NL 42 **Evolution® RL dilatatorsheathset met rotatieregeling**
Gebruiksaanwijzing
- NO 46 **Evolution® RL dilatatorhylsesett med kontrollert rotasjon**
Bruksanvisning
- PL 50 **Zestaw koszułek rozszerzających Evolution® RL z regulacją obrotu**
Instrukcja użycia
- PT 54 **Conjunto de bainhas de dilatação de rotação controlada Evolution® RL**
Instruções de utilização
- RU 58 **Комплект расширителей-интродьюсеров с регулируемым вращением Evolution® RL**
Инструкция по применению
- SV 62 **Evolution® RL dilatatorhylsset med rotationskontroll**
Bruksanvisning

Evolution® RL Controlled-Rotation Dilator Sheath Set**Instructions for Use**

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Evolution RL Controlled-Rotation Dilator Sheath Set is comprised of an inner polymer sheath connected to a handle capable of mechanically rotating the sheath and an outer telescoping polymer sheath. The inner sheath has a stainless steel tip at its distal end.

INTENDED USE

The Evolution RL Controlled-Rotation Dilator Sheath Set is intended for use in patients requiring the percutaneous dilation of tissue surrounding cardiac leads, indwelling catheters and foreign objects.

Sheath Size	REF Sheath Set
9 Fr	LR-EVN-9.0-RL
11 Fr	LR-EVN-11.0-RL
13 Fr	LR-EVN-13.0-RL

Use With Other Devices

The Evolution RL Controlled-Rotation Dilator Sheath Set may be used in conjunction with the following catheter/lead extraction devices from Cook:

- Locking Stylet
- Lead Extender

Be sure to closely follow the "Suggested Instructions for Use" for each device used.

CONTRAINDICATIONS

None known

WARNINGS

- **When using sheaths, do not insert sheaths over more than one lead at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.**
- **When using a Locking Stylet:**
 - **Do not abandon a catheter/lead in a patient with a Locking Stylet still in place inside the catheter/lead. Severe vessel or endocardial wall damage may result from the stiffened catheter/lead or from fracture or migration of the abandoned stylet wire.**
 - **Do not apply weighted traction to an inserted Locking Stylet as myocardial avulsion, hypotension, or venous wall tearing may result.**
 - **Be aware that a lead that has a J-shape retention wire that occupies its inner lumen (rather than being outside the coil) may not be compatible with the Locking Stylet. Insertion of the Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.**
- **Weigh the relative risks and benefits of intravascular catheter/lead removal procedures in cases when:**
 - **The item to be removed is of a dangerous shape or configuration,**
 - **The likelihood of catheter/lead disintegration resulting in fragment embolism is high or vegetations are attached directly to the catheter/lead body.**
- **Catheter/lead removal devices should be used only at institutions with thoracic surgical capabilities.**
- **Catheter/lead removal devices should be used only by physicians knowledgeable in the techniques and devices for catheter/lead removal.**

PRECAUTIONS

- Prior to the procedure, consider the size of the catheter/lead in relation to the size of the LEAD EXTRACTION™ devices to determine possible incompatibility.
- If selectively removing catheters/leads with the intent to leave one or more chronic catheters/leads implanted intact, the non-targeted catheters/leads must be subsequently tested to ensure that they were not damaged or dislodged during the extraction procedure.

As a result of the risk of complications, doctors highly experienced in this procedure have advised the following:

PREPARATIONS

- Obtain a thorough patient history, including patient blood type. Appropriate blood products should be rapidly available.
- Ascertain the manufacturer, model number, and implant date of the catheter/lead to be removed.
- Perform radiographic/echocardiographic evaluation of catheter/lead condition, type, and position. Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, thoracotomy tray, and pericardiocentesis tray.
- Echocardiography and cardiothoracic surgery should be rapidly available.
- Prep/drape the patient's chest for possible thoracotomy; prep/drape the patient's groin for femoral approach or possible femoral approach.
- Establish back up pacing as necessary.
- Have available an extensive collection of sheaths, Lead Control Devices (Locking Stylet and Lead Extender), stylets to unscrew active fixation leads, snares, and accessory equipment.

PROCEDURE

- Use fluoroscopic monitoring during ALL catheter/lead and sheath manipulations. Monitor ECG and arterial blood pressure continuously throughout the procedure and during recovery.
- If using sheaths or sheath sets, including the Evolution or Evolution RL Controlled-Rotation Dilator Sheath Set, the following precautions should be followed:
 - Prior to using sheaths including the Evolution or Evolution RL Controlled-Rotation Dilator Sheath Set it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.
 - The Evolution Shortie RL Controlled-Rotation Dilator Sheath Set or Byrd Telescoping Stainless Steel Dilator Sheath Set should be used only to minimally enter the vessel.
 - Be aware that insertion of a stainless steel dilator sheath over a Byrd Dilator Sheath Set Polypropylene, the Evolution or Evolution RL Controlled-Rotation Dilator Sheath Set may damage the inner sheath.
 - When advancing sheaths including the Evolution or Evolution RL Controlled-Rotation Dilator Sheath Set, use proper sheath technique and maintain adequate tension on the catheter/lead (via a Locking Stylet or directly) to avoid damage to vessel walls.
 - If excessive scar tissue or calcification prevents safe advancement of sheaths, consider an alternate approach.
 - Excessive force with sheaths, including the Evolution or Evolution RL Controlled-Rotation Dilator Sheath Set used intravascularly may result in damage to the vascular system requiring surgical repair.
 - If the catheter/lead breaks, evaluate fragment; retrieve as indicated.
 - If hypotension develops, rapidly evaluate; treat as appropriate.
- Due to rapidly evolving catheter/lead technology, this device may not be suitable for the removal of all types of catheters/leads. If there are questions or concerns regarding compatibility of this device with particular catheters/leads, contact the catheter/lead manufacturer.

POTENTIAL ADVERSE EVENTS

Potential adverse events related to the procedure of intravascular extraction of catheters/leads include (listed in order of increasing potential effect):

- dislodging or damaging nontargeted catheter/lead
- chest wall hematoma
- thrombosis
- arrhythmias
- acute bacteremia
- acute hypotension
- pneumothorax
- stroke
- migrating fragment from catheter/object
- pulmonary embolism
- laceration or tearing of vascular structures or the myocardium
- hemopericardium/pericardial effusion
- cardiac tamponade
- hemothorax
- cardiac arrest
- death

INSTRUCTIONS FOR USE

Suggested Instructions for Use

Evolution RL Controlled-Rotation Dilator Sheath Set

WARNING: When using dilator sheaths or sheath sets do not insert more than one sheath set into a vein at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.

SUPERIOR APPROACH

1. Surgically expose the proximal end of the indwelling catheter/lead and remove the catheter/lead from its connections (if connected). Remove all suture and tie-down materials.

2. Cut off all proximal fittings, if present, using clippers or other cutters. It is important to cut the catheter/lead very close to the connector (but past any crimp joints) leaving as long a portion of the indwelling catheter/lead to work with as possible. Avoid closing off the interior lumen (or coil) of the catheter/lead when cutting it.

3. Consider passing a Locking Stylet through the inner lumen of the catheter/lead to stabilize the catheter/lead during dilation of surrounding tissues. Closely follow the "Suggested Instructions for Use" for the Locking Stylet to:

- A. Expose the inner coil of the catheter/lead
- B. Check patency of the coil lumen
- C. Determine the appropriate size of Locking Stylet based on the inner diameter of the catheter/lead
- D. Advance the Locking Stylet to the distal end of the catheter/lead
- E. Lock the Locking Stylet in place

4. Unless the catheter/lead insulation is damaged, degraded or too thin, tie a ligature or use a compression coil (Cook) at the proximal end of the catheter/lead, compressing the insulation against the coil and Locking Stylet to help prevent the coil and insulation from stretching. The ligature can be tied to the loop handle or to the suture tie loop.

NOTE: If a Locking Stylet has not been used, be aware that damage to the catheter/lead caused by pulling on it may prevent subsequent passage of a Locking Stylet through the lumen and/or make dilation of scar tissue more difficult.

5. For an active fixation catheter/lead, attempt to unscrew the catheter/lead by rotating the catheter/lead and Locking Stylet counterclockwise if appropriate.

6. Gently pull back on the catheter/lead to see if it is still engaged in tissue. If the catheter/lead is sufficiently loose in the tissue, gently pull on the Locking Stylet and catheter/lead to remove it.

NOTE: If removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.

7. If the catheter/lead is not removed from the vessel with gentle pulling then using dilator sheaths, including the Evolution or Evolution RL Controlled-Rotation Dilator Sheath Set may help separate the catheter/lead from tissue encapsulation.

8. With the inner Evolution RL Controlled-Rotation Dilator Sheath placed inside the appropriate outer sheath for telescopic action, insert the proximal free end of the catheter/lead with Locking Stylet in place, into the distal end of the inner Evolution RL Controlled-Rotation Dilator Sheath. Advance the catheter/lead with Locking Stylet, until it completely exits the opposite (proximal) end of the sheath set.

9. Apply adequate retracting pressure or tension on the catheter/lead via its Locking Stylet. **This is critical to safe passage of the sheath set over the catheter/lead. If tension is inadequate, the catheter/lead may buckle, precluding the sheath set from advancement along the appropriate path.**

10. With the catheter/lead in tension, and under fluoroscopic guidance, advance the inner sheath along the length of the catheter/lead and into the blood vessel. Rotate the inner Evolution® RL Controlled-Rotation Dilator Sheath, by squeezing and releasing the trigger activation handle.

For Bi-Rotational motion

Squeeze trigger to activate sheath rotation. A return of the trigger to its forward most position will mechanically change the rotational direction of the sheath to the direction opposite of the last rotation. Repeat as appropriate.

For Uni-Rotational motion

Squeeze trigger to activate sheath rotation and repeat without allowing trigger to return to its forward most position. Repeat as appropriate.

Advance the outer sheath over the inner sheath, keeping it within the vessel.

11. In a telescopic manner, alternately advance the inner, and then outer, sheaths while maintaining adequate tension on the catheter/lead and/or its wire guide or Locking Stylet, until the catheter/lead is free of tissue restriction. In general, the inner sheath should not be advanced more than 4 cm beyond the outer sheath.

NOTE: Always monitor fluoroscopically while advancing the sheaths to avoid shearing the catheter/lead or rupturing the vessel. Continue to maintain adequate tension on the catheter/lead and Locking Stylet. Adjust the angle of entry to keep the sheaths aligned with the catheter/lead in the vessel and to minimize sheath curvature during advancement. Rotating the inner Evolution RL Controlled-Rotation Dilator Sheath Set using the trigger activation handle and its outer sheath during advancement may facilitate progress through exuberant scar tissue.

NOTE: When the sheath set is tracking around a bend along the targeted lead/catheter, keep the outer sheath's beveled tip (if applicable) to the inside radius of the bend.

12. For cardiac leads, if the lead has not been freed by the time the sheaths near the myocardium, position the outer sheath so that the blunt end is directed toward the myocardium. Pull the inner sheath back several centimeters.

13. Apply countertraction with the outer sheath as follows:

Firmly hold the outer sheath one centimeter from the heart wall and steadily pull the Locking Stylet back; the pacing tip will be pulled into the sheath.

14. When the catheter/lead has become unentangled and freed from tissue attachment, carefully remove from the patient, the catheter/lead, Locking Stylet and Evolution or Evolution RL Controlled-Rotation Dilator Sheath Set.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Vandergrift, PA 15690 U.S.A.
COOK MEDICAL EUROPE LTD