- EN 2 **Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set**  
Instructions for Use
- CS 6 **Souprava dilatačního sheathu Evolution® Shortie RL s ovládáním rotace**  
Návod k použití
- DA 10 **Evolution® Shortie RL dilatatorsheathsæt med kontrolleret rotation**  
Brugsanvisning
- DE 14 **Evolution® Shortie RL Dilatatorschleusen-Set mit kontrollierter Rotation**  
Gebrauchsanweisung
- EL 18 **Σετ θηκαριών διαστολέα ελεγχόμενης περιστροφής**  
**Evolution® Shortie RL**  
Οδηγίες χρήσης
- ES 22 **Equipo de vainas dilatadoras de rotación controlada**  
**Evolution® Shortie RL**  
Instrucciones de uso
- FR 26 **Set de gaines de dilatation Evolution® Shortie RL à rotation contrôlée**  
Mode d'emploi
- HU 30 **Evolution® Shortie RL kontrollált forgatású dilatátorhüvely-készlet**  
Használati utasítás
- IT 34 **Set di guaine dilatatrici a rotazione controllata Evolution® Shortie RL**  
Istruzioni per l'uso
- LT 38 **„Evolution® Shortie RL“ kontroliuojamojo sukimo plėtiklio vamzdelių rinkinys**  
Naudojimo nurodymai
- NL 42 **Evolution® Shortie RL gecontroleerde-rotatiedilatatorsheath set**  
Gebruiksaanwijzing
- NO 46 **Evolution® Shortie RL dilatatorhysesett med kontrollert rotasjon**  
Bruksanvisning
- PL 50 **Zestaw koszulek rozszerzających Evolution® Shortie RL z regulacją obrotu**  
Instrukcja użycia
- PT 54 **Conjunto de bainha de dilatação com rotação controlada**  
**Evolution® Shortie RL**  
Instruções de utilização
- RU 58 **Набор расширителей-интродьюсеров Evolution® Shortie RL с контролируемым вращением**  
Инструкция по применению
- SV 62 **Evolution® Shortie RL dilatatorhylsset för styrd rotation**  
Bruksanvisning

**Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set**

Federal law restricts this device to sale by or on the order of a physician.

**DEVICE DESCRIPTION**

The Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set is comprised of an inner polymer sheath connected to a handle capable of mechanically rotating the sheath and an outer telescoping polymer sheath. The inner sheath has a stainless steel tip at its distal end.

**INTENDED USE**

Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set is intended for use in patients requiring the percutaneous dilation of tissue surrounding cardiac leads, indwelling catheters and foreign objects.

Sheath Size	REF Sheath Set
9 Fr	LR-EVN-SH-9.0-RL
11 Fr	LR-EVN-SH-11.0-RL

**Use With Other Devices**

The Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set may be used in conjunction with the following catheter/lead extraction devices from Cook:

- Locking Stylet
- Lead Extender

Be sure to closely follow the “Suggested Instructions for Use” for each device used.

**CONTRAINDICATIONS**

None known

**WARNINGS**

- **Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set should only be used to minimally enter the vessel.**
- **Do not attempt to negotiate the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set past a bend in the vessel as vessel wall or cardiac lead insulation damage may occur.**
- **When using dilator sheaths or sheath sets, do not insert more than one sheath set into a vein at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.**
- **When using a Locking Stylet:**
  - **Do not abandon a catheter/lead in a patient with a Locking Stylet still in place inside the catheter/lead. Severe vessel or endocardial wall damage may result from the stiffened catheter/lead or from fracture or migration of the abandoned stylet wire.**
  - **Do not apply weighted traction to an inserted Locking Stylet as myocardial avulsion, hypotension, or venous wall tearing may result.**
  - **Be aware that a lead that has a J-shape retention wire that occupies its inner lumen (rather than being outside the coil) may not be compatible with the Locking Stylet. Insertion of the Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.**
- **Weigh the relative risks and benefits of intravascular catheter/lead removal procedures in cases when:**
  - **The item to be removed is of a dangerous shape or configuration,**
  - **The likelihood of catheter/lead disintegration resulting in fragment embolism is high or vegetations are attached directly to the catheter/lead body.**
- **Catheter/lead removal devices should be used only at institutions with thoracic surgical capabilities.**
- **Catheter/lead removal devices should be used only by physicians knowledgeable in the techniques and devices for catheter/lead removal.**

**PRECAUTIONS**

- When activating the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set be careful not to cut or ensnarl lead suture-ties that are attached to the locking stylet targeted lead or to companion leads.
- Position the distal tip of the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set so that it does not engage companion leads during activation.
- Prior to the procedure, consider the size of the catheter/lead in relation to the size of the LEAD EXTRACTION™ devices to determine possible incompatibility.
- If selectively removing catheters/leads with the intent to leave one or more chronic catheters/leads implanted intact, the non-targeted catheters/leads must be subsequently tested to ensure that they were not damaged or dislodged during the extraction procedure.

As a result of the risk of complications, doctors highly experienced in this procedure have advised the following:

### **PREPARATIONS**

- Obtain a thorough patient history, including patient blood type. Appropriate blood products should be rapidly available.
- Ascertain the manufacturer, model number, and implant date of the catheter/lead to be removed.
- Perform radiographic/echocardiographic evaluation of catheter/lead condition, type, and position. Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, thoracotomy tray, and pericardiocentesis tray.
- Echocardiography and cardiothoracic surgery should be rapidly available.
- Prep/drape the patient's chest for possible thoracotomy; prep/drape the patient's groin for femoral approach or possible femoral approach.
- Establish back up pacing as necessary.
- Have available an extensive collection of sheaths, Lead Control Devices (Locking Stylet and Lead Extender), stylets to unscrew active fixation leads, snares, and accessory equipment.

### **PROCEDURE**

- Use fluoroscopic monitoring during ALL catheter/lead and sheath manipulations. Monitor ECG and arterial blood pressure continuously throughout the procedure and during recovery.
- If using sheaths or sheath sets, including the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set, the following precautions should be followed:
  - Prior to using sheaths including the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.
  - When advancing sheaths including the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set, use proper sheath technique and maintain adequate tension on the catheter/lead (via a Locking Stylet or directly) to prevent damage to vessel walls.
  - If excessive scar tissue or calcification prevents safe advancement of sheaths, consider an alternate approach.
  - Excessive force with sheaths used intravascularly may result in damage to the vascular system requiring surgical repair.
  - If the catheter/lead breaks, evaluate fragment; retrieve as indicated.
  - If hypotension develops, rapidly evaluate; treat as appropriate.
  - Due to rapidly evolving catheter/lead technology, this device may not be suitable for the removal of all types of catheters/leads. If there are questions or concerns regarding compatibility of this device with particular catheters/leads, contact the catheter/lead manufacturer.

### **POTENTIAL ADVERSE EVENTS**

Potential adverse events related to the procedure of intravascular extraction of catheters/leads include (listed in order of increasing potential effect):

dislodging or damaging nontargeted catheter/lead

chest wall hematoma

thrombosis

arrhythmias

acute bacteremia

acute hypotension

pneumothorax

stroke

migrating fragment from catheter/object

pulmonary embolism

laceration or tearing of vascular structures or the myocardium

hemopericardium/pericardial effusion

cardiac tamponade

hemothorax

cardiac arrest

death

## SUGGESTED INSTRUCTIONS FOR USE

### Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set

**WARNING: When using dilator sheaths or sheath sets, do not insert more than one sheath set into a vein at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.**

**WARNING: Do not attempt to negotiate the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set past a bend in the vessel as vessel wall or cardiac lead insulation damage may occur.**

1. Surgically expose the proximal end of the indwelling catheter/lead and remove the catheter/lead from its connections (if connected). Remove all suture and tie-down materials.
2. Cut off all proximal fittings, if present, using clippers or other cutters. It is important to cut the catheter/lead very close to the connector (but past any crimp joints) leaving as long a portion of the indwelling catheter/lead to work with as possible. Avoid closing off the interior lumen (or coil) of the catheter/lead when cutting it.
3. Consider passing a Locking Stylet through the inner lumen of the catheter/lead to stabilize the catheter/lead during dilation of surrounding tissues. Closely follow the Suggested Instructions for Use for the Locking Stylet to:
  - A. Expose the inner coil of the catheter/lead
  - B. Check patency of the coil lumen
  - C. Determine the appropriate size of Locking Stylet based on the inner diameter of the catheter/lead
  - D. Advance the Locking Stylet to the distal end of the catheter/lead
  - E. Lock the Locking Stylet in place
4. Unless the catheter/lead insulation is damaged, degraded or too thin, tie a ligature or use a compression coil (Cook) at the proximal end of the catheter/lead, compressing the insulation against the coil and Locking Stylet to help prevent the coil and insulation from stretching. The ligature can be tied to the loop handle or to the suture tie loop.

NOTE: If a Locking Stylet has not been used, be aware that damage to the catheter/lead caused by pulling on it may prevent subsequent passage of a Locking Stylet through the lumen and/or make dilation of scar tissue more difficult.
5. For an active fixation catheter/lead, attempt to unscrew the catheter/lead by rotating the catheter/lead and Locking Stylet counterclockwise if appropriate.
6. Gently pull back on the catheter/lead to see if it is still engaged in tissue. If the catheter/lead is sufficiently loose in the tissue, gently pull on the Locking Stylet and catheter/lead to remove it.

NOTE: If removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.
7. If the catheter/lead is not removed from the vessel with gentle pulling then using dilator sheaths, including the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set may help separate the catheter/lead from tissue encapsulation.
8. With the inner Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set placed inside the appropriate outer sheath for telescopic action, insert the proximal free end of the catheter/lead into the distal end of the inner Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set. Advance the catheter/lead until it completely exits the opposite (proximal) end of the sheath set.
9. Apply adequate retracting pressure or tension on the catheter/lead and Locking Stylet. **This is critical to safe passage of the sheath set over the catheter/lead. If tension is inadequate, the catheter/lead may buckle, precluding the sheath set from advancement along the appropriate path.**

10. With the catheter/lead in tension and under fluoroscopic guidance advance the inner and outer Evolution® Shortie RL Controlled-Rotation Dilator Sheaths along the catheter/lead activating the device as necessary until vessel entry is obtained.

To activate the inner Evolution® Shortie RL Controlled-Rotation Dilator Sheath, slowly squeeze and release the trigger activation handle (one squeeze every one to two seconds max):

**For Bi-Rotational motion**

Squeeze trigger to activate sheath rotation. A return of the trigger to its forward most position will mechanically change the rotational direction of the sheath to the direction opposite of the last rotation. Repeat as appropriate.

**For Uni-Rotational motion**

Squeeze trigger to activate sheath rotation and repeat without allowing trigger to return to its forward most position. Repeat as appropriate.

11. Immediately end the use of the device when the first of either, minimal entry into the vessel has been achieved or fluoroscopic monitoring confirms that a mid-clavicular location along the cardiac lead has been reached. Exchange the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set for the appropriately sized Evolution® RL Controlled-Rotation Dilator Sheath Set.

**WARNING: The Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set should be used only to minimally enter the vessel.**

NOTE: If using Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set when minimal entry into the vessel has been verified, exchange the Evolution® Shortie RL Controlled-Rotation Dilator Sheath Set for the Evolution® RL Controlled-Rotation Dilator Sheath Set as described in the "Suggested Instructions for Use".

**HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Vandergrift, PA 15690 U.S.A.
COOK MEDICAL EUROPE LTD