EN
3
CS
6
DA
10
DE
14
EL
18
ES
22
FI
26
FR
30
HU
34
IT
38
LT
42
LV
46
NL
50
NO
54
PL
58
PT
62
SK
66
SL
70
SV
74
TR
78**Pin Vise Polymer**

Instructions for Use

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Pin vise polymer

Brugsanvisning

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


3

PIN VISE POLYMER - LR-PVP001

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION**Performance Characteristics**

A chuck-like device that can be advanced over locking stylets. Tightening down of the rotational mechanism at its distal end causes the collet to close inward and grasp the proximal end of the outer diameter of a sheath for easier manipulation.

Device Compatibility**Compatible with the Pin Vise Polymer: Cook Byrd Dilator Sheath Sets Polypropylene**

Regular	LR-PPLBES-7.0	LR-PPLBES-8.5	LR-PPLBES-10.0	LR-PPLBES-11.5	LR-PPLBES-13.0
Extra Long	LR-PPLBES-7.0-XL	LR-PPLBES-8.5-XL	LR-PPLBES-10.0-XL	LR-PPLBES-11.5-XL	LR-PPLBES-13.0-XL

Be sure to closely follow the "Instructions for Use" for each device used. Use only with Cook devices.

Patient Population

The Pin Vise Polymer is used in adult patients who have cardiac implantable electronic devices (permanent pacemakers or implantable cardioverter-defibrillators) and indwelling cardiac leads that require lead extraction.

Intended User

The Pin Vise Polymer is intended for use by licensed physicians experienced in techniques used for lead removal and located at institutions with appropriate thoracic surgical capabilities. Physicians should possess expertise in heart disease as well as an understanding of unusual lead positions.

INTENDED USE

The Pin Vise Polymer is intended to aid the gripping of a dilator sheath by attaching to the proximal end of the sheath.

INDICATIONS FOR USE

The Pin Vise Polymer is designed as an auxiliary tool indicated for use in patients requiring percutaneous retrieval of cardiac leads.

Cardiac leads may need to be removed/extracted for a variety of medical indications, including:

- Cardiac implantable electronic device-related infection;
- Venous stenosis/occlusion (including thrombosis);
- Mechanical lead failure;
- Lead-related events (e.g., cardiac perforation, thromboembolism, superior vena cava syndrome, lead-lead interaction, lead misplacement);
- Lead advisory/recall (for potential lead malfunction);
- Removal of abandoned leads;
- Chronic pain after cardiac implantable electronic device implantation; or
- Cardiac implantable electronic device, cardiac lead, or system upgrade

CLINICAL BENEFITS

These devices provide assistance in the percutaneous retrieval of cardiac leads.

CONTRAINDICATIONS

None.

WARNINGS

This single use device is not designed for re-use. Risks of device re-use include exposure to contaminated device possibly resulting in systemic infection requiring medical intervention. Attempts to reprocess (re-sterilize) and/or to re-use may lead to contamination, device failure, and/or transmission of disease.

Be aware of the presence of vegetations larger than 2 centimeters in diameter. Before proceeding with the procedure, determine if the benefits of continuing outweigh the risks.

PRECAUTIONS

Prior to the procedure, ensure that the Pin Vise Polymer is securely fastened onto the Byrd Polypropylene Dilator Sheath.

General precautions for lead extraction procedures should be taken, including the following:

Prior to the procedure, consider the size of the lead in relation to the size of the lead extraction devices to determine possible incompatibility.

Ensure sufficient medical personnel and equipment are available.

Prep/drape the patient's chest for possible thoracotomy.

Before removing a lead with an extravascular portion, it is essential to carefully inspect the extravascular lead tract to ensure removal of all suture sleeves, sutures, tie-down materials, and Cook One-Tie Compression Coil, if used.

When using sheaths, do not insert sheaths over more than one lead at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.

When using a Liberator Beacon Tip Locking Stylet (referred to hereafter as Locking Stylet):

- Do not abandon a lead in a patient with a Locking Stylet still in place inside the lead. Severe vessel or endocardial wall damage may result from the stiffened lead or from fracture or migration of the abandoned stylet wire.

- Do not apply weighted traction to an inserted Locking Stylet as myocardial avulsion, hypotension, or venous wall tearing may result.

- Be aware that a lead that has a J-shape retention wire that occupies its inner lumen (rather than being outside the coil) may not be compatible with the Locking Stylet. Insertion of the Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.

Weigh the relative risks and benefits of intravascular lead removal procedures in cases when:

- the item to be removed is of a dangerous shape or configuration,
- the likelihood of lead disintegration resulting in fragment embolism is high,
- or vegetations are attached directly to the lead body.

- If selectively removing leads with the intent to leave one or more chronic leads implanted intact, the nontargeted leads must be subsequently tested to ensure that they were not damaged or dislodged during the extraction procedure.

The use of the Pin Vise Polymer in conjunction with the Byrd Dilator Sheath Set Polypropylene devices in pediatric and pregnant patients has not been well reported or studied. The safety and effectiveness of the Pin Vise Polymer in conjunction with the Byrd Dilator Sheath Set Polypropylene devices in these patients has not been established.

POTENTIAL ADVERSE EVENTS

Potential adverse events related to the procedure of intravascular extraction of leads include (listed in English alphabetical order):

- Acute bacteremia
- Acute hypotension
- Arrhythmias
- AV fistula that requires intervention
- Bleeding that requires blood transfusion
- Cardiac arrest
- Cardiac tamponade
- Cerebrovascular accident/event (e.g., stroke)
- Chest wall hematoma

Coronary sinus dissection

Death

Dislodging or damaging nontargeted lead

Hematoma that requires evacuation

Hemothorax that requires intervention

Laceration or tearing of vascular structures or the myocardium (e.g., cardiac avulsion, vascular laceration, need for vascular repair, coronary sinus dissection)

Migrated lead fragment without sequelae

Pericardial effusion that requires intervention, pericardial effusion without intervention

Pneumothorax that requires a chest tube

Pulmonary embolism

Thromboembolism that requires intervention

Tricuspid valve regurgitation (or damage)

Venous thrombosis that requires medical intervention

STERILIZATION AND HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. This product is single use and shall not be reused. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

STORAGE

Keep the device dry and away from sunlight.

INSPECTION OF DEVICE

Visually inspect the device thoroughly, including packaging, to verify that there is no damage prior to use. Visually inspect and confirm that the sterile barrier has not been compromised in any way. Do not use if the sterile packaging is damaged or unintentionally opened before use.

INSTRUCTIONS FOR USE

Pin Vise Polymer

1. Adjust opening of Pin Vise Polymer to accommodate Byrd Polypropylene Dilator Sheath
See illustration 1 on page 2
2. Place Byrd Polypropylene Dilator Sheath into Pin Vise Polymer
See illustration 2 on page 2
3. Tighten Pin Vise Polymer until it firmly grasps Byrd Polypropylene Dilator Sheath
See illustration 3 on page 2

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines. Risks of improper safe disposal include microbial and systemic infections, and airborne pathogens.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, and the competent authority and/or regulatory authority of the country where the device was used.

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Cook Vandergrift Inc.
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd