EN

AcuJect Flex[®] Injection Needle
Instructions for Use

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ENGLISH

INDICATIONS

Endoscopic injection needles are used for the endoscopic injection of legally marketed sclerotherapy agents and dyes into esophageal or colonic varices.

CONTRAINDICATIONS

Contraindications for this device are those applicable to injection therapy and include, but may not be limited to, those patients allergic to dyes, sclerotherapy or vasoconstriction agents and patients with lesions inappropriate for injection therapy with sclerotherapy or vasoconstriction agents.

WARNING AND PRECAUTIONS

The device is sharp. Use same care as for any other surgical device to avoid punctures or cuts to hands, gloves, etc.

When inserting or withdrawing device from endoscope, ensure the endoscope is not in a retroflexed position and the needle is fully retracted to avoid damaging the working channel.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure, infection, and/or transmission of disease.

Do not use if the sterile packaging is compromised, any abnormalities are found or if the instrument is not operating properly.

INSTRUCTIONS

1. Peel open the pouch and remove the device.
2. Visually inspect the device for kinks, loose or broken parts.
3. Uncoil the device and actuate the handle to ensure that the needle is operating properly.
4. Attach the syringe with injection media to the needle via the Luer-Lock and purge the needle of air.
5. Retract the needle before inserting the device into the endoscope.
6. Advance the device in small increments until the tip emerges from end of the endoscope. If the endoscope is equipped with an elevator, load the device past the elevator before raising fully.
7. Deploy the needle until you feel the handle "click" into place and proceed with the injection(s).
8. Following the injection(s), press the release button on the handle to fully retract the needle.
9. Upon completion of the procedure remove the device from the endoscope. If the endoscope has an elevator, lower the elevator before withdrawing the device.

NOTES

Use of this device restricted to a trained healthcare professional.

In the event of a malfunction, discontinue use and save for return and contact your local representative.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

DISPOSAL

After use, this product may be a potential biohazard. Dispose of it in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

STORAGE

Store sterile instruments at room temperature. Do not expose to solvents, radiation, UV light or any other substance or environment that could compromise the sterility of the unit.

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A symbol glossary can be found at cookmedical.com/symbol-glossary

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If symbol appears on product label, X = maximum insertion portion width

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If symbol appears on product label, X = minimum accessory channel

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EU REP

Advena Ltd.

Tower Business Centre, 2nd Flr,
Tower Street, Swatar, BKR 4013 Malta

Distributed By:
Cook Medical, LLC
400 N. Daniels Way
Bloomington, IN 47404 USA

cookmedical.com

Manufactured By:
Endo-Therapeutics, Inc.
15251 Roosevelt Blvd. Suite 204
Clearwater, FL 33760 USA

©2025 Endo-Therapeutics, Inc. All rights reserved. 190801 Rev. C Sept 2025

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ACRA Regulatory Services Pty Ltd