
EN

Caesar[®] Alligator Rat Tooth Grasping Forceps

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ENGLISH

INTENDED USE

This product is used in combination with endoscopes to remove tissue, stones, stents, or foreign bodies from the digestive tract.

INDICATIONS

This product is used in combination with endoscopes to remove tissue, stones, stents, or foreign bodies from the digestive tract.

CLINICAL BENEFIT

Removing tissue, stones, stents, or foreign bodies from the gastrointestinal tract.

CONTRAINDICATIONS

Contraindications include, but may not be limited to the following:

Grasping Forceps are contraindicated in patients undergoing anticoagulant therapy;

Grasping Forceps are contraindicated in children;

All contraindications for an endoscopic operation apply to the Grasping Forceps.

This device has not been designed, nor is it intended to be used with electrosurgical devices or their accessories.

POTENTIAL COMPLICATIONS

Procedural bleeding and delayed bleeding; perforation; infection; sepsis and abdominal pain.

WARNINGS

Do not use if the product is damaged, presents a defect or is not functional.

The product is only intended for adult populations.

MATERIAL

The distal jaw of the Grasping Forceps is made of medical grade stainless steel.

PACKAGING AND STERILIZATION

Grasping Forceps are individually packaged and sterilized with ethylene oxide.

Grasping Forceps cannot be re-sterilized.

Do not use if the sterile primary packaging is damaged.

Shelf life: 3 years (performance and sterility).

Do not use if expired.

FREQUENCY OF USE

Single use only. Use on different patients is forbidden.

PRECAUTIONS

Grasping Forceps must be used exclusively by physicians and nurses, under the supervision of a physician, who are properly trained in therapeutic endoscopy for retrieval of tissue, stones, stents and foreign bodies.

During use, gently advance the grasping forceps to the retrieval site through the endoscope channel. To avoid damage to the endoscope, do not apply force if resistance is encountered while advancing

Confirm that the endoscopic view is clear before use. Do not insert the instrument into the endoscope unless there is a clear endoscopic field of view. Insertion without a clear endoscopic field of view could cause patient injury, such as perforation, hemorrhage or mucous membrane damage. Damage to the endoscope and/or the instrument may also occur.

Maintain traction on the device during retrieval to prevent the retrieved object from becoming loose.

ENDOSCOPE COMPATIBILITY

Jaw diameter (mm)	Min working channel diameter of compatible endoscope (mm)
2.3	2.8

INSTRUCTIONS FOR USE

Before using the Grasping Forceps, perform a diagnostic endoscopy to visualize the area that requires sampling.

1. Open the primary packaging.
2. Visually inspect the Grasping Forceps to detect any abnormality with particular attention to kinks, bends and breaks.
3. Repeatedly activate the forceps handle by inserting a thumb into the designated ring and making the slip ring move back and forth to confirm the proper opening and closing of the forceps.
4. Before inserting the device into the endoscope, ensure that the Grasping Jaws are in the closed position. Place the Grasping Forceps into the endoscope's working channel, and gently advance to the target site under endoscopic view.
5. Once the object is visible through the endoscope, gently advance the device down the accessory channel to avoid bending the outer tube.
6. If resistance is encountered while inserting the Grasping Forceps into the endoscope, do not force the device. Instead, reduce the angulation of the endoscope or lower the forceps elevator in the endoscope until the device can pass smoothly.
7. Carefully observe the object for retrieval and move the Grasping Jaws to the targeted object. Optimal position can be achieved by grabbing the knob at the front end of the handle and rotating it to rotate the Grasping Jaws.
8. Open the Grasping Jaws and grasp the object gently while maintaining continuous traction. Do not use excessive force on the handle. Do not coil the outer tube outside of the endoscope.
9. Keep the distal end of the device endoscopically visible to ensure you do not lose sight of the object during removal.
10. Once Grasping Jaws grasp the object, slowly withdraw the object together with the endoscope from the patient. Do not retract the device and the object into the endoscope; this could cause damage to the endoscope and the device.
11. Once the endoscope and the device have been withdrawn from the patient, open the Grasping Jaws, and remove the object.

12. Remove the device from the endoscope with the Grasping Jaws in the closed position.
13. If removing additional objects from the patient, confirm that the Grasping Jaws open and close smoothly. Repeat steps 4 -13 above.

STORAGE AND TRANSPORTATION

Store at room temperature between 50 °F and 104 °F (10 °C and 40 °C), away from sunlight, non-corrosive gases, and in a well-ventilated area.

WORKING ENVIRONMENT

Surgery room of general hospitals, ambient temperature: 50-104 °F (10-40 °C). Relative humidity: 30%-85% and atmospheric pressure between 70-106 kpa.

PRODUCT DISPOSAL

Used devices must be disposed of according to institutional guidelines and local regulations.

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Cook Medical, LLC
Technology Co.,Ltd