EN
3**VITAL-PORT® Vascular Access Systems**

Suggested Instructions For Use

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













Rekommenderad bruksanvisning

10

The Vital-Port® Vascular Access System

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

INTRODUCTION

The Vital-Port Vascular Access System is a completely implantable subcutaneous system for patients who require repeated intermittent vascular access for simultaneous or separate infusion and/or withdrawal of fluids. Power Injection of contrast media may be performed on specified models when accessed with an infusion set indicated for power injection. The Vital-Port System consists of either one or two independent injection port chambers, each with a self-sealing septum accessed by percutaneous needle puncture using a non-coring needle, and a catheter.

The titanium or polysulfone port body has suture points that permit anchoring. (**Fig. 1**) In titanium models, these suture points are silicone-filled. A silicone or polyurethane catheter, with a catheter lock, is supplied either pre-attached or detached. A non-coring stainless steel needle is included with all models, and some systems also include a vascular access introducer set, tunneler and non-coring infusion set.

Intended Use

The **Vital-Port Vascular Access System** is indicated for patient therapy requiring repeated vascular access for simultaneous or separate injection or infusion therapy and/or blood sampling.

The **Vital-Port Vascular Access System Power Injectable Port** is intended for patient therapy requiring repeated vascular access for infusion therapy, power injected diagnostic techniques using contrast media, and blood infusion/withdrawal. For the power injection of contrast media, the maximum recommended infusion rate is 5mL/s using media with a maximum viscosity of 11.8 cP.

The following Vital-Port devices are indicated for power injection of contrast media when used with a power injectable infusion set. For further information on using these models for contrast power injection, please refer to the **Contrast Power Injection using Vital-Port System** section in these Instructions for Use.

- IP-7110, IP-S7010, IP-S7110, IP-S9010, IP-S9010P, IP-S9110, IP-S9110P

CONTRAINDICATIONS


The Vital-Port System should not be implanted in the presence of known or suspected infections, bacteremia, septicemia, or peritonitis, or in patients who have demonstrated prior intolerance to the materials of construction.

WARNINGS

- Introduction of the catheter into the subclavian vein (**Fig. 2, reference A**) using standard percutaneous techniques may subject the catheter to periodic compression forces within the narrow costoclavicular space between the clavicle (**Fig. 2, reference B**) and first rib (**Fig. 2, reference C**). Reported complications from repeated subclavian compression include catheter pinch-off syndrome, catheter fracture, and catheter shear followed by embolization of the distal portion. The likelihood of catheter occlusion and damage can be greatly reduced by choosing an access site lateral to the clavicle (preferably at or lateral to the midclavicular line) (**Fig. 2, reference D**) and prior to the vein entering the thorax at the costoclavicular space.
 - RECOMMENDED: Central Venous Catheter introduced at or lateral to midclavicular line (**Fig. 2 reference E**).
 - AVOIDED: Central Venous Catheter introduced medial to costoclavicular space (**Fig. 2, reference F**).
- Prevent air embolism by applying finger pressure over the opening of the introducer sheath.
- Failure to adequately anchor the port to the fascia increases the risk of catheter fracture and/or disconnection which could result in catheter migration.

Peripheral Implanted Port

- To avoid catheter damage to a vascular access system implanted peripherally in the upper arm, position the catheter so that it enters the venous system in the lower third of the upper arm and below the vein's passage through the deep or brachial fascia. Following insertion, confirm that the catheter follows a smooth course into the vein and that no abrupt directional changes are present.

- Abrupt directional changes may result in catheter fracture and subsequent possible embolization. A radiographic confirmation of the catheter insertion should be made to ensure that the catheter is not being pinched.

For Warnings related to models indicated for contrast power injection, please refer to the **Contrast Power Injection using Vital-Port System** section in these Instructions for Use.

Power Injection

- The indication for power injection of contrast media implies the port system's ability to withstand the procedure, but does not imply appropriateness of the procedure for a particular patient. A trained clinician is responsible to evaluate the health status of the patient, the integrity of the port system and the ability of staff to properly access the system for power injection.
- Verify that the port system is power injection compatible prior to using it for such a procedure. Check the patient's medical records and/or patient identification materials.
- If the catheter and lock are connected and then disconnected the catheter must be re-trimmed to ensure a secure re-connection.
- The Vital-Port Vascular Access Systems are only indicated for power injection when using the supplied port in conjunction with the catheter included in the tray.
- Do not place sutures on catheter, as catheter may be compromised.
- The Vital-Port Vascular Access Systems are only indicated for power injection when accessed with a power injectable infusion set.
- Failure to use an infusion set indicated for power injection may result in needle component failure.
- Do not exceed the maximum flow rate; doing so may result in septum leakage, catheter displacement and/or system failure.
- Do not exceed the maximum pressure; doing so may result in septum leakage, catheter displacement and/or system failure.
- Failure to ensure patency of the catheter lumen prior to injection may result in catheter failure.
- Do not use a port system that exhibits signs of compression or obstruction, as it may result in port system failure.
- Do not force injection or withdrawal of fluids through an impeded lumen. Notify attending physician immediately.
- Stop the injection immediately if local pain, swelling or sign of extravasation are noted.
- Flow rates were achieved using room temperature contrast at 11.8 centipoise. A change in temperature or viscosity of the contrast media used will result in a change in achievable flow rates.
- It is recommended that a 10 ml or larger syringe be used for injecting or infusing fluids through the system by hand.

PRECAUTIONS

General

- Consult the "Suggested Instructions for Use" before implanting or accessing the product.
- Use aseptic technique when implanting or accessing the device.
- Use talc-free gloves for handling the device.
- To avoid the potential problems associated with unnecessary or improperly maintained long-term device placement, consider removing the entire port system following cessation of therapy. A post-operative x-ray is suggested to verify system removal.

Implant Considerations

- Withdrawal or manipulation of distal spring coil portion of wire guide through needle tip may result in breakage.
- Avoid positioning the port septum directly below the incision line since the septum will be punctured repeatedly.
- Avoid contacting the catheter with sharp objects which may lead to catheter shearing.
- Avoid excessive clamp force which could damage the catheter. Clamp the catheter only on a section that will be cut off prior to implant. If the catheter was previously cut to length, keep the distal end elevated rather than clamping.
- Avoid trimming the catheter at an angle: an angled catheter tip can be sealed off against the side of a vessel.
- Improper catheter connection may result in catheter damage, leakage, or possible disconnection.
- Prior practice is recommended to ensure dexterity in connecting the catheter to the port. However, do **not** practice making a port-catheter connection with a catheter intended for implantation; doing so may result in catheter damage.

- Verify that the catheter lock and reinforcing sleeve are retrieved, along with the port body and total length of the catheter, upon removal of the Vital-Port System. A post operative x-ray is suggested for verification.

Accessing the Port

- Follow all manufacturer guidelines for infusates.
- Do not exert excessive force during needle insertion. The needle point may become damaged and could tear the septum when the needle is removed.
- Do not exceed 40 psi pressure when administering fluid through the system. Use of a 10 mL or larger syringe is recommended for all injection and flushing procedures. Pressures in excess of 40 psi can be generated with most syringes; the smaller the syringe, the higher the pressure that can be generated. For Vital-Port system models that are labeled as a POWER INJECTABLE PORT, refer to the **Contrast Power Injection using Vital-Port System** section in these Instructions for Use.
- Inspect and replace the infusion set and sterile dressing according to CDC guidelines, or at intervals according to hospital policy or current practice.

POTENTIAL ADVERSE EVENTS

Use of the Vital-Port System involves potential risks normally associated with the insertion or use of any implanted device or indwelling catheter. These include, but not limited to: air embolism; bacteremia, infection; sepsis; brachial plexus injury; cardiac arrhythmia, puncture, or tamponade; catheter disconnection, fracture, fragmentation, malposition, pinch-off, rupture, or shearing; embolization, extravasation, or formation of calcium-like deposits; fibrin sheath formation; hemothorax; hematoma; port/catheter erosion, migration, occlusion, rejection, or rotation; pneumothorax; thoracic duct injury; thrombosis; thromboembolism; vessel erosion, laceration, perforation, or trauma.

PATIENT COUNSELING INFORMATION

Each Vital-Port System is supplied with a patient information card, containing the port model and lot number as well as manufacturer contact information. Please complete the requested information on the card and instruct the patient to carry the card with him or her.

A "Patient Information" booklet is also supplied with each port or available separately.

INSTRUCTIONS FOR USE: PLACEMENT

System Preparation

Prior to placement, all air must be flushed from the chamber(s) and catheter lumen(s) of the Vital-Port system using a non-coring needle, a 10mL or larger syringe, and heparinized saline (100 IU/mL).

CAUTION: Use aseptic techniques and extreme care when handling the catheter. Avoid both contacting the catheter with sharp objects (which may lead to catheter shearing) and excessive clamp force (which could damage the catheter). Clamp the catheter only on a section that will be cut off prior to implant. If the catheter was previously cut to length, keep the distal end elevated rather than clamping. When trimming catheter to size, perform a perpendicular cut in one motion.

Port Placement Considerations

- The injection port reservoir should be placed where it will be supported by an underlying bony structure. The location should be convenient and comfortable for the patient. Patient mobility should also be considered.
- The incision should not be significantly larger than the length of the port base. Avoid positioning the port septum directly below the incision line since the septum will be punctured repeatedly.
- The port should not be placed too deeply to be easily located, or too superficially, causing erosion due to pressure necrosis. A depth between 0.5cm and 2cm under the skin surface is recommended. The port should be easily located by palpation for access.

Catheter Placement Considerations

- The catheter tip should be placed in an area of high blood flow (usually the SVC) **Figure 3**.
- To allow for body movement, the catheter should not be implanted with tension. Also avoid excessive catheter slack in the port pocket to prevent possible kinking.

WARNING

• Introduction of the catheter into the subclavian vein (**Fig. 2, reference A**) using standard percutaneous techniques may subject the catheter to periodic compression forces within the narrow costoclavicular space between the clavicle (**Fig. 2, reference B**) and first rib (**Fig. 2, reference C**). Reported complications from repeated subclavian compression include catheter pinch-off syndrome, catheter fracture, and catheter shear followed by embolization of the distal portion. The likelihood of catheter occlusion and damage can be greatly reduced by choosing an access site lateral to the clavicle (preferably at or lateral to the midclavicular line) (**Fig. 2, reference D**) and prior to the vein entering the thorax at the costoclavicular space.

RECOMMENDED: Central Venous Catheter introduced at or lateral to mid-clavicular line (**Fig. 2, reference E**).

AVOIDED: Central Venous Catheter introduced medial to costoclavicular space (**Fig. 2, reference F**).

For additional information and references on this subject, a booklet entitled **Understanding and Avoiding Subclavian Compression of Central Venous Catheters** is available from Cook.

Catheter Placement

For percutaneous access:

- Choose an appropriately sized introduction system and attach the introducer needle to a syringe.
- Introduce the needle into the desired vessel, gently aspirating during introduction.
- Remove the syringe from the needle, being sure to cover the needle opening to prevent air embolism.
- Insert the wire guide into the needle using wire guide straightener if included.
- Advance the wire guide to the appropriate position. Confirm position using radiographic technology.
- Remove the needle and straightener.
- Advance the introducer dilator/sheath over the wire guide.
- Remove the dilator. **WARNING: Prevent air embolism by applying finger pressure over the opening of the introducer sheath.**
- Insert the catheter into the sheath. Position the distal end of the catheter at the desired location using appropriate imaging technologies (**Fig. 3**)
- Peel away the sheath while withdrawing from the vessel.
- Verify correct catheter tip placement using fluoroscopy or other appropriate imaging technology.

For cutdown access:

- Make a small incision to expose the entry vessel of choice.
- Isolate and stabilize the vessel. Perform vessel incision.
- Introduce the catheter and advance the tip to the desired location. (**Fig. 3**)
- Verify correct catheter tip position using fluoroscopy or other appropriate imaging technology.

Pre-Attached Catheter Models

- Prior to catheter placement select the site for the port pocket.
- Make an incision at the edge of the desired anatomic location for the port pocket.
- Create the subcutaneous pocket using blunt dissection
- Measure the catheter length. Provide enough slack to allow for body movement.
- Trim the excess catheter by cutting the distal end squarely.
- Verify that the port chamber(s) and catheter are filled with heparinized saline (100 IU/mL).
- Gain access to the vascular system and place the catheter using either percutaneous or cutdown approach as described above. Use a standard tunneling device when required.
- Secure port body in the subcutaneous pocket.

Detached Catheter Models

- Verify that the catheter is filled with heparinized saline (100 IU/mL).
- Gain access to the vascular system using standard percutaneous techniques or surgical cutdown.
- Make an incision at the edge of the desired anatomic location for the port pocket.
- Create the subcutaneous pocket using blunt dissection
- Position the distal end of the catheter in the vessel at the desired location. Advance the catheter from the selected access vessel back to the port pocket. Use a standard tunneling device when required.

- After placement, confirm correct position of the catheter using radiographic techniques. **NOTE:** To allow for body movement, the catheter should not be implanted with tension; however, be sure to avoid excessive catheter slack in the port pocket to prevent possible kinking.
- Trim the excess catheter by cutting the proximal end squarely. **CAUTION:** Cutting squarely will aid in the proper sealing of the locking mechanism.
- Verify that the port chamber(s) and outlet tube(s) are filled with heparinized saline (100 IU/mL), then attach the catheter to the port outlet tube(s) as indicated in the SYSTEM ASSEMBLY section below.
- Secure port body in the subcutaneous pocket.

Mini Titanium 5116W models catheter placement

The Vital-Port Mini model 5116W is supplied with a detached catheter having a Slip-Coat™ Wire Guide Obturator pre-inserted for facilitating peripheral placement.

The Vital-Port Mini Titanium 5116W models should be used only after careful consideration and with extreme caution for patients who have had extensive trauma to the veins or have anatomical irregularities which may hamper proper placement of the system.

- Connect a syringe with heparinized saline (100 IU/mL) to the stopcock and flush the catheter. Close the stopcock by rotating the handle 90° (**Fig. 4, reference A**).
- The catheter should be placed via a percutaneous puncture technique or a cutdown procedure.
- Position the distal end of the catheter in the vessel at the desired location.
- With the catheter in the desired position, remove the obturator from the catheter by removing the stopcock (**Fig. 5, reference A**) and gently withdrawing the obturator (**Fig. 5, reference B**).
- Next, flush the catheter with heparinized saline (100 IU/mL) via the Luer connector on the end of the catheter.
- Clamp the catheter closed at the proximal end to be connected to the port, clamping the catheter only on the section that will be cut off prior to implant. Cut the Luer connector off of the catheter.

SYSTEM ASSEMBLY

For single-chamber models:

- Position the catheter lock over the proximal end of the catheter (**Fig. 6, reference A and Fig. 7, reference A**) and push the catheter over the port outlet (**Fig. 6, reference B and Fig. 7, reference B**). To facilitate advancement over the ring, rotate the catheter while pushing forward. **NOTE:** On Standard size models, the end of the catheter should be advanced halfway (**Fig. 6, reference D**) between the ring (**Fig. 6, reference C**) and shoulder of the connecting tube (**Fig. 6, reference E**). On petite and mini models, the end of the catheter should be advanced halfway (**Fig. 7, reference D**) between the ring (**Fig. 7, reference C**) and port body (**Fig. 7, reference E**).
- Slide the catheter lock over the port outlet tube and catheter until it is engaged in the recess of the port body.
- When properly positioned, the end of the catheter lock should not be visible. The end of the catheter should be visible through the clear portion of the catheter lock (**Fig. 8**).

For dual-chamber models:

- Position the catheter lock over the proximal end of the catheter and push the catheter over the port outlet tubes (**Fig. 9, reference A**), aligning the catheter lumens with the port outlet tubes.
- The end of the catheter should be positioned against the port body, completely covering the port outlet tubes. This can be visualized through the window in the bottom of the port body (**Fig. 9, reference B**).
- Slide the catheter lock over the port outlet tubes and catheter until engaging the threads in the port body. Turn the catheter lock clockwise until rotation is stopped by the locator pin. When properly positioned, the wings of the catheter lock will be parallel with the base of the port.
- The end of the catheter should now appear as a small "O" ring or donut within the port body window (**Fig. 9, reference C**).

CAUTION: Improper catheter connection may result in catheter damage, leakage, or possible disconnection. Prior practice is recommended to ensure dexterity in connection the catheter to the port. However, do not practice with a catheter intended for implantation; doing so may result in damage to the catheter.

Flow Verification

Before closing the port pocket, verify that there is unimpeded fluid flow. To confirm two-way fluid flow, use aspiration followed by a flush with 20 mL of saline to cleanse the catheter lumen(s) and port chamber(s) of blood. Flush each chamber and corresponding catheter lumen with 5 mL of heparinized saline (100 IU/mL; 3 mL for Mini systems). For dual models, check both chambers and corresponding catheter lumens. Use a non-coring needle and a 10mL or larger syringe with saline to confirm that the flow is not obstructed and that there are no leaks. After flow verification, establish a heparin lock with 5mL of heparinized saline (100 IU/mL; 3 mL for Mini systems). Continue to inject fluid and depress the port with the fingers while withdrawing the needle to avoid reflux.

Port Suturing and Site Closure

Anchor the port to the fascia with nonabsorbable sutures. Three sutures, with at least one at each end of the port, are recommended. Close the site, making sure that the port does not lie directly beneath the incision. Access the port, to ensure a proper heparinized saline lock, before sending the patient to the recovery area.

WARNING: Failure to adequately anchor the port to the fascia increases the risk of catheter fracture and/or disconnection which could result in catheter migration.

INSTRUCTIONS FOR USE: INJECTION, INFUSION AND BLOOD SAMPLING

General Considerations

- Always use aseptic technique before accessing the system.
- Locate the port by palpation.
- Observe for wound hematoma or accumulation of serous fluid around the implant site.
- Use 10 mL syringes or larger for all flushing or injection procedures.
- Use only 19 to 22 gage, single-use, non-coring needles. Non-coring needles preserve the life of the self-sealing septum.

Needle Insertion

Needles should be inserted perpendicular to the implanted port after locating the port septum.

1. Puncture the skin immediately over the septum and advance the needle point slowly through the septum until it makes contact with the bottom of the port chamber (**Fig. 10**). **NOTE:** Incomplete needle insertion may result in needle obstruction by the septum. Do not rock or tilt the needle.
2. Stabilize the needle in the port septum and cover with sterile dressing. **NOTE:** Never leave the system open to air.

Saline Flushes

- Prior to drug administration, flush the system with saline solution.
- If more than one drug is administered during the treatment, flush the system with saline solution between drugs.
- After completing patient treatment, always flush the system with saline to cleanse the catheter and port chamber.

Heparin Lock

To keep the Vital-Port system patent, establish a heparin lock by flushing each port chamber and catheter lumen with 5 mL heparinized saline (100 IU/mL; 3 mL for Mini systems) after each patient treatment and at least once every 4 weeks.

Injection/Infusion

1. Before accessing the port chamber, connect the syringe or IV set to a non-coring needle, and fill the infusion set with saline to remove air.
2. Access the port as previously described and begin infusion or injection. Examine the puncture site closely for any signs of subcutaneous irritation.
3. Flush the system as previously described.

Blood Sampling

1. Flush the port with saline to confirm system patency.
2. Withdraw 3-4 mL blood sample and discard.
3. Withdraw the required blood volume.
4. Immediately follow the sampling with a 20 mL saline flush to cleanse the port chamber and catheter.
5. If no further treatments are required, establish a heparin lock as previously described.

Contrast Power Injection using Vital-Port System

Only the following Vital-Port models are indicated for power injection of contrast media:

IP-7110, IP-S7010, IP-S7110, IP-S9010, IP-S9010P, IP-S9110 and IP-S9110P.

1. Access the port with a power-injectable infusion set. The needle should be inserted perpendicular to the implanted port and should be advanced slowly until it makes contact with the bottom of the port chamber.

WARNING: The Vital-Port Vascular Access Systems are only indicated for power injection when accessed with a power-injectable infusion set. Failure to use a power-injectable infusion set may result in needle component failure.

2. Attach a 10 mL or larger syringe filled with sterile, normal saline.
3. Place the patient in the position that he or she will be in during the power injection.
4. Aspirate for adequate blood return and vigorously flush the system with at least 10 mL of sterile, normal saline.

WARNING: Failure to ensure patency of the catheter lumen prior to injection may result in catheter failure.

- During the saline flush, observe for swelling or signs of extravasation. Ask the patient if there is any burning, pain or discomfort.

WARNING: Stop the injection immediately if local pain, swelling or sign of extravasation are noted.

- Detach the syringe.
- Warm contrast media to body temperature prior to injection as per the contrast manufacturer's instructions for use.
- Attach the power injection device to the power-injectable infusion set and ensure that the connection is secure.
- Check the indicated flow rate for the infusion set and confirm power injection device settings for flow and pressure.

Infusion Set Gauge Size	19 Ga	20 Ga
Max Flow Rate Setting	5 mL/sec	5 mL/sec
Max Pressure Setting	325 psi	

WARNING: Do not exceed the maximum flow rate or pressure; doing so may result in septum leakage, catheter displacement and/or system failure.

- Instruct the patient to notify staff immediately if there is any pain or a change in feeling occurs during the injection.
- Complete the power injection study.
- Disconnect the power injection device.
- Flush the Vital-Port System with 10 mL of sterile saline using a 10 mL or larger syringe.
- Establish a heparin lock, clamping while injecting the last of the flush solution.
- Remove the power injectable infusion set and follow the Instructions for Use for any special needle disposal considerations.

Power Injection Information

Catheter and Size	Maximum Indicated CT Flow Rate	Average Pressure in the Port Reservoir	Average Port-Catheter Burst Pressure
9.5 Fr Silicone	5 mL/s	42 psi	169 psi
7.5 Fr Polyurethane	5 mL/s	54 psi	200 psi
7.5 Fr Silicone	5 mL/s	52 psi	171 psi
Note: CT injector pressure should be set at a maximum of 325 psi. Flow rates less than 5 mL/s and/or lower viscosity contrast will generate lower pressures in the port and catheter. Catheter length of at least 70 cm used for all testing. The Vital-Port System testing included 40 power injection cycles using 11.8 cP viscosity media and a flow rate of 5 mL/s.			

1cP = 1mPa·s

1psi = 6.894kPa

DISCONTINUING USE OF THE VITAL-PORT SYSTEM

In vitro testing has shown sealing integrity after 1200 and 2500 punctures for 20 and 22 ga, respectively (400 and 1000 punctures for the Mini). Septum lifetime may vary, depending on clinical use.

It is suggested that the clinician consider removing the Vital-Port System upon completion of treatment.

Upon removal of the Vital-Port System, verify that the catheter lock and reinforcing sleeve are retrieved along with the port body and total length of the catheter. A post operative x-ray is suggested for verification.

Ports and catheters that are removed because the clinician believes there has been a system malfunction should be returned to Cook for evaluation. Please contact customer service for a returned goods authorization number and instructions.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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EN MR-Conditional

Testing has demonstrated that the Vital-Port® Vascular Access Systems are MR-conditional. A patient with the implant can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less - Spatial gradient magnetic field of 720 Gauss/cm or less
- Maximum MR system reported whole body-average specific absorption rate (SAR) of 3-W/kg for 15 minutes of scanning

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