EN
2
CS
6
DA
10
DE
14
EL
19
ES
24
FR
29
IT
34
LT
39
NL
43
PL
48
PT
53
RU
58
SV
63

Liberator® Beacon® Tip Locking Stylet

Instructions for Use

Uzamykací stilet Liberator® s hrotem Beacon®

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Liberator® låsestilet med Beacon® spids

Brugsanvisning

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Bruksanvisning

Liberator® Beacon® Tip Locking Stylet

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Liberator® Beacon® Tip Locking Stylet is comprised of a handle and a stylet wire having a fixation mechanism attached to the distal end.

INDICATIONS

The Liberator Beacon Tip Locking Stylet is intended for use in patients requiring the percutaneous removal of cardiac leads, indwelling catheters and foreign objects, with a central lumen.

Use With Other Devices

The Liberator Beacon Tip Locking Stylet may be used in conjunction with the following catheter/lead extraction devices from COOK Incorporated:

- Bulldog™ Lead Extender
- Byrd Telescoping Stainless Steel Dilator Sheath Set
- Byrd Dilator Sheath Set - Polypropylene or Teflon
- Evolution® Mechanical Dilator Sheath Set
- Evolution® Shortie Mechanical Dilator Sheath Set
- Byrd Workstation™
- Needle's Eye Snare®

Be sure to closely follow the "Suggested Instructions for Use" for each device used.

CONTRAINDICATIONS

None known

WARNINGS

When using sheaths, do not insert sheaths over more than one lead at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.

When using the Liberator Beacon Tip Locking Stylet:

Do not abandon a catheter/lead in a patient with the Liberator Beacon Tip Locking Stylet still in place inside the catheter/lead. Severe vessel or endocardial wall damage may result from the stiffened catheter/lead or from fracture or migration of the abandoned stylet wire.

Do not apply weighted traction to an inserted Liberator Beacon Tip Locking Stylet as myocardial avulsion, hypotension, or venous wall tearing may result.

Be aware that a lead that has a **J-shape retention wire** that **occupies its inner lumen** (rather than being outside the coil) may not be compatible with the Liberator Beacon Tip Locking Stylet. Insertion of the Liberator Beacon Tip Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.

Weigh the relative risks and benefits of intravascular catheter/lead removal procedures in cases when:

- the item to be removed is of a dangerous shape or configuration,
- the likelihood of catheter/lead disintegration resulting in fragment embolism is high,
- or vegetations are attached directly to the catheter/lead body.

Catheter/lead removal devices should be used only at institutions with thoracic surgical capabilities.

Catheter/lead removal devices should be used only by physicians knowledgeable in the techniques and devices for catheter/lead removal.

PRECAUTIONS

Prior to the procedure, consider the size of the catheter/lead in relation to the size of the LEAD EXTRACTION™ devices to determine possible incompatibility.

If selectively removing catheters/leads with the intent to leave one or more chronic catheters/leads implanted intact, the nontargeted catheters/leads must be subsequently tested to ensure that they were not damaged or dislodged during the extraction procedure.

As a result of the risk of complications, doctors highly experienced in this procedure have advised the following precautions:

Preparations

Obtain a thorough patient history, including patient blood type. Appropriate blood products should be rapidly available.

Ascertain the manufacturer, model number, and implant date of the catheter/lead to be removed.

Perform radiographic/echocardiographic evaluation of catheter/lead condition, type, and position.

Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, thoracotomy tray, and pericardiocentesis tray. Echocardiography and cardiothoracic surgery should be rapidly available.

Prep/drape the patient's chest for possible thoracotomy; prep/drape the patient's groin for femoral approach or possible femoral approach.

Establish back-up pacing as necessary.

Have available an extensive collection of sheaths, Locking Stylets, stylets to unscrew active fixation leads, snares, and accessory equipment.

Procedure

Use fluoroscopic monitoring during ALL catheter/lead and sheath manipulations. Monitor ECG and arterial blood pressure continuously throughout the procedure and during recovery.

If using sheaths or sheath sets, the following precautions should be followed:

Prior to using sheaths, it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

The Evolution Shortie or Byrd Telescoping Stainless Steel Dilator Sheath Set should only be used to minimally enter the vessel.

When advancing sheaths, use proper sheath technique and maintain adequate tension on the catheter/lead (via a Locking Stylet or directly) to prevent damage to vessel walls.

If excessive scar tissue or calcification prevents safe advancement of dilator sheaths, consider an alternate approach.

Excessive force with sheaths used intravascularly may result in damage to the vascular system requiring surgical repair.

If the catheter/lead breaks, evaluate fragment; retrieve as indicated.

If hypotension develops, rapidly evaluate; treat as appropriate.

Due to rapidly evolving catheter/lead technology, this device may not be suitable for the removal of all types of catheters/leads. If there are questions or concerns regarding compatibility of this device with particular catheters/leads, contact the catheter/lead manufacturer.

ADVERSE EVENTS

Potential Adverse Events

Potential adverse events related to the procedure of intravascular extraction of catheters/leads include (listed in order of increasing potential effect):

dislodging or damaging nontargeted catheter/lead

thrombosis

acute bacteremia

pneumothorax

migrating fragment from catheter/object

laceration or tearing of vascular structures or the myocardium

cardiac tamponade

death

chest wall hematoma

arrhythmias

acute hypotension

stroke

pulmonary embolism

hemopericardium

hemothorax

INSTRUCTIONS FOR USE

Suggested Instructions for Use: Liberator Beacon Tip Locking Stylet

WARNING: Do not abandon a catheter/lead in a patient with the Liberator Beacon Tip Locking Stylet still in place inside the catheter/lead. Severe vessel or endocardial wall damage may result from the stiffened catheter/lead or from fracture or migration of the abandoned stylet wire.

WARNING: Be aware that a lead that has a **J-shape retention wire** that **occupies its inner lumen** (rather than being outside the coil) may not be compatible with the Liberator Beacon Tip Locking Stylet. Insertion of the Liberator Beacon Tip Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.

WARNING: Do not apply weighted traction to an inserted Liberator Beacon Tip Locking Stylet as myocardial avulsion, hypotension, or venous wall tearing may result.

NOTE: For targeted active fixation leads with a retractable distal helix, consider retracting the helix (by following lead manufacturers instruction) before prepping the lead.

1. After the proximal end of the catheter/lead has been exposed, sutures and tie-down materials removed, and the proximal fittings (if present) cut off, expose the inner coil of the catheter/lead.

NOTE: For a unipolar pacing lead, if the lead is long enough, use a scalpel to circumscribe the insulation two centimeters from the cut proximal end of the lead (do this carefully to avoid deforming the coil), then pull the insulation off to expose the inner coil. For a bipolar lead with coaxial coils, remove the outer coil and inner insulation to expose the inner coil and prevent it from being pushed deeply into the outer coil.

NOTE: Examine the lumen to be sure the interior coil is not flattened and there are no burrs that would inhibit passage of the Liberator Beacon Tip Locking Stylet into the lumen.

NOTE: Because implant procedures vary dramatically, catheters/leads may be secured in a variety of ways at the time of implant. Carefully inspect the site to ensure all suture sleeves, sutures, and tie down materials which may be hidden or blind have been removed.

NOTE: It is important to cut the catheter/lead very close to the connector (but past any crimp joints) leaving as long a portion of the indwelling catheter/lead to work with as possible. Avoid closing off the interior lumen (or coil) of the catheter/lead when cutting it.

NOTE: Do not pull on the catheter/lead because it may stretch, distort, or break, making subsequent removal more difficult. Damage to the lead may prevent passage of a Liberator Beacon Tip Locking Stylet through the lumen and/or make dilation of scar tissue more difficult.

2. Check patency of the coil lumen. Pass a standard wire guide or stylet wire through the inner lumen of the catheter/lead to clear any blood, coagulation, or blockage in the lumen prior to insertion of the Liberator Beacon Tip Locking Stylet. After confirmation of a clear lumen, remove the standard stylet wire.

NOTE: It may be helpful to mark the depth of penetration of a normal stylet wire into the catheter/lead by placing a mosquito clamp at that point where the stylet wire exits the cut end of the catheter/lead. Knowing the depth of penetration will be useful later in monitoring progress during insertion of the Liberator Beacon Tip Locking Stylet.

3. Advance the Liberator Beacon Tip Locking Stylet to the distal end of the catheter/lead, with fluoroscopic monitoring.

4. To lock the Liberator Beacon Tip Locking Stylet in place:

Figure 1

After the Liberator Beacon Tip Locking Stylet is advanced to the tip of the lead, unlatch the tubular handle of the locking mechanism by removing the latchpin (Figure 1).

Figure 2

Lock the Liberator Beacon Tip Locking Stylet in the lead by firmly advancing the tubular handle forward while holding integral wire extender in place (Figure 2)

5. Unless the catheter/lead insulation is damaged, degraded or too thin, tie a ligature at the proximal end of the catheter/lead, compressing the insulation against the coil and the Liberator Beacon Tip Locking Stylet to help prevent the coil and insulation from stretching. The ligature can be tied to the loop handle.

6. For an active fixation catheter/lead with an exposed distal helix, attempt to unscrew the catheter/lead by rotating the catheter/lead (and Locking Stylet if used) counterclockwise.

7. Gently pull back on the catheter/lead to see if it is still engaged in tissue. If the catheter/lead is sufficiently loose in the tissue, gently pull on the Liberator Beacon Tip Locking Stylet and catheter/lead to remove it.

NOTE: If a Locking Stylet has not been used, be aware that damage to the catheter/lead caused by pulling on it may prevent subsequent passage of a Locking Stylet through the lumen and/or make dilation of scar tissue more difficult.

NOTE: If removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.

8. If the catheter/lead is not removed from the vessel with gentle pulling, dilator sheaths or other retrieval devices may help separate the catheter/lead from any tissue encapsulation as follows:

Advance dilator sheaths (or other retrieval device) over proximal end of the Liberator Beacon Tip Locking Stylet.

9. If for some reason removal is unsuccessful or becomes medically contraindicated, removal of the Locking Stylet may be facilitated as described below.

Pull on the tubular handle to retract it toward the loop handle. If the Liberator Beacon Tip Locking Stylet is freed, remove it.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.