EN
2

Byrd Workstation™

Instructions for Use

CS
6

Byrd Workstation™

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DA
10

Byrd Workstation™

Brugsanvisning

DE
14

Byrd Workstation™

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EL
19

Byrd Workstation™

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Bruksanvisning

Byrd Workstation™

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Byrd Workstation™ consists of a 16 Fr plastic introducer sheath, a retrieval basket within a 12 Fr plastic sheath, and a 0.045" diameter 125 cm long tip deflecting wire guide.

INDICATIONS

The Byrd Workstation™ is intended for use in patients requiring the percutaneous retrieval of cardiac leads, indwelling catheters, fragments of catheter tubing or wire guides, and other foreign objects.

Use With Other Devices

The Byrd Workstation™ may be used in conjunction with the following catheter/lead extraction devices from Cook Incorporated:

- Locking Stylet
- Byrd Telescoping Stainless Steel Dilator Sheath Set
- Byrd Dilator Sheath Set - Polypropylene or Teflon
- Metal-Reinforced Flexible Dilator Sheath
- Needle's Eye® Snare

Be sure to closely follow the "Suggested Instructions for Use" for each device used.

CONTRAINDICATIONS

None known

WARNINGS

Weigh the relative risks and benefits of intravascular catheter/lead removal procedures in cases when:

- the item to be removed is of a dangerous shape or configuration, the likelihood of catheter/lead disintegration resulting in fragment embolism is high,
- or vegetations are attached directly to the catheter/lead body.

Catheter/lead removal devices should be used only at institutions with thoracic surgical capabilities.

Catheter/lead removal devices should be used only by physicians knowledgeable in the techniques and devices for catheter/lead removal.

PRECAUTIONS

Prior to the procedure, consider the size of the catheter/lead in relation to the size of the LEAD EXTRACTION™ devices to determine possible incompatibility. If selectively removing catheters/leads with the intent to leave one or more chronic catheters/leads implanted intact, the nontargeted catheters/leads must be subsequently tested to ensure that they were not damaged or dislodged during the extraction procedure.

Competency with the Seldinger technique is necessary. Care should be taken to avoid the artery and to enter the vein cleanly and at a shallow angle.

Do not use excessive force to grasp and pull large portions of doubled-over lead/catheter into the 12 Fr or 16 Fr sheath, as it could become lodged, preventing further movement of the sheaths or snare(s).

As a result of the risk of complications, doctors highly experienced in this procedure have advised the following precautions:

Preparations

Obtain a thorough patient history, including patient blood type.

Appropriate blood products should be rapidly available.

Ascertain the manufacturer, model number, and implant date of the catheter/lead to be removed. Perform radiographic/echocardiographic evaluation of catheter/lead condition, type, and position.

Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, thoracotomy tray, and pericardiocentesis tray. Echocardiography and cardiothoracic surgery should be rapidly available.

Prep/drape the patient's chest for possible thoracotomy; prep/drape the patient's groin for femoral approach.

Establish back-up pacing as necessary.

Have available an extensive collection of sheaths, Locking Stylets, stylets to unscrew active fixation leads, snares, and accessory equipment.

Procedure

Use fluoroscopic monitoring during ALL sheath, catheter/lead, or snare manipulations. Monitor ECG and arterial blood pressure continuously throughout the procedure and during recovery.

Prior to removing a catheter/lead that has an extravascular portion, it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

When advancing sheath(s) over the catheter/lead, maintain adequate tension on the catheter/lead (via the basket snare or Needle's Eye® Snare) to guide the sheath(s) and prevent damage to vessel walls.

If excessive scar tissue or calcification prevents safe advancement of sheath(s), consider an alternate approach.

Excessive force with sheath(s) used intravascularly may result in damage to the vascular system requiring surgical repair.

If the catheter/lead breaks, evaluate fragment; retrieve as indicated.

If hypotension develops, rapidly evaluate; treat as appropriate.

Due to rapidly evolving catheter/lead technology, this device may not be suitable for the removal of all types of catheters/leads. If there are questions or concerns regarding compatibility of this device with particular catheters/leads, contact the catheter/lead manufacturer.

ADVERSE EVENTS

Potential Adverse Events

Potential adverse events related to the procedure of intravascular extraction of catheters/leads include (listed in order of increasing potential effect):

- dislodging or damaging nontargeted catheter/lead
- chest wall hematoma
- thrombosis
- arrhythmias
- acute bacteremia
- acute hypotension
- pneumothorax
- stroke
- migrating fragment from catheter/object
- pulmonary embolism
- laceration or tearing of vascular structures or the myocardium
- hemopericardium
- cardiac tamponade
- hemothorax
- death

INSTRUCTIONS FOR USE

Suggested Instructions for Use:

Byrd Workstation™

NOTE: In addition to closely following the “Suggested Instructions for Use” for this device, refer to the “Clinical Considerations for Removing Leads, Catheters, or Other Indwelling Objects: Femoral Approach” shown in the next section of this insert.


NOTE: The Tip Deflecting Wire Guide is intended for use only to position the indwelling object for snaring. Excessive force may damage the Tip Deflecting Wire Guide.

CAUTION: Competency with the Seldinger technique is necessary. Care should be taken to avoid the artery and to enter the vein cleanly and at a shallow angle.

CAUTION: Use fluoroscopic monitoring during all sheath, catheter/lead, or snare manipulations.

CAUTION: Prior to removing a catheter/lead that has an extravascular portion, it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

1. After obtaining access to the vessel using the Seldinger technique and inserting the 16 Fr introducer sheath under fluoroscopic guidance, inspect the pre-loaded 12 Fr sheath. The Tip Deflecting Wire Guide should cross through or be alongside the Dotter Basket. Note their relative positions and make sure the Tip Deflecting Wire Guide extends beyond the Dotter Basket. Retract the Tip Deflecting Wire Guide and Dotter Basket into the 12 Fr sheath.
2. Attach the Deflecting Handle to the Tip Deflecting Wire Guide.
3. Under fluoroscopic guidance, introduce the pre-loaded 12 Fr sheath through the 16 Fr sheath and advance it beyond the end of the 16 Fr sheath approximately 5 to 6 centimeters below the right atrium.
4. Follow Option A or B to snare the indwelling object to be removed, then proceed to step 5. Option B may be preferable for larger leads or defibrillator leads. It also allows the ability to gradually pull the subclavian end of the lead down to the atrium while restricting the tension below the level that would apply traction to the myocardium.

Option A

Figure 1


Deploy the basket and Tip Deflecting Wire Guide into the right atrium by carefully advancing each from the 12 Fr femoral sheath. Deflect the Tip Deflecting Wire Guide around the indwelling object within the atrium, encircling the object (A). When the Tip Deflecting Wire Guide has encircled the lead, lock the Tip Deflecting Wire Guide in place by tightening the locking screw on the handle.

Use the Tip Deflecting Wire Guide to stabilize the indwelling object. Slide the helical loop basket forward so that a portion of the object's body is positioned within the confines of the helical loop basket (B).

To close the basket, do not withdraw it. Instead, while maintaining the position of the indwelling object, basket, and Tip Deflecting Wire Guide, advance the 12 Fr sheath forward to close the basket around the object (C). Releasing and slightly withdrawing the Tip Deflecting Wire Guide may facilitate advancing the sheath over the basket.

If free end of the indwelling object is still extravascular, cut it as short as possible.

When the lead is captured, cautiously pull the entire Workstation™ downward so that the proximal portion of the lead is drawn into the vasculature (D).

Option B

Figure 2

After deploying the helical loop basket and Tip Deflecting Wire Guide, deflect the Tip Deflecting Wire Guide to encircle the indwelling object, positioning the distal end of the Tip Deflecting Wire Guide past and perpendicular to the shaft of the Tip Deflecting Wire Guide (A). Lock the lip Deflecting Wire Guide in position by tightening the single locking screw on the handle.

Slide the helical loop basket forward and capture the distal end of the lip Deflecting Wire Guide in the basket (B).

Close the basket over the Tip Deflecting Wire Guide by carefully advancing the 12 Fr sheath, creating a strong, reversible loop around the body of the indwelling object (C).

If free end of the indwelling object is still extravascular, cut it as short as possible.

Slide the complete femoral Workstation™ downward to draw the proximal portion of the indwelling object into the vasculature and into the inferior vena cava (D).

With the proximal end of the indwelling object in the right atrium or inferior vena cava, release the loop by 1) loosening the locking screw on the handle of the Tip Deflecting Wire Guide and 2) opening the basket by pushing it forward. Position the opened basket below the proximal end of the indwelling object (E). Position the proximal end of the indwelling object in the basket by slowly advancing the basket (F). Slight rotation of the basket may be helpful.

Once a portion of the indwelling object passes through the basket, advance the 12 Fr sheath forward to close the basket around the indwelling object (G).

5. After the indwelling object has been snared, it is freed from entrapments by advancing the sheath(s) over the object while maintaining adequate tension on the snared object. Refer to Part III of the **"Clinical Considerations for Removing Leads, Catheters, or Other**

Indwelling Objects: Femoral Approach” shown in the next section of this insert.

Clinical Considerations for Removing Leads, Catheters or Other Indwelling Objects:

FEMORAL APPROACH

Clinical experience in the removal of leads has identified several considerations for lead removal techniques using superior or femoral approaches. Physicians highly experienced with lead removal techniques have suggested the following considerations for removal of leads, catheters, or other indwelling objects using a femoral approach.

Part I. Accessing the Vessel

CAUTION: Competency with the Seldinger technique is necessary. Care should be taken to avoid the artery and to enter the vein cleanly and at a shallow angle.

CAUTION: Use fluoroscopic monitoring during all sheath, catheter/lead, or snare manipulations.

CAUTION: Prior to removing a catheter/lead that has an extravascular portion, it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

1. Insert the dilator into the 16 Fr sheath (outer sheath). Use the Seldinger technique for sterile insertion of the syringemounted needle and wire guide into the femoral vein. Being certain that the dilator leads the sheath, insert the 16 Fr introducer sheath and dilator under fluoroscopic guidance. Remove the dilator and advance the sheath until it extends to the middle third of the inferior vena cava.

Part II. Snaring the Indwelling Object

2. Follow the “Suggested Instructions for Use” to: advance the 12 Fr sheath containing the snare to be used (Tip Deflecting Wire Guide with Dotter Basket, or the Needle’s Eye® Snare) inside the 16 Fr introducer sheath deploy the snare, and snare the catheter/lead or other indwelling object.

Part III. Removing the Indwelling Object

3. If removing an active fixation lead, rotate the snare (helical loop basket or Needle’s Eye® Snare) counterclockwise to unscrew the fixation mechanism.

4. While maintaining adequate tension on the indwelling object via the closed snare (basket or Needle’s Eye®Snare), prepare the indwelling object for removal by advancing the 16 Fr sheath over the object to encapsulate it. Cautiously disrupt fibrotic adhesions encountered during sheath advancement. (Where possible, the 12 Fr sheath can also be advanced to aid in dilation).

5. For a cardiac lead, if the lead has not been freed by the time the sheath nears the myocardium, position the 16 Fr sheath about one centimeter from the myocardium, making sure that the 12 Fr sheath has been withdrawn into the 16 Fr sheath at least 6 cm. As the 16 Fr sheath is held in position by stabilizing the sheath at its venous entry site in the groin, apply countertraction by adding tension to the snare (basket or Needle’s Eye®Snare), thereby pulling the lead tip into the sheath and drawing the endocardium against the 16 Fr sheath’s distal tip. Steadily pull on the lead until the lead tip is freed from the myocardium and is pulled into the sheath.

6. When the indwelling object has been freed from tissue entrapment and retracted into the 16 Fr sheath, withdraw the sheath back into the inferior vena cava. Remove the entire ensemble while maintaining hemostasis.

NOTE: If removing multiple catheters/leads, leave the 16 Fr sheath in place until all catheters/leads have been extracted. If using the Tip Deflecting Wire Guide and Dotter Basket, a new Tip Deflecting Wire Guide, Dotter Basket, and 12 Fr sheath should be used to extract subsequent leads.

NOTE:To prevent potential complications, carefully examine any broken catheters/leads. A long exposed lead coil could lacerate the heart, or a loose fragment could migrate, possibly into the pulmonary artery. A pigtail catheter or closed loop snare may be used to retrieve a catheter/lead from the pulmonary artery, repositioning it in the ventricle for removal.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.