COOK®

MEDICAL

CE
0088

EN
2

**Trigeminal Ganglion
Microcompression**

Instructions for Use

DA
14

**Mikrokompression
af ganglion trigeminale**

Brugsanvisning

DE
26

**Mikrokompression des Ganglion
trigeminale**

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EL
38

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F M _ 2 7 4 9 _ R E V 0

MULLAN PERCUTANEOUS TRIGEMINAL GANGLION MICROCOMPRESSION SET

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Mullan Percutaneous Trigeminal Ganglion Microcompression Set (MTNS-) consists of an access needle with stylet, a blunt obturator and a microcompression balloon catheter.

The set is available with the following order number suffixes: -MN: containing a beveled access needle and stylet, -GN: containing a blunt access needle with pencil point stylet, -BN: containing a blunt access needle with beveled stylet and small probing stylets.

All sets also include a blunt obturator, scalpel, syringe, adhesive bandage, 22 ga. needle, angled needle gauge and gauze pad.

Refer to label for further information.

All Tables and References are located in document FM-2741 which accompanies this Instruction for Use.

INTENDED USE

The Mullan Percutaneous Trigeminal Ganglion Microcompression Set is intended for use in the percutaneous treatment of trigeminal neuralgia. This set is supplied sterile in a peel-open package and is intended for one-time use.

PATIENT CONSIDERATIONS

- Percutaneous surgical treatment of symptomatic Trigeminal Neuralgia (TN) may be considered for patients who are intolerant of or refractory to pharmaceutical therapy such as carbamazepine and/or diphenhydantoin.
- Postoperative numbness is to be expected and is usually well tolerated. Patients who are unwilling to accept facial sensory loss should consider posterior fossa surgery for microvascular decompression.
- Technical difficulties in needle placement and subsequent balloon compression have been reported in patients with impaired access to the foramen ovale (References 1,3,5,6,8) such as those with Paget's disease or those who have undergone a previous treatment of alcohol injection. Alternative treatment may be considered if impaired access to the foramen ovale is suspected.
- Clinical studies have shown that 90 to 100 percent of patients treated with balloon compression experience initial pain relief. (References 1, 2, 3, 5, 6, 8)
- Greatest success with balloon compression has been reported for patients who have no prior surgical intervention for TN. (Reference 8)
- Patients with recurrent TN after glycerolization may respond well to balloon compression if access to the foramen ovale is not impaired. (Reference 8)
- Recurrent TN after balloon compression may be treated with additional balloon compression, however, some studies suggest that treatment success rates may be lower. (Reference 8)
- Recurrence of pain over time has been reported with this and other surgical treatments for TN. (References 1, 2, 3, 5, 6, 8) Table 1 provides a summary of clinical data published concerning the outcome of balloon compression for TN for the described patient population. (References 1, 2, 4, 6, 8)
- Balloon compression has been reported as an effective nonselective therapy for TN affecting the 1st and 2nd divisions. (References 1, 2, 4, 5, 6, 8)
- For patients with isolated 3rd division trigeminal pain, radiofrequency thermocoagulation has been proposed as a preferred percutaneous treatment. (Reference 8)

Table 2 provides a comparative summary of published clinical results for the treatment of TN with radiofrequency thermocoagulation, glycerolization, and balloon compression. (Reference 8)

CONTRAINDICATIONS

There are no known contraindications of this device in the treatment of TN.

WARNINGS

- Contrast agents should not be used in patients with known or suspected contrast agent allergies or other adverse reactions.
- Manipulation of products requires fluoroscopic control.
- The possibility of balloon rupture must be taken into account when considering the risks involved in the use of this device in the treatment of TN.
- Upon engagement of the needle in the foramen ovale and during balloon inflation, bradycardia of short duration with associated hypotension may occur. In patients with heart disorders it is advisable to guard against bradycardia by administration of atropine. An external pacemaker has also been cited as useful in the temporary treatment of bradycardia in some patients. (References 1, 2, 4, 5, 6, 8)
- Air should not be used for balloon inflation.
- Hypertension has been reported in some patients (Reference 4) at the time of balloon inflation. Blood pressure changes and possible cardiac arrhythmias during the hypertensive response may be controlled by administration of antihypertensive drugs. Upon cessation of balloon compression, cardiovascular changes have been reported to return to normal within 3 to 10 minutes. (Reference 4)
- Pulsatile blood from the needle may indicate penetration of a significant vascular structure. The procedure should be aborted immediately.
- Arterial bleeding may arise from penetration of the carotid artery or an accessory meningeal artery traversing the foramen ovale. Care must be taken to screen patients with possible anatomical abnormalities that place the carotid artery in or near the foramen ovale (Reference 4). Patients with primitive foramen lacerum medius may be at greater risk for such complications. (Reference 4)
- Sinus arrest has been reported in one patient (Reference 2) upon engagement of the needle in the foramen ovale. Reversal without complications was achieved by administration of atropine in conjunction with chest compressions.

PRECAUTIONS

- Manipulation of products requires fluoroscopic control.
- Care must be taken during the placement of the needle to ensure that the needle does not penetrate through the inside of the cheek and into the mouth, as the risk of infection is significantly increased.
- Care must be taken to ensure that the balloon catheter is not extended further than necessary as entry into the posterior fossa is possible.
- Longer inflation times have been reported to reduce the recurrence rate of TN; however, longer inflation times have also been reported to increase the risk of significant sensory loss. (References 1, 2, 4, 5, 6, 8)
- Care must be taken to ensure that the balloon catheter is not damaged by the needle. Once the balloon catheter has been placed within the needle, do not attempt to pull the catheter back through the needle. If repositioning is required, remove the needle and catheter as a unit, then carefully remove the catheter from the needle under direct vision and verify via balloon inflation and visual examination that the catheter is not damaged.
- No single inflation time has been reported as optimal for all patients. Some degree of permanent nerve damage has been reported for compression times of 4 to 10 minutes, while others report that shorter inflation times of approximately 1 minute yield acceptable results with fewer complications. (References 1, 4, 5, 6, 8)
- As with any surgical procedure, care must be taken to reduce the risk of infection via strict adherence to aseptic procedures

- Use of non-diluted contrast agents, catheter kinking, or over tightening of the Tuohy-Borst adapter may result in partial or total occlusion of the catheter lumen preventing balloon inflation or deflation. In the event of slow (or no) balloon deflation, replacement of the balloon catheter stylet may reopen the catheter lumen.
- No complications as a result of balloon rupture during the course of compression have been reported; however, premature balloon deflation may result in poor treatment outcome requiring additional treatment.
- Reports indicate that, regardless of the treatment procedure (radiofrequency thermocoagulation, glycerolization, or balloon compression) some sensory loss is likely, a balance between greater sensory loss and lower recurrence rate must be determined.

Table 3 has been added to describe reported complications associated with balloon compression for TN and information to minimize these events.

ADVERSE EFFECTS

- Clinical observations noted in literature which represent possible adverse reactions and complications include both minor and severe numbness, and minor dysesthesia in some patients. (References 1, 2, 4, 5, 6) Most patients readily adjust to minor sensory deficits. Constant, severe dysesthesias in an anesthetic or analgesic zone (anesthesia dolorosa) rarely occur with this procedure.
- Eye complications have rarely been reported. Literature states that neurolytic keratitis occurs in less than 2% of patients. However, care must be taken to protect against corneal ulceration in these patients.
- As with any surgical procedure, venous or arterial bleeding is possible during needle placement. Hematoma with transient discomfort have been noted in some patients.
- Transient diplopia has also been noted as a rare complication which resolved in all patients within several weeks.

POTENTIAL ADVERSE EFFECTS

- Other possible adverse effects which may occur as a result of this procedure include hemorrhage, meningitis, loss of corneal sensation, paresis of the muscles innervated by the motor root of the trigeminal nerve including weakness of the masseter, temporalis, and pterygoid muscles. Muscle weakness could result in jaw deviation and loss of chewing power. Difficulty in hearing has also been reported in similar (radiofrequency) treatments related to roaring and popping sounds due to paresis of the small muscles about the eustachian tube and tympanic membrane. Difficulty with mastication without motor paresis is also possible.

SUGGESTED INSTRUCTIONS FOR USE

Patient and Equipment Preparation

1. The patient is typically prepared under general anesthesia.
2. Intubation is recommended to ensure maintenance of airway.
3. Use of fluoroscopy is recommended.
4. Prepare about 3 ml of nonionic contrast agent. A 50 percent dilution with sterile saline will reduce viscosity for easier balloon filling and emptying.
5. With the patient in supine position, induce anesthesia and prepare the cheek of the patient. Care should be taken to ensure that the intended puncture site, typically 2.5 cm external to and slightly above the lateral angle of the mouth, is prepared.
6. Ensure that balloon catheter air space is evacuated prior to use. This is easily accomplished by replacing the air with contrast agent. After the air in the balloon and catheter has been replaced with contrast agent, close the stopcock on the balloon catheter.


Fig. 1

1. Use of the scalpel rather than direct needle puncture is optional. If using the scalpel included with the set, make a small nick in the skin of the patient's cheek at the site of entry of the needle approximately 2.5 cm external to and slightly above the lateral angle of the mouth (Fig. 1).


Fig. 2

2. Insert the access needle and engaged sharp stylet through the skin and tissues of the cheek under fluoroscopic control until engaging the foramen ovale. This is done by inserting the needle laterally into the mobile cheek for a short distance before directing it to the foramen magnum. The patient may require repositioning along with repositioning of the fluoroscopy unit (typically C-Arm) to ensure proper trajectory and engagement of the foramen ovale. Care must be taken to ensure that the needle does not penetrate through the inside of the cheek and into the mouth, increasing the risk of infection (Fig. 2).
3. The needle should not penetrate beyond the foramen, but should remain engaged in the foramen.
4. After insertion of the needle through the cheek and facial tissues, the sharp stylet may be removed and replaced with the blunt stylet included in the set. The blunt stylet may then be used to probe for the foramen.

WARNING: Upon engagement of the foramen, bradycardia of short duration may occur.

5. Upon engagement of the foramen and fluoroscopic verification, advance the balloon catheter through the needle.
6. The balloon catheter should be advanced until the balloon extends just beyond the tip of the needle within Meckel's Cave. The two circumferential marks on the balloon catheter, when aligned with the top of the needle female Luer fitting, indicate when the tip of the balloon is at the tip of the needle (most distal mark) and when the balloon is completely outside the needle (most proximal mark) (Fig. 3).

CAUTION: Care must be taken to ensure that the balloon catheter is not extended further than necessary as entry into the posterior fossa is possible.

Fig. 3


Fig. 4



Fig 5

7. To prevent the balloon from entering further than necessary, the catheter is fitted with a Tuohy-Borst adapter that may be placed at the desired location on the catheter and tightened to provide a positive stop against the needle hub. Tighten the adapter by grasping the knurled portion in one hand and turning the female luer adapter clockwise. Do not over tighten as occlusion of the catheter lumen can occur (Fig. 4).

8. Verify placement of the balloon by visualizing the catheter under fluoroscopy.

CAUTION: If repositioning is required, remove the catheter and needle as a unit. Do not pull the catheter back through the needle as damage to the balloon or catheter may occur.

9. Slowly inflate the balloon under fluoroscopic control using the 1 cc syringe included with the set.

CAUTION: Care must be taken to slowly inflate the balloon and to not exceed 1cc of liquid volume. Typically 0.7 to 1 ml of contrast/saline solution is required.

WARNING: During balloon inflation, bradycardia of short duration may occur. The balloon typically takes on a pear shape as Meckel's Cave is filled and the trigeminal nerve is compressed. (Figure 5)

Note: Repeated balloon failures may occur in certain patients due to varying anatomy of the sphenoid bone, foramen ovale or within Meckel's Cave. Alternative treatments should be explored in these cases.

The pear shape is indicative of an appropriate amount of compression of the nerve.

10. The balloon is typically inflated for 1 minute. The fluid volume may be maintained by use of the catheter stopcock or the 1 cc syringe. Patients with recurrent or difficult to control symptoms may require longer inflation times to achieve pain control.

CAUTION: Longer inflation times may increase the risk of severe numbness or dysesthesia.

11. Deflate the balloon opening the catheter stopcock and withdrawing the syringe plunger.
12. Once the balloon has deflated, withdraw the balloon catheter and needle together.

CAUTION: Do not attempt to withdraw the balloon through the needle as the needle may damage the catheter.

13. Firm digital pressure should be applied at the site of needle entry to prevent hematoma.

Post-operative Care

1. Upon completion of the procedure, the patient's cheek should be cleansed. The small adhesive bandage included with the set may be placed over the site of the needle entry.
2. Upon recovery from general anesthesia a complete neurological examination should be performed to ensure that the patient is free from pain or complications.
3. Specific medications for control of TN are usually discontinued; however, if long-term high dose therapy has been employed the medications should be tapered to avoid withdrawal symptoms.
4. Postoperative headache is not uncommon. Mild analgesics are routinely used to control pain.
5. Neurological examination with special care taken to evaluate the extent of pain relief and numbness should be undertaken.
6. In the event of reduced corneal sensation, meticulous eye care including the use of artificial tears may be instituted.
7. In the event of reduced oral sensation and/or motor function a temporary soft diet may be instituted.
8. Sensory reductions are not uncommon and are expected as a result of balloon compression. Patients should be made aware of this aspect of the treatment and informed that sensory changes are reported to resolve or lessen over time. (References 1, 2, 4, 5, 6, 8) However, some sensory reductions may persist indefinitely.
9. Patients with suspected reduced corneal sensation should be advised to observe the eye for redness and avoid rubbing or touching.
10. Postoperative use of ice packs and/or pressure to the jaw have been reported to reduce pain, swelling, and hematoma formation. (References 1, 2, 4, 5, 6, 8)

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

If device packaging is damaged, do not use and discard the device.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

MULLAN PERCUTANEOUS TRIGEMINAL GANGLION MICROCOMPRESSION SET

Order Numbers: MTNS-4.0-40-MN
MTNS-4.0-40-GN
MTNS-4.0-40-BN

COMPONENTS


TABLE 1 - CLINICAL DATA SUMMARY

FACTOR	MULLAN (1983)
Study Duration [yrs]	5
# Pts	50
Age [yrs]	16-88
# Males	22
# Females	28
% w/Previous Treatments	48%
% w/Initial Relief	98%
Overall Complications	
Cardiac Arrest	---
Sinus Bradycardia	7 ^a
Hypertension	---
Cardiac Arrhythmia	---
Difficulty in Access	10%
Dysesthesia	12%
Hypesthesia	8%
Anesthesia	---
Anesthesia Dolorosa	---
Transient Motor Weakness	7 ^a
Persistent Motor Weakness	---
Transiently Diminished Corneal Reflex	---
Bleeding	---
Herpetic Eruption	---
Diplopia	---
Cheek Hematoma	---
Headache	---
Aseptic Meningitis	---
Transient Palsy	2%
Length of Follow-Up [yrs]	0.5-4.5
% w/Recurrence of Pain	12%
Subsequent Treatment Resulting in Lasting Relief	
Repeat Balloon Compression	6%
Partial Posterior Fossa Section	2%
Posterior Fossa Decompression	2%
Low-Dose Carbamazepine	2%
RF Coagulation	---
Microvascular Decompression	---

--- *No occurrence noted*

^a *Occurrence noted, but number/incidence not specified*

^b *4 additional patients had pain recurrence after a second balloon compression. For these patients, as well as for the remaining 4 patients having recurrence, further treatment resulting in relief of pain was not noted.*

TABLE 2 - CLINICAL DATA SUMMARY

FACTOR
Pts
Avg Length of Follow-Up [yrs]
w/Pain Relief
w/Recurrence of Pain
Subsequent Treatment Resulting in Lasting Relief
RF Thermocoagulation
Glycerolization
Balloon Compression
Post-operative Complications
Hypesthesia (Unwanted)
Analgesia or Anesthesia w/disappearance of corneal reflex
Keratitis
Anesthesia Dolorosa
Paresthesia Requiring Treatment
Transient Masticatory Weakness
Persistent Masticatory Weakness
Transient Nerve Palsy

--- *No occurrence noted*

FOR BALLOON COMPRESSION OF THE TRIGEMINAL NERVE IN THE TREATMENT OF TRIGEMINAL NEURALGIA			
BELBER (1987)	LOBATO (1990)	BROWN (1993)	FRAIOLI (1989)
7	4.5	7.5	10
25	144	50	159
48-86	83.3% > 50	35-86	95% > 60
15	58	7 ^A	7 ^A
10	86	7 ^A	7 ^A
60%	43%	52%	8%
100%	93.8%	94%	89.9%
---	---	---	---
4%	---	---	---
4%	---	7 ^A	---
---	7 ^A	---	---
---	7 ^A	---	---
---	20%	---	0.6%
6%	19%	20%	6.9%
---	4.1%	---	---
---	---	---	2.5%
---	---	---	0.6%
48%	12%	7 ^A	6.9%
---	---	---	3.1%
---	7 ^A	---	---
---	13.9%	---	---
48%	11%	---	---
---	2.8%	---	---
16%	3.5%	---	---
---	2.8%	---	---
---	---	6%	---
---	---	2%	---
0.5-7	0.5-4.5	0.75-7	Avg. 3.5 yrs.
24%	9.7%	26%	9.8%
---	---	---	---
---	7.6%	8% ^A	5%
---	---	---	---
---	---	---	---
24%	1.4%	---	---
---	0.7%	---	1.8%
---	---	2%	---

(FRAIOLI, ET AL. 1989) FOR COMPARATIVE STUDY OF RF THERMOCOAGULATION, GLYCEROLIZATION, AND BALLOON COMPRESSION FOR TREATING TRIGEMINAL NEURALGIA		
RF Thermocoagulation	PROCEDURE Glycerolization	Balloon Compression
533	32	159
6.5	5	3.5
519 (97.4%)	23 (71.9%)	143 (89.9%)
52 (10%)	13 (56.5%)	14 (9.8%)
---	---	---
52	11	3
---	---	---
---	10	8
---	---	---
92 (17.3%)	---	---
16 (3%)	1 (3.1%)	4 (2.5%)
10 (1.9%)	1 (3.1%)	---
8 (1.5%)	---	1 (0.6%)
81 (15.2%)	3 (9.4%)	11 (6.9%)
---	---	11 (6.9%)
16 (3%)	---	5 (3.1%)
1 (0.2%)	---	---

TABLE 3**COMPLICATIONS (INCIDENCE)**

Cardiac Arrest (4%)

Hypertension (Not Specified)

Cardiac Arrhythmia (4%)

Difficulty in Access (0.8-20%)

Dysesthesia (6-20%)

Hypesthesia (2.5-8%)

Transient Motor Weakness w/Asymmetrical Mastication (6.9-48%)

Persistent Motor Weakness (3.1%)

Transiently Diminished Corneal Reflex (3.1%)

Bleeding (13.9%)

Herpetic Eruption (11-48%)

Diplopia (2.8%)

Cheek Hematoma (3.5-16%)

Headache (2.8%)

Aseptic Meningitis (6%)

Transient Palsy (2%)

TABLE 4**ANALYSIS OF REPORTED RESULTS
OF BALLOON COMPRESSION
(REFERENCES 1, 2, 4, 6, 8)**

Excellent	64.0%
Good	16.9%
Fair	12.9%
Poor	2.0%
Failure	4.9%

Note: Results in Table 4 are mean values determined from the above references. Therefore, total percentages add to slightly over 100%.

Excellent results are defined as no TN pain and no postoperative complications of dysesthesia, hypesthesia, anesthesia, motor weakness (transient or persistent), reduced corneal reflex, anesthesia dolorosa, transient palsy, or diplopia.

Good results are defined as no TN pain but complications of minor hypesthesia, anesthesia, and transient motor weakness.

Fair results are defined as no TN pain but complications of paresthesia requiring treatment, reduced corneal reflex, and dysesthesia.

Poor results are defined as no TN pain but complications of persistent motor weakness, anesthesia dolorosa, diplopia, and transient palsy.

Failure is defined as patients who receive no relief of TN immediately post-treatment.

TABLE 5**THERAPY RESULTING IN LASTING
RELIEF OF TN AFTER RECURRENCE
(REFERENCES 1, 2, 4, 6, 8)**

Repeat Balloon Compression	42.0%
Drug Therapy	26.2%
Microvascular Decompression	8.2%
Radiofrequency Thermocoagulation	5.1%
Unreported	18.4%

REPORTED BALLOON COMPRESSION COMPLICATIONS

COMMENT

Continuously monitor heart rate throughout procedure, and administer atropine if necessary. An external pacemaker may be used.

Hypertension may occur in response to compression, and last throughout the inflation duration. Continuously monitor blood pressure, maintain appropriate anesthetic level, and administer anti-hypertensive agents if necessary.

Continuously monitor heart rate throughout procedure. An external pacemaker may be used.

Only experienced, skilled physicians should perform this procedure. If impaired access to the foramen ovale is suspected (such as in patients with Paget's disease, or those who have undergone a previous treatment of alcohol injection), consider alternative treatment.

Sensory reductions are not uncommon, and may resolve or lessen over time. Determine balance between decreased sensory loss and increased rate of recurrence and adjust balloon compression time accordingly.

Sensory reductions are not uncommon, and may resolve or lessen over time. Determine balance between decreased sensory loss and increased rate of recurrence and adjust balloon compression time accordingly.

A soft diet is recommended.

A soft diet is recommended.

Meticulous eye care, including the use of artificial tears, is recommended. Advise patient to observe eye for redness and avoid rubbing or touching.

Careful technique is required. Screen patients with possible anatomical abnormalities that place the carotid artery in or near the foramen ovale. Needle repositioning may be necessary.

Careful technique is required. Placing a finger inside the mouth may facilitate placement.

Diplopia usually resolves over time. Its occurrence may be minimized by ensuring the inflated balloon is kept in Meckel's Cave.

Careful technique is required. Use of ice packs and/or pressure may reduce hematoma formation.

When post-operative headaches occur mild analgesics may be used to control pain.

Careful technique is required. Placing a finger inside the mouth may facilitate placement, and prevent needle penetration through the oral mucosa.

May be related to dysesthesia. Determine balance between decreased sensory loss and increased rate of recurrence and adjust balloon compression time accordingly.

REFERENCES

1. S. Mullan, T. Lichtor: "Percutaneous Microcompression of the Trigeminal Ganglion for Trigeminal Neuralgia," *Journal of Neurosurgery*, 59 (1983), 1007-1012.
2. C. J. Belber, R. A. Rak: "Balloon Compression Rhizolysis in the Surgical Management of Trigeminal Neuralgia," *Neurosurgery*, 20 (1987), 908-913.
3. A. M. Gerber: "Needle for Use During Percutaneous Compression of the Gasserian Ganglion for Trigeminal Neuralgia," *Journal of Neurosurgery*, 71 (1989), 455.
4. R. D. Lobato, J. J. Rivas, R. Sarabia, E. Lamas: "Percutaneous Microcompression of the Gasserian Ganglion for Trigeminal Neuralgia," *Journal of Neurosurgery*, 72 (1990), 546-553.
5. T. Lichtor, J. F. Mullan: "A 10-Year Follow-up Review of Percutaneous Microcompression of the Trigeminal Ganglion," *Journal of Neurosurgery*, 72 (1990), 49-54.
6. J. A. Brown, M. D. McDaniel, M. T. Weaver: "Percutaneous Trigeminal Nerve Compression for Treatment of Trigeminal Neuralgia: Results in 50 Patients," *Neurosurgery*, 32 (1993), 570-573.
7. A. M. Gerber: "Improved Visualization of the Foramen Ovale for Percutaneous Approaches to the Gasserian Ganglion," *Journal of Neurosurgery*, 80 (1994), 156-159.
8. B. Fraioli, V. Esposito, B. Guidetti: "Treatment of Trigeminal Neuralgia by Thermocoagulation, Glycerolization, and Percutaneous Compression of the Gasserian Ganglion and/or Retrogasserian Rootlets: Long-term Results and Therapeutic Protocol," *Neurosurgery*, 24 (1989), 239-244.

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1. S. Mullan, T. Lichtor: "Percutaneous Microcompression of the Trigeminal Ganglion for Trigeminal Neuralgia", *Journal of Neurosurgery*, 59 (1983), 1007-1012.
2. C. J. Belber, R. A. Rak: "Balloon Compression Rhizolysis in the Surgical Management of Trigeminal Neuralgia", *Neurosurgery*, 20 (1987), 908-913.
3. A. M. Gerber: "Needle for Use During Percutaneous Compression of the Gasserian Ganglion for Trigeminal Neuralgia", *Journal of Neurosurgery*, 71 (1989), 455.
4. R. D. Lobato, J. J. Rivas, R. Sarabia, E. Lamas: "Percutaneous Microcompression of the Gasserian Ganglion for Trigeminal Neuralgia", *Journal of Neurosurgery*, 72 (1990), 546-553.
5. T. Lichtor, J. F. Mullan: "A 10-Year Follow-up Review of Percutaneous Microcompression of the Trigeminal Ganglion", *Journal of Neurosurgery*, 72 (1990), 49-54.
6. J. A. Brown, M. D. McDaniel, M. T. Weaver: "Percutaneous Trigeminal Nerve Compression for Treatment of Trigeminal Neuralgia: Results in 50 Patients", *Neurosurgery*, 32 (1993), 570-573.
7. A. M. Gerber: "Improved Visualization of the Foramen Ovale for Percutaneous Approaches to the Gasserian Ganglion", *Journal of Neurosurgery*, 80 (1994), 156-159.
8. B. Fraioli, V. Esposito, B. Guidetti: "Treatment of Trigeminal Neuralgia by Thermocoagulation, Glycerolization, and Percutaneous Compression of the Gasserian Ganglion and/or Retrogasserian Rootlets: Long-term Results and Therapeutic Protocol", *Neurosurgery*, 24 (1989), 239-244.

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Bloomington, IN 47404 U.S.A.